

DEG602 BACK MOTOR GEAR INSTALLATION INSTRUCTIONS

- **1.** Raise chair and seat motor to highest position. Disconnect AC power to chair.
- 2. Remove screws holding bellows top cover to bellows hangers and lower cover. Remove plastic motor skirt cover.
- **3.** Remove back link (bent-bar) from crank by unscrewing round-head shoulder screw at back motor bracket. Hold back of chair from falling.
- **4.** Disconnect motor from chair junction box.
- **5.** Remove socket-head cap screw, lock washer and nut from rear of back motor bracket.
- **6.** Remove two cap screws holding back motor and bracket assembly to cradle assembly.

CAUTION When removing cap screws, support back motor and bracket assembly.


- 7. Remove bracket assembly from cradle.
- **8.** Loosen set screws in yoke of back motor bracket. Remove motor mounting trunion pins.
- 9. Remove snap rings from screw tube retaining pin. Remove pin.

CAUTION Do not move motor adjustment screws.

- **10.** Remove cap screws at each corner of limit switch cover on back motor. Remove cover.
- **11.** Remove screws from gear housing cover. Pull motor shaft slightly, separating cover from housing.
- **12.** Unscrew motor screw shaft and gear assembly from actuator, pulling motor gear shaft assembly out of housing.


***NOTE REGARDING ACTUATOR**

To determine which actuator you need, measure thread size and direction of thread. For actuators to fit the DEN-TAL-EZ back motor, please see the RPI catalog.



Motor Screw Shaft and Gear Assembly Alignment

- **13.** Remove roll pin from gear shaft assembly, pull gear off shaft.
- **14.** Place spacer on shaft before installing replacement gear (RPI Part #DEG602). Install gear on shaft. Line up pin holes and install roll pin.
- **15.** Place spacer on top of gear. Install roller bearing then spacer. Install complete screw shaft and gear assembly into gear housing.
- **16.** Screw fine-threaded end of shaft into actuator. Stop when actuator is between limit switches.
- **17.** Push screw shaft into gear housing, twisting shaft to make gear mesh with worm drive gear. Installed properly, gear housing cover will be firmly seated against gear housing.
- **18.** Replace screws in gear housing cover and limit switch cover.
- **19.** Connect motor to chair junction box. Connect AC power to chair.
- **20.** Press top of the Back Control switch so screw tube will turn counter-clockwise and run motor to its full "DOWN" position. Disconnect AC power to chair.
- **21.** Allow 7/16["] maximum clearance between edge of screw tube and edge of gear housing. Adjust screw tube as necessary.



- 22. Replace motor in bracket assembly with motor pins and screw tube pin.
- **23**. Secure motor and bracket assembly to cradle assembly with screws.

CAUTION

The socket-head screw cap should be installed from under the bracket assembly.

- **24**. Attach the back link (bent-bar) to bracket crank with round-head shoulder screw.
- **25**. Connect back motor to chair junction box.
- **26**. Replace plastic motor skirt with screws.
- **27**. Raise bellows cover and fasten it to bellows hangers with screws.
- **28**. Connect chair to AC power and check for proper operation.