



replacement parts industries, inc.

CALL 800-221-9723 or 818-882-8611 • FAX 818-882-7028  
E-MAIL: techsupport@rpiparts.com • WEBSITE: www.rpiparts.com

## SCK020 AIR COMPRESSOR KIT INSTALLATION INSTRUCTIONS

### **SERVICE TIP: REPLACE THE CHECK VALVE WHEN REPLACING THE COMPRESSOR!**

A leaky check valve can cause moisture or steam to seep inside a compressor thus causing the compressor to fail. We recommend replacing the check valve at the same time you replace the compressor so, RPI includes a Check Valve Kit (RPI Part #SCK011) with the AIR Compressor Kit (RPI Part #SCK020). It's the RPI Advantage!

- 1) Turn off and unplug the autoclave.
- 2) Remove the tubing connecting the check valve to the compressor.
- 3) Use a 7/16" wrench to secure fitting on steam generator/boiler and remove the check valve using a 1/2" wrench.
- 4) After the check valve is removed, clean any debris from the steam generator/boiler fitting. **Warning: Incomplete removal of debris from the steam generator/boiler fitting before installing new check valve may result in valve failure.**
- 5) Apply approximately 3-4 drops of the enclosed Threadlocker 545 (RPI Part #RPA369) to completely cover the threads on the steam generator/boiler fitting.
- 6) Install new check valve hand-tight onto boiler fitting. Use a 7/16" wrench to secure the boiler fitting and tighten the new check valve with a 1/2" wrench (approximately 3-4 turns).

### **CAUTION:**

Care must be exercised to ensure that fittings on the steam generator/boiler are properly wrenched when loosening or tightening mating fittings. Failure to do so may result in damage to the fitting on the steam generator/boiler. If damage occurs, it will be an expensive repair that can only be made at the OEM's repair facility. There are no "field serviceable" options for this type of damage.


- 7) Note the routing and colors of the compressor wires. Remove cable ties securing wires to chassis and remove wires from main PC board terminal block.
- 8) Remove and discard the four mounting screws and washers.
- 9) Remove ground wire from old compressor and reinstall onto the new compressor using the ground wire screw previously removed from the old compressor.
- 10) Mount the new compressor using the screws provided in this kit.
- 11) Reroute the wires and attach them to the main PC board terminal block (White wire to J1-10, Black wire to J1-9). Secure wires to chassis at the same points the old compressor wires were secured.
- 12) Cut the new transparent Tubing (RPI Part #RPT373) to size and secure to the barb fitting on the check valve and air compressor using the Hi-Temp Cable Ties (RPI Part #RPT480) provided.
- 13) Run a test cycle or two to ensure autoclave runs properly.

### STATIM® 2000 Model Shown

#### TRANSPARENT SILICONE TUBING (RPI Part #RPT373)

Checking for a leaky Check Valve just got easier. The RPI tubing used in this kit is transparent, making the leak easy to spot! Early detection of a leaky Check Valve can prevent a costly repair. - It's the RPI Advantage!

#### CAUTION: NEVER OIL COMPRESSOR



This MEDO Pump draws air through the Compressor Filter (RPI Part #SCF002) located in the front of the compressor. Air is then pushed through the Biological Filter (RPI Part #SCF032) (if present) and through the Check Valve into the Steam Generator/Boiler.