

Models: Century V116, V120, V136, V148 & V160



Parts to fit AMSCO/STERIS® Century Sterilizers

New Door Motor, Capacitor, and Power Door Assembly Parts!

RPI Part #AMM279
OEM Part #136809-709




DOOR MOTOR

- 115 VAC, 60Hz; 6.5 RPM
- Includes: Mounting hardware
- Also available: Capacitor (3µF, 370 VAC) (RPI Part #AMC280)

Fits: Door Motor Housing

Models: Century V116, V120, V136, V148 & V160



RPI Part #AMK241

OEM Part #129367-563 (Turnbuckle),
129367-562 (Nut, knurled), 150824-374 (Nut, hex)



TURNBUCKLE KIT

- Includes: Turnbuckle, knurled nut, and (1 pc) Hex Nut (#10-24) (RPI Part #RPH112)

Fits: Power Door Assembly

Models: Century V116 & V120



RPI Part #AMC236
OEM Part #93911-961




CABLE (DOOR)

- Size: 31.5" lg.;
- 1/4-28 stud x 5/16" eyelet

Fits: Power Door Assembly

Model: Century V116



RPI Part #AMC235
OEM Part #93918-039




CABLE (SHEAVE)

- Size: 56.25" lg.;
- #10-24 stud x #10 eyelet

Fits: Power Door Assembly

Model: Century V116



RPI Part #AMC238
OEM Part #93918-357




CABLE (DOOR DRIVE)

- Size: 64.25" lg.;
- #10-24 stud x 3/16" catch ball

Fits: Power Door Assembly

Model: Century V120



RPI Part #AMC237
OEM Part #93918-358




CABLE (COUNTER WEIGHT)

- Size: 37.75" lg.;
- 1/4-28 stud x 3/8" eyelet

Fits: Power Door Assembly

Model: Century V120



RPI Part #AMC280
OEM Part #93918-050



CAPACITOR (3µF, 370 VAC)

- Size: 2.06" x 1.22" (oval case) x 2.00" lg.
- Quad 1/4" tabs

Fits: Door Motor Housing

Models: Century V116, V120, V136, V148 & V160



RPI Part #AMP239
OEM Part #93911-952

PULLEY (16")

Fits: Power Door Assembly
Model: Century V116



RPI Part #AMP240
OEM Part #93918-990

PULLEY (20")

Fits: Power Door Assembly
Model: Century V120

* All Cables are stranded stainless steel cable with clear PVC jacket.



replacement parts industries, inc.


CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@riparts.com • WEBSITE www.riparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. Flyer Code 2006AM Rd2 SHOW 07/21

New Parts to fit AMSCO/STERIS® Century Sterilizers ...



Door Motor



Capacitor



Pulleys (16" & 20")



New Parts to fit AMSCO/STERIS® Century Sterilizers

ADDRESS SERVICES REQUESTED

PRESENTED STANDARD
U.S. POSTAGE PAID
VAN NUYS, CA
PERMIT NO. 2399

replacement parts industries, inc.
625 Cochran Street • Simi Valley, California • 93065




Also Available to fit AMSCO/STERIS® Century Sterilizers ...



Turnbuckle Kit



Cables



Door Cable • Sheave Cable • Door Drive Cable • Counterweight Cable