

Air Shields•Hill-Rom•Drager Models: Isolette® C2000, Resuscitaire® & Versalet® 7700 and Ohmeda Models: 190, Care Plus, Giraffe® Omnibed®, IC, & IWS

**INFANT WARMER & INCUBATOR SKIN TEMPERATURE PROBES**




**Ohmeda**  
Fits Models: Care Plus, IC & IWS

**Skin Temperature Probe (Disposable)**  
RPI Part #OMP006 (10 count)  
OEM Part #6600-0208-700

- Includes (10) probes and (10 pieces) Reflective Probe Cover (RPI Part #RPC619)
- \*Note: Also available by the case (RPI Part #OMP007 - 50 count)
- OEM Part #6600-0196-700

**Fits:** Controller



**Skin Temperature Probe (Reusable)**  
RPI Part #OMP008 (1 count)  
New OEM Part #6600-0628-700  
Old OEM Part #0208-0697-700

**Fits:** Controller




**Ohmeda**  
Fits Models: Giraffe® Omnibed®

**Skin Temperature Probe (Disposable)**  
RPI Part #OMP009 (10 count)  
OEM Part #6600-0873-700

- Includes (10) probes and (10 pieces) Reflective Probe Cover (RPI Part #RPC619)
- \*Note: Also available by the case (RPI Part #OMP010 - 50 count)
- OEM Part #6600-0874-700

**Fits:** Probe Panel



**Skin Temperature Probe (Reusable)**  
RPI Part #OMP011 (1 count)  
OEM Part #6600-0875-700

**Fits:** Probe Panel




**Air Shields•Hill-Rom•Drager**  
Fits Models: Isolette® C2000 & Versalet® 7700

**Skin Temperature Probe (Disposable)**  
RPI Part #AIP133 (10 count)  
OEM Part #83 104 15/MU12551

- Includes (10) probes and (10 pieces) Reflective Probe Cover (RPI Part #RPC619)
- \*Note: Also available by the case (RPI Part #AIP134 - 100 count) (OEM Part #Not Available)

**Fits:** Controller



**Skin Temperature Probe (Reusable)**  
RPI Part #AIP135 (1 count)  
OEM Part #83 102 41/MU12533

**Fits:** Controller




**Air Shields•Hill-Rom•Drager**  
Fits Model: Resuscitaire® (RW82-1)

**Skin Temperature Probe (Disposable)**  
RPI Part #AIP130 (10 count)  
OEM Part #81 300 06 & 81 300 08/MU11228

- Includes (10) probes and (10 pieces) Reflective Probe Cover (RPI Part #RPC619)
- \*Note: Also available by the case (RPI Part #AIP131 - 100 count)
- OEM Part #81 300 09/MU11229

**Fits:** Controller



**Skin Temperature Probe (Reusable)**  
RPI Part #AIP132 (1 count)  
OEM Part #81 300 05/MU11225



**Ohmeda Fits Model: 190**

**Skin Temperature Probe (Reusable)**  
RPI Part #OMP012 (1 count)  
OEM #208-0161-300

**Fits:** Controller



**Fits All Infant Warmer and Incubator Skin Temperature Probes**

**Reflective Probe Cover**  
RPI Part #RPC619 (100 count)  
OEM Part #0203-1980-300/68 209 48/MU06944

- Single use only
- Latex free

**Fits:** Skin Temperature Probe




**SWITCHES TO FIT AIR SHIELDS•HILL-ROM•DRAGER**

*Stabilet™ 200, 2000A/B & 2200*

**Power Switch**  
RPI Part #AIS138  
OEM Part #BDB154

- 10A, 125/250 VAC

**Fits:** Controller




*Includes Green Power Switch Lens (RPI Part #AIL139) and (2) Midget Flange Lamps (RPI Part #LMP024)*

*Stabilet™ 200, 1250 & 2000A/B/C*

**Lamp Switch**  
RPI Part #AIS136  
OEM Part #BCZ956/BCZ010

- 6A, 125/250 VAC

**Fits:** Controller




*Includes Blue Lamp Switch Lens (RPI Part #AIL137) and Bi-Pin Lamp (RPI Part #LMP023)*

**Power Switch Lens (Green)**  
RPI Part #AIL139  
OEM Part #BDB155

- Green plastic

**Fits:** Power Switch



**Lamp (Midget Flange)**  
RPI Part #LMP024  
OEM Part #BDB156

- 12VDC, 1.0 W
- Base: Midget flange, T-1-3/4; Incandescent

**Fits:** Power Switch



**Lamp Switch Lens (Blue)**  
RPI Part #AIL137  
(OEM Part # Not Available)

- Blue plastic

**Fits:** Lamp Switch



**Lamp (Bi-Pin)**  
RPI Part #LMP023  
(OEM Part # Not Available)

- 14 VDC, 1.12 W
- Base: Bi-Pin, T-1-3/4; Incandescent

**Fits:** Lamp Switch



\* SEE PAGE 209B FOR PRICING



replacement parts industries, inc.  
CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Ohmeda/GE & Air Shields•Hill-Rom•Drager equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Website 04/06

