
















## DISPOSABLE: DIN/UNSHIELDED LEADWIRE SETS

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END	NOTES
LDM301	3 LEADWIRE SET DISPOSABLE 24" DIN/SNAP	3	~24"	Din/Unshielded; .060" Safety Socket 	Snap  WHITE RED BLACK	<ul style="list-style-type: none"> <li>• 10 sets per package</li> <li><b>Also available:</b> <ul style="list-style-type: none"> <li>• Reusable Leadwire Set (RPI Part #LKM005)</li> <li>• 3 Lead Organizer (RPI Part #LAR004)</li> </ul> </li> </ul>
LDM302	3 LEADWIRE SET DISPOSABLE 24" DIN/SNAP	3	~24"	Din/Unshielded; .060" Safety Socket 	Snap  WHITE GREEN BLACK	<ul style="list-style-type: none"> <li>• 10 sets per package</li> <li><b>Also available:</b> <ul style="list-style-type: none"> <li>• Reusable Leadwire Set (RPI Part #LKM011)</li> <li>• 3 Lead Organizer (RPI Part #LAR004)</li> </ul> </li> </ul>
LDM303	5 LEADWIRE SET DISPOSABLE 24" DIN/SNAP	5	~24"	Din/Unshielded; .060" Safety Socket 	Snap  RED WHITE GREEN BLACK BROWN	<ul style="list-style-type: none"> <li>• 10 sets per package</li> <li><b>Also available:</b> <ul style="list-style-type: none"> <li>• Reusable 24" Leadwire Set (RPI Part #LKM029)</li> <li>• Reusable 40" Leadwire Set (RPI Part #LKM030)</li> <li>• 5 Lead Organizer (RPI Part #LAR005)</li> </ul> </li> </ul>
LDM304	7 LEADWIRE SET DISPOSABLE 24" DIN/SNAP	7	~24"	Din/Unshielded; .060" Safety Socket 	Snap  ORANGE BLUE BROWN BLACK GREEN WHITE RED	<ul style="list-style-type: none"> <li>• 10 sets per package</li> <li><b>Also available:</b> <ul style="list-style-type: none"> <li>• Reusable Leadwire Set (RPI Part #LKM041)</li> </ul> </li> </ul>
LDM308	3 LEADWIRE SET DISPOSABLE 24" DIN/PINCH	3	~24"	Din/Unshielded; .060" Safety Socket 	Pinch  WHITE RED BLACK	<ul style="list-style-type: none"> <li>• 10 sets per package</li> <li><b>Also available:</b> <ul style="list-style-type: none"> <li>• Reusable Leadwire Set (RPI Part #LKM002)</li> <li>• 3 Lead Organizer (RPI Part #LAR004)</li> </ul> </li> </ul>
LDM309	3 LEADWIRE SET DISPOSABLE 24" DIN/PINCH	3	~24"	Din/Unshielded; .060" Safety Socket 	Pinch  WHITE GREEN BLACK	<ul style="list-style-type: none"> <li>• 10 sets per package</li> <li><b>Also available:</b> <ul style="list-style-type: none"> <li>• Reusable Leadwire Set (RPI Part #LKM008)</li> <li>• 3 Lead Organizer (RPI Part #LAR004)</li> </ul> </li> </ul>
LDM310	5 LEADWIRE SET DISPOSABLE 24" DIN/PINCH	5	~24"	Din/Unshielded; .060" Safety Socket 	Pinch  RED WHITE GREEN BLACK BROWN	<ul style="list-style-type: none"> <li>• 10 sets per package</li> <li><b>Also available:</b> <ul style="list-style-type: none"> <li>• Reusable Leadwire Set (RPI Part #LKM026)</li> <li>• 5 Lead Organizer (RPI Part #LAR005)</li> </ul> </li> </ul>

# DISPOSABLE: DIN/UNSHIELDED LEADWIRE SETS *(Continued)*

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END	NOTES
LDM311	7 LEADWIRE SET DISPOSABLE 24" DIN/PINCH	7	~24"	Din/Unshielded; .060" Safety Socket	Pinch 	<ul style="list-style-type: none"> <li>• 10 sets per package</li> </ul> <p><b>Also available:</b></p> <ul style="list-style-type: none"> <li>• Reusable Leadwire Set (RPI Part #LKM038)</li> </ul>