





















DUAL/SHIELDED PATIENT CABLES

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END	SUGGESTED LEADWIRE KITS
KCC003	PATIENT CABLE - 3 LEAD DUAL	3	~10 ft.	6-Pin; 1K Ω Resistors Jumper E to C 	Dual/Shielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM014 & LKM017 Disposable Leadwire Sets RPI Part #'s LDM305 & LDM312 Veterinary Leadwire Set RPI Part #LVM705
KCC009	PATIENT CABLE - 3 LEAD DUAL	3	~10 ft.	6-Pin; Jumper E to B 	Dual/Shielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM014 & LKM017 Disposable Leadwire Sets RPI Part #'s LDM305 & LDM312 Veterinary Leadwire Set RPI Part #LVM705
KCC012	PATIENT CABLE - 3 LEAD DUAL	3	~10 ft.	12-Pin; 1K Ω Resistors 	Dual/Shielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM014 & LKM017 Disposable Leadwire Sets RPI Part #'s LDM305 & LDM312 Veterinary Leadwire Set RPI Part #LVM705
KCC021	PATIENT CABLE - 3 LEAD DUAL	3	~10 ft.	11-Pin; Jumper 5 to 6 with 120M Ω Resistor 	Dual/Shielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM014 & LKM017 Disposable Leadwire Sets RPI Part #'s LDM305 & LDM312 Veterinary Leadwire Set RPI Part #LVM705
KCC034	PATIENT CABLE - 3 LEAD DUAL	3	~10 ft.	10-Socket; 4.7K Ω Resistors 	Dual/Shielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM014 & LKM017 Disposable Leadwire Sets RPI Part #'s LDM305 & LDM312 Veterinary Leadwire Set RPI Part #LVM705 NOTES: • OEM Part #545302
KEC003	PATIENT CABLE - 5 LEAD DUAL	5	~10 ft.	11-Pin 	Dual/Shielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM032 & LKM035 Disposable Leadwire Sets RPI Part #'s LDM307 & LDM314 Veterinary Leadwire Set RPI Part #LVM707
KEC007	PATIENT CABLE - 5 LEAD DUAL	5	~10 ft.	12-Pin; 1K Ω Resistors 	Dual/Shielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM032 & LKM035 Disposable Leadwire Sets RPI Part #'s LDM307 & LDM314 Veterinary Leadwire Set RPI Part #LVM707

DUAL/SHIELDED PATIENT CABLES *(Continued)*

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END	SUGGESTED LEADWIRE KITS
KEC010	PATIENT CABLE - 5 LEAD DUAL	5	~10 ft.	6-Pin; 1K Ω Resistors 	Dual/Shielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM032 & LKM035 Disposable Leadwire Sets RPI Part #'s LDM307 & LDM314 Veterinary Leadwire Set RPI Part #LVM707
KEC016	PATIENT CABLE - 5 LEAD DUAL	5	~10 ft.	17-Pin; 1K Ω Resistors Jumper 1 to 3 with 16K Ω Resistor 	Dual/Shielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM032 & LKM035 Disposable Leadwire Sets RPI Part #'s LDM307 & LDM314 Veterinary Leadwire Set RPI Part #LVM707
KEC022	PATIENT CABLE - 5 LEAD DUAL	5	~10 ft.	10-Socket; 4.7K Ω Resistors 	Dual/Shielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM032 & LKM035 Disposable Leadwire Sets RPI Part #'s LDM307 & LDM314 Veterinary Leadwire Set RPI Part #LVM707 NOTES: • OEM Part #545303