





REUSABLE: BONDED LEADWIRE SETS TO FIT GE MULTI-LINK®

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END	NOTES
LKM052	3 LEADWIRE SET - 30" MULTI-LINK/ SNAP Fits: GE Multi-Link® Cable, Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021) and Patient Cable - 3 Lead Multi-Link® (RPI Part #KCD033)	3 Bonded	~30"	Dual Shielded; Jumper 1 to 5 with 130M Ω Resistor; Compatible with GE Multi-Link® 	Snap 	<ul style="list-style-type: none"> OEM Part #41203-001 Includes: 3 Lead Organizer (RPI Part #LAR004) Also available: Bonded Leadwire Sets (RPI Part #'s LKM053, LKM054 & LKM055) Combined Leadwire Sets (RPI Part #'s LKM056 & LKM057)
LKM053	3 LEADWIRE SET - 30" MULTI-LINK/ PINCH Fits: GE Multi-Link® Cable, Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021) and Patient Cable - 3 Lead Multi-Link® (RPI Part #KCD033)	3 Bonded	~30"	Dual Shielded; Jumper 1 to 5 with 130M Ω Resistor; Compatible with GE Multi-Link® 	Pinch 	<ul style="list-style-type: none"> OEM Part #412682-001 Includes: 3 Lead Organizer (RPI Part #LAR004) Also available: Bonded Leadwire Sets (RPI Part #'s LKM052, LKM054 & LKM055) Combined Leadwire Sets (RPI Part #'s LKM056 & LKM057)
LKM054	5 LEADWIRE SET - 30" MULTI-LINK/ SNAP Fits: GE Multi-Link® Cable, Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021)	5 Bonded	~30"	Dual Shielded; Compatible with GE Multi-Link® 	Snap 	<ul style="list-style-type: none"> OEM Part #41202-001 Includes: 5 Lead Organizer (RPI Part #LAR005) Also available: Bonded Leadwire Sets (RPI Part #'s LKM052, LKM053 & LKM055) Combined Leadwire Sets (RPI Part #'s LKM056 & LKM057)
LKM055	5 LEADWIRE SET - 30" MULTI-LINK/ PINCH Fits: GE Multi-Link® Cable, Patient Cable - 5 Lead Multi-Link® (RPI Part #KED021)	5 Bonded	~30"	Dual Shielded; Compatible with GE Multi-Link® 	Pinch 	<ul style="list-style-type: none"> OEM Part #412681-001 Includes: 5 Lead Organizer (RPI Part #LAR005) Also available: Bonded Leadwire Sets (RPI Part #'s LKM052, LKM053 & LKM054) Combined Leadwire Sets (RPI Part #'s LKM056 & LKM057)