







**REUSABLE: LEADWIRE SETS TO FIT GE MULTI-LINK®**

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END	NOTES
LKM049	10 LEADWIRE SET - CAM 14/BANANA  Fits: CAM 14 Module	10	(6) ~27" (4) ~37"	Dual Shielded; Compatible with GE Multi-Link®  	Banana Plug  	<ul style="list-style-type: none"> <li>OEM Part #420101-002</li> <li>Includes: (2) 5 Lead Organizer (RPI Part #LAR005)</li> <li>Length: ~27": Red (V1), Yellow (V2), Green (V3), Blue (V4), Orange (V5), Violet (V6); and, ~37": Green (RL), White (RA), Black (LA) &amp; Red (LL)</li> </ul> <p><b>Also available:</b></p> <ul style="list-style-type: none"> <li>10 Leadwire Sets (RPI Part #'s LKM050 &amp; LKM051)</li> </ul>
LKM050	10 LEADWIRE SET - CAM 14/PINCH  Fits: CAM 14 Module	10	(6) ~27" (4) ~37"	Dual Shielded; Compatible with GE Multi-Link®  	Pinch  	<ul style="list-style-type: none"> <li>OEM Part #420101-002</li> <li>Includes: (2) 5 Lead Organizer (RPI Part #LAR005)</li> <li>Length: ~27": Red (V1), Yellow (V2), Green (V3), Blue (V4), Orange (V5), Violet (V6); and, ~37": Green (RL), White (RA), Black (LA) &amp; Red (LL)</li> </ul> <p><b>Also available:</b></p> <ul style="list-style-type: none"> <li>10 Leadwire Sets (RPI Part #'s LKM049 &amp; LKM051)</li> </ul>
LKM051	10 LEADWIRE SET - CAM 14/SNAP  Fits: CAM 14 Module	10	(6) ~27" (4) ~37"	Dual Shielded; Compatible with GE Multi-Link®  	Snap  	<ul style="list-style-type: none"> <li>OEM Part #420101-002</li> <li>Includes: (2) 5 Lead Organizer (RPI Part #LAR005)</li> <li>Length: ~27": Red (V1), Yellow (V2), Green (V3), Blue (V4), Orange (V5), Violet (V6); and, ~37": Green (RL), White (RA), Black (LA) &amp; Red (LL)</li> </ul> <p><b>Also available:</b></p> <ul style="list-style-type: none"> <li>10 Leadwire Sets (RPI Part #'s LKM049 &amp; LKM050)</li> </ul>