

Models: 100 Series, 400 Series, 500 Series, 700 Series, 3522, 3523, 3525, 3533, 6300 Series, 6400 Series, 6900 Series, & MC3600

# RPI is Your Best Source for Door Gaskets to fit Getinge/Castle Bulk Sterilizers!


**GETINGE/CASTLE**



**RPI Part #CSG102**  
**OEM Part #61301607691**

**DOOR GASKET (17.50" x 17.50")**

- Material: Blue silicone
  - Custom formulated for greater lubricity!
  - .65" OD cross section
- Fits:** Chamber  
**Models:** 400 Series




**RPI Part #CSG002**  
**OEM Part #61301605017 & 58397**

**DOOR GASKET (16" X 16")**

- Material: Blue silicone
- New material with greater lubricity!
- .43" OD cross section


**Fits:** 16" x 16" Door  
**Models:** 3522, 3525 & 3533



**RPI Part #CSG103**  
**OEM Part #61301607692**

**DOOR GASKET (21" x 21")**

- Material: Blue silicone
  - Custom formulated for greater lubricity!
  - .65" OD cross section
- Fits:** 21.50" x 21.50" Door  
**Models:** 500 Series




**RPI Part #CSG030**  
**OEM Part #513966**

**DOOR GASKET (16" x 16")**

- Material: Blue silicone
- New material with greater lubricity!
- .40" x .40" "U" shaped cross section


**Fits:** 16" x 16" Door  
**Models:** 3522, 3525 & 3533



**RPI Part #CSG104**  
**OEM Part #61301609906**

**DOOR GASKET (26.50" x 36")**

- Material: Blue silicone
  - Custom formulated for greater lubricity!
  - .65" OD cross section
- Fits:** 26.5" x 36.00" Door  
**Models:** 700 Series




**RPI Part #CSG001**  
**OEM Part #61301605016 & 58398**

**DOOR GASKET (20" X 20")**

- Material: Blue silicone
- New material with greater lubricity!
- .43" OD cross section

**Fits:** 20" x 20" Door  
**Models:** 100 Series & 3523




**RPI Part #CSG031**  
**OEM Part #513967**

**DOOR GASKET (20" x 20")**

- Material: Blue silicone
- New material with greater lubricity!
- .40" x .40" "U" shaped cross section

**Fits:** 20" x 20" Door  
**Models:** 100 Series & 3523



**See reverse side for more!**


**Warning: These gaskets are not for use in Ethylene Oxide (E.O. Gas) sterilizers**

**rpi** replacement parts industries, inc.  
CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@ripiparts.com](mailto:order@ripiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Getinge/Castle equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. 1109CS Show 05/17

**SECTION 10**  
**PAGE 310**  
**UPDATED 05/17**

# RPI is Your Best Source for Door Gaskets to fit Getinge/Castle Bulk Sterilizers!



RPI is pleased to introduce 3 new Door Gaskets to fit the following Getinge/Castle Bulk Sterilizers:

- 400 Series
- 500 Series
- 700 Series

Plus, RPI also offers Door Gaskets that are now custom formulated for greater lubricity to fit the following Models:

- 100 Series
- 6300 Series
- 6900 Series
- 3522, 3523, 3525 & 3533
- 6400 Series
- MC3600

RPI is Your Best Source for Door Gaskets to fit Getinge/Castle Bulk Sterilizers!

ADDRESS SERVICES REQUESTED

PREPARED STANDARD  
U.S. POSTAGE PAID  
VAN NUYS, CA  
PERMIT NO. 2399


**replacement parts industries, inc.**   
20338 Corisco Street • Chatsworth, California • 91311



**RPI Part #CSG003**  
**OEM Part #40441**  
**DOOR GASKET (24" x 36")**

- Material: Blue silicone
- New material with greater lubricity!
- .43" OD cross section

**Fits:** 24" x 36" Door  
**Models:** 6300 Series, 6400 Series & 6900 Series




  
CSG003  
CROSS SECTION  
NOT SHOWN ACTUAL SIZE



**RPI Part #CSG032**  
**OEM Part #536651**  
**DOOR GASKET (24" x 36")**

- Material: Blue silicone
- New material with greater lubricity!
- 1" OD cross section

**Fits:** 24" x 36" Door  
**Models:** MC3600



  
CSG032  
CROSS SECTION  
NOT SHOWN ACTUAL SIZE

**Warning: These gaskets are not for use in Ethylene Oxide (E.O. Gas) sterilizers**