


Models: 75L (411), 105, 106, 107, 111, 112, 113, 114, 116, 117, 119, 204, 222, 223, 230, 305, 306, 307, 311, 317, 318, 319, 405, 411, 411 (75L), 413, 414, 415, 417, 419, 604, 622, 623, 630, 641, 646, 647 & Evolution 75

New
Parts
to fit

MIDMARK® • RITTER®

Exam Tables & Chairs




RPI Part #MIA272
OEM Part #002-0564-08

ACTUATOR (FOOT)

- 2.7A @ 115VAC
- Also available: Brake & Coupler Kit (RPI Part #MIK167), Capacitor (RPI Part #MIC182), and Limit Switch (RPI Part #MIS175)

Fits: Foot Section

Model: See RPI website for Serial # breakdown of the model this part fits: 419




RPI Part #MIA218
OEM Part #002-0498-00 & 002-0498-02

ACTUATOR (TILT)

- 2.7A @ 115VAC
- Also available: Brake & Coupler Kit (RPI Part #MIK167) and Capacitor (RPI Part #MIC219)

Fits: Seat Section

Model: See RPI website for Serial # breakdown of the models this part fits: 75L (411), 411 (75L), 411, 419 & Evolution 75




RPI Part #MIA216
OEM Part #002-0496-00 & 002-0564-06

ACTUATOR (BACK)

- 2.7A @ 115VAC
- Also available: Brake & Coupler Kit (RPI Part #MIK167) and Capacitor (RPI Part #MIC183)

Fits: Back Section

Models: See RPI website for Serial # breakdown of the models this part fits: 75L (411), 411 (75L), 411, 419 & Evolution 75



RPI Part #MIK167
OEM Part #016-0237-00 (Brake only)

BRAKE & COUPLER KIT

- This kit is for **Motion Systems version** motor only. (For **Fasco version** motor use **Motor Coupler** (RPI Part #MIC166)).
- Included in: Back & Foot Actuator (RPI Part #MIA164), Base & Tilt Actuator (RPI Part #MIA165), Actuator (Foot) (RPI Part #MIA215), Actuator (Back) (RPI Part #MIA216), Actuator (Tilt) (RPI Part #MIA218), and Actuator (Foot) (RPI Part #MIA272)

Fits: Only fits Back & Foot Actuator and Base & Tilt Actuator manufactured by **Motion Systems**


Models: See RPI website for Serial # breakdown of the models this part fits: 75L (411), 411 (75L), 411, 419 & Evolution 75



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@ripiparts.com • WEBSITE www.ripiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Midmark®/Ritter® equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Flyer Code 1807MI Show 03/19




RPI Part #MIC298
OEM Part #016-0356-00

GAS CYLINDER

- See **Model** below for quantities required for each model serial number

Fits: Power Base Assembly

Model: 414-001: All Serial #'s, Serial #BM1000 to BM1050 requires (1) Cylinder, Serial #BM1051 to present requires (2) Cylinders




RPI Part #MIC299
OEM Part #016-1625-00

LOCKING GAS CYLINDER

Fits: Back Assembly

Models: See RPI website for Serial # breakdown of the models this part fits: 222 & 622




RPI Part #MIC300
OEM Part #016-1017-00

LOCKING GAS CYLINDER

Fits: Back Assembly

Models: See RPI website for Serial # breakdown of the models this part fits: 204 & 604



RPI Part #MIC310
OEM Part #002-1065-00

GAS CYLINDER SET

- Includes: (2) Gas Cylinders and (4) Safety Clips

Fits: Scissor Assembly

Models: See RPI website for Serial # breakdown of the models this part fits: 230, 630, 641, 647 & Evolution 75

RPI Part #MIC176


OEM Part #002-0213-00, 016-0236-00 & 016-0236-01 ("Double Gas Cylinder" design)

GAS CYLINDER

- **Caution** - Use the Support Bar (RPI Part #MIT177) to support the Table Top and the Pry Bar (RPI Part #MIT178) to safely compress the Gas Cylinder during the removal and installation process.
- This Gas Cylinder replaces the "Double Gas Cylinder" design used in older versions of these tables

Fits: Middle Member and Chassis Assembly

Models: See RPI website for Serial # breakdown of the models this part fits: 411, 413 & 419




RPI Part #RXB035
OEM Part #016-0131-21

FLANGE BEARING

- 12 per package
- .375" ID x .501" OD x .69" dia. x .062" Flange
- Material: Oil embedded bronze

Fits: Foot, Seat, Back and Column Assemblies

Models: See RPI website for Serial # breakdown of the models this part fits: 230, 630 & Evolution 75



RPI Part #RXF036
OEM Part #015-0346-40

FUSE (6.3A, 250V)
5mm x 20mm TIME DELAY

- 5 per package
- Type: Ceramic
- Size: 5mm dia. x 20mm lg.

Fits: Power Inlet Module



Models: See RPI website for Serial # breakdown of the models this part fits: 230, 630, 641, 647 & Evolution 75

RPI Part #MIS255
OEM Part #015-1055-00 & 62977400

LIMIT SWITCH

- Momentary (on - off)
- 0.1A @ 125VAC, SPDT
- 3/16" tabs
- Panel cutout: 1.125" x .550"

Models: See RPI website for Serial # breakdown of the models this part fits: 75L (411), 222, 223, 230, 411, 411 (75L), 622, 623, 630, 641, 646, 647 & Evolution 75




RPI Part #MIS175
OEM Part #015-0476-00 & 015-0430-00

LIMIT SWITCH

- 0.1 Amp @ 125VAC
- 3/16" Tabs

Fits: Column Assembly and Base Actuator

Models: See RPI website for Serial # breakdown of the models this part fits: 75L (411), 411, 411 (75L), 413, 415, 419 & Evolution 75




RPI Part #MIS016
OEM Part #015-0421-00 & 015-0381-00

RETURN LIMIT SWITCH

- 20A@125/250VAC; DPDT
- 1/4" Tabs

Fits: Column Assembly and Base Actuator

Models: See RPI website for Serial # breakdown of the models this part fits: 111, 112, 113, 114, 117 & 317



RPI Part #MIS015
OEM Part #002-0101-00 & 015-0072-00

MICROSWITCH

- 15A @ 125/250VAC; 1/2A @ 125 VDC; 1/4A @ 250 VDC
- Includes: Switch and (2 pcs) 3/16" QD terminal (RPI Part #RPT257)
- 3/16" Tabs

Fits: Foot Control and Front Panel

Category: Exam Tables & Chairs

Models: 1 See RPI website for Serial # breakdown of the models this part fits: 105, 106, 107, 117, 222, 305, 306, 307, 405, 413, 415 & 622


RPI Part #TUS009
OEM Part #015-0650-03

POWER SWITCH

- Rocker switch (On-Off), Legend: I/O; Green lens
- 20A @ 125/250VAC; DPST
- Panel Cutout: 1.185" x .874"; 1/4" Tabs

Fits: Main Housing

Models: See RPI website for Serial # breakdown of the models this part fits: 222, 223, 622 & 623



RPI Part #MIC284
OEM Part #002-0911-04




HAND CONTROL (BASIC)

- Includes all parts as shown
- See **Note** below

Fits: Front or rear Connection Port

Models: See RPI website for Serial # breakdown of the models this part fits: 230 & Evolution 75



Just Need The Panel?

RPI Part #MIP285
OEM Part #015-1866-00

MEMBRANE PANEL

- Size: 2.47" x 5.05"; Adhesive backing



RPI Part #MIC288
OEM Part #002-0911-01




HAND CONTROL (NON-PROGRAMMABLE)

- Includes all parts as shown
- See **Note** below

Fits: Front or Rear Connection Port

Models: See RPI website for Serial # breakdown of the models this part fits: 630, 641 & Evolution 75



Just Need The Panel?

RPI Part #MIP289
OEM Part #015-1677-00

MEMBRANE PANEL

- Size: 2.47" x 5.05"; Adhesive backing



RPI Part #MIC286
OEM Part #002-0911-02




HAND CONTROL (PROGRAMMABLE)

- Includes all parts as shown
- See **Note** below

Fits: Front or Rear Connection Port

Models: See RPI website for Serial # breakdown of the models this part fits: 630, 641 & Evolution 75



Just Need The Panel?

RPI Part #MIP287
OEM Part #015-1678-00

MEMBRANE PANEL

- Size: 2.47" x 5.05"; Adhesive backing



RPI Part #MIC290
OEM Part #002-0911-05




HAND CONTROL (PROGRAMMABLE WITH HEATER)

- Includes all parts as shown
- See **Note** below

Fits: Front or Rear Connection Port

Model: See RPI website for Serial # breakdown of the model this part fits: 630



Just Need The Panel?

RPI Part #MIP291
OEM Part #015-1765-00

MEMBRANE PANEL

- Size: 2.47" x 5.05"; Adhesive backing



RPI Part #MIC282
OEM Part #002-0911-00/
002-0911-06




HAND CONTROL

- Includes all parts as shown
- See **Note** below

Fits: Front or Rear Connection Port

Models: See RPI website for Serial # breakdown of the models this part fits: 223 & 623



Just Need The Panel?

RPI Part #MIP283
OEM Part #015-1684-00

MEMBRANE PANEL

- Size: 2.47" x 5.05"; Adhesive backing



RPI Part #MIC307
OEM Part #002-0624-00




HAND CONTROL

- Includes all parts as shown
- See **Note** below

Fits: Front or Rear Connection Port

Model: See RPI website for Serial # breakdown of the model this part fits: 405



Just Need The Panel?

RPI Part #MIP308
OEM Part #015-1243-00

MEMBRANE PANEL

- Size: 2.47" x 5.05"; Adhesive backing




Don't Forget the The Hand Control Holster Kit (RPI Part #MIK252)!

RPI Part #MIK252
OEM Part #002-1474-00

HAND CONTROL HOLSTER KIT

Models: See RPI website for Serial # breakdown of the models this part fits: 223, 230, 405, 411, 419, 623, 630, 641 & Evolution 75

Allow your customer to securely store, and easily access any of the hand controls shown above. Simply install an RPI Hand Control Holster Kit on the side rail of the exam chair, or on a wall near the chair. This Kit includes everything you need to quickly mount the holster and/or the side rail bracket for your customer's convenience.



RPI Part #MIC304
OEM Part #029-0681-00
 (Velcro Version)/002-0409-00
 (Clip On Version)




HAND CONTROL

- Includes options for both clip on and Velcro mounting.
- Includes all parts as shown

Fits: Hand Control

Model: See RPI website for Serial # breakdown of the model this part fits: 405



Just Need The Panel?

RPI Part #MIP305
OEM Part #015-0528-00



HAND CONTROL PANEL

- Size: 2.16" W x 8.28" lg.
 - Metal backing plate
- Fits:** Hand Control




RPI Part #MIC277
OEM Part #002-0499-00



HAND CONTROL (NON-PROGRAMMABLE)

- Includes all parts as shown

Fits: Chassis Connector
Models: See RPI website for Serial # breakdown of the models this part fits: 411, 411 (75L), 75L (411) & Evolution 75




Just Need The Panel?

RPI Part #MIA278
OEM Part #015-0995-00



CONTROL PANEL ASSEMBLY

- Size: 2.12" x 5.11"
- Fits:** Hand Control Assembly




RPI Part #MIC280
OEM Part #002-0578-00



HAND CONTROL (PROGRAMMABLE)

- Includes all parts as shown

Fits: Chassis Connector
Models: See RPI website for Serial # breakdown of the models this part fits: 411 (75L), 75L (411) & Evolution 75



Just Need The Panel?

RPI Part #MIP281
OEM Part #015-1127-00




HAND CONTROL PANEL

- Size: 2.16" x 5.56"
 - Metal backing plate
- Fits:** Hand Control



RPI Part #MIB250
OEM Part #015-1066-00




PC BOARD

- Includes: Mounting hardware
- Included in: Hand Control (Programmable) (RPI Part #'s MIC248 and MIC280)
- Also available: Coiled Cord (RPI Part #MIC251)

Fits: Hand Control Case

Models: See RPI website for Serial # breakdown of the models this part fits: 75L (411), 411 (75L) & Evolution 75

RPI Part #MIB306
OEM Part #002-0347-06




HAND CONTROL PC BOARD

- Included in: Hand Control (RPI Part #MIC304)
- Also available: Coil Cord (RPI Part #MIC174) and Membrane Panel (RPI Part #MIP305)

Fits: Hand Control

Model: See RPI website for Serial # breakdown of the model this part fits: 405

RPI Part #MIB309
OEM Part #015-1173-00



CONTROL INLET PCB

- Also available: Coiled Cord (RPI Part #MIC251) and Coiled Cord (Foot Control) (RPI Part #MIC293)

Fits: Chassis

Models: See RPI website for Serial # breakdown of the models this part fits: 75L (411), 411, 411 (75L) & Evolution 75


RPI Part #MIC174
OEM Part #015-0505-03

COILED CORD

- Length: ~ 12 ft. extended
- Included in: Hand Control (RPI Part #'s MIC168 and MIC304)

Fits: Hand Control

Models: See RPI website for Serial # breakdown of the models this part fits: 405, 411 & 417




RPI Part #MIC251
OEM Part #015-0839-00


COILED CORD

- Length: ~ 25 ft. extended
- Included in: Hand Control (Programmable) (RPI Part #'s MIC248 and MIC280), Hand Control (Non-Programmable) (RPI Part #'s MIC262 and MIC277), and Hand Control (RPI Part #MIC307)

Fits: Hand Control and Foot Control

Models: See RPI website for Serial # breakdown of the models this part fits: 75L (411), 319, 405, 411, 411 (75L), 419, & Evolution 75






RPI Part #MIC256
OEM Part #015-1675-00

COILED CORD

- Length: ~ 14 ft. extended
- Included in: Hand Control (Basic) (RPI Part #MIC284), Hand Control (Programmable) (RPI Part #MIC286), Hand Control (Non-Programmable) (RPI Part #MIC288), and Hand Control (Programmable With Heater) (RPI Part #MIC290)
- Also available: Connection Port Harness (Short) (RPI Part #MIH257) and Connection Port Harness (Long) (RPI Part #MIH258)

Fits: Hand Control and Foot Control

Models: See RPI website for Serial # breakdown of the models this part fits: 230, 630, 641, 646, 647 & Evolution 75




RPI Part #MIC292
OEM Part #015-1675-04

COILED CORD

- Length: ~ 25 ft. extended
- Included in: Hand Control (RPI Part #MIC282)

Fits: Hand Control and Foot Control

Models: See RPI website for Serial # breakdown of the models this part fits: 223 & 623




RPI Part #MIC293
OEM Part #015-1785-00

COILED CORD (FOOT CONTROL)

- Length: ~ 25 ft. extended
- One end of coiled cord is 12" long to easily fit through the foot switch case

Fits: Foot Control

Models: See RPI website for Serial # breakdown of the models this part fits: 223 & 623




RPI Part #MIC273
OEM Part #002-0048-00

FOOT CONTROL CORD

- Sold by the foot (reel not included)
- 8 conductor #18 AWG 19/30 strands
- Colors: Black, Red, Orange, Green, Blue, White, White/Black and Red/Black
- Grey TPR jacket, UL 20234, 176°F (80°C), 600V; .54" OD cord

Fits: Foot Control and Base Assembly

Models: See RPI website for Serial # breakdown of the models this part fits: 105, 106, 107, 111, 112, 113, 114, 116, 117, 119, 305, 306, 307, 311, 317, 318, 319 & Evolution 75




RPI Part #MIH257
OEM Part #002-1110-00

CONNECTION PORT HARNESS (SHORT)

- Material: 18 conductor cable; black braided sleeve; (2) 18 Pin connectors
- Length: 13"
- Also available: Connection Port Harness (Long) (RPI Part #MIH258)

Fits: Hand Control and Foot Control

Models: See RPI website for Serial # breakdown of the models this part fits: 230, 630, 641 & Evolution 75




RPI Part #MIH258
OEM Part #002-1110-01

CONNECTION PORT HARNESS (LONG)

- Material: 18 conductor cable; black braided sleeve; (2) 18 Pin connectors
- Length: 26"
- Also available: Connection Port Harness (Short) (RPI Part #MIH257)

Fits: Hand Control and Foot Control

Models: 230, 630, 641, 647 & Evolution 75




RPI Part #MIC171
OEM Part #053-0412-00

TOP END CAP (CLIP ON)

- Includes: Mounting hardware
- Included in: Hand Control (RPI Part #'s MIC168 and MIC304), Hand Control (Non-Programmable) (RPI Part #MIC277), and Hand Control (Programmable) (RPI Part #MIC280)
- Also available: Top End Cap (Velcro) (RPI Part #MIC172) and Bottom End Cap (RPI Part #MIC173)

Fits: Hand Control

Models: See RPI website for Serial # breakdown of the models this part fits: 75L (411), 405, 411, 411 (75L) & Evolution 75




RPI Part #MIC172
OEM Part #053-0257-00

TOP END CAP (VELCRO)

- Includes: Mounting hardware
- Included in: Hand Control (RPI Part #'s MIC168 and MIC304)
- Also available: Top End Cap (Clip On) (RPI Part #MIC171) and Bottom End Cap (RPI Part #MIC173)

Fits: Hand Control

Models: See RPI website for Serial # breakdown of the models this part fits: 75L (411), 319, 405, 411, 411 (75L), 417 & Evolution 75




RPI Part #MIC173
OEM Part #053-0256-00

BOTTOM END CAP

- Includes: Mounting hardware
- Included in: Hand Control (RPI Part #'s MIC168 and MIC304), Hand Control (Non-Programmable) (RPI Part #MIC277), and Hand Control (Programmable) (RPI Part #MIC280)
- Also available: Top End Cap (Clip On) (RPI Part #MIC171) and Top End Cap (Velcro) (RPI Part #MIC172)

Fits: Hand Control

Models: See RPI website for Serial # breakdown of the models this part fits: 75L (411), 319, 405, 411, 411 (75L) 417 & Evolution 75




RPI Part #MIA302
OEM Part #015-1758-00

**SWITCH ASSEMBLY
 (UPHOLSTERY HEATER)**

- Momentary 3A @ 125/250VAC; SPDT with 9" lg. cable and (4) 1/4" male fully insulated QC terminals
- Illuminated green lens; 12 VDC lamp
- Panel cut-out: 7/8" dia.
- Includes: Mounting hardware

Fits: Chassis

Models: See RPI website for Serial # breakdown of the models
 this part fits: 604, 622 & 623



RPI Part #MIA303
OEM Part #015-1759-00

**SWITCH ASSEMBLY
 (DRAWER HEATER)**

- 3A @ 125/250 VAC; DPST with 20-1/2" lg. cable with (2) 1/4" male fully insulated QC terminals, and (2) 1/4" female fully insulated flag QC terminals
- Illuminated white lens; 12 VDC lamp
- Panel cut-out: 7/8" dia.
- Includes: Mounting hardware

Fits: Chassis


Models: See RPI website for Serial # breakdown of the models
 this part fits: 204, 222, 223, 604, 622 & 623

New Parts to fit Midmark®•Ritter® Exam Tables and Chairs!

ADDRESS SERVICES REQUESTED

PERMIT NO. 2399
 VAN NUYS, CA
 U.S. POSTAGE PAID
 PRESORTED STANDARD

r.p.i. replacement parts industries, inc.
 625 Cochran Street • Simi Valley, California • 93065




RPI Part #MIC295
OEM Part #053-0848-00

STEP COVER (GREY)

- Material: Grey slip-resistant tread with PSA backing
- Size: 12" x 17"

Fits: Step Assembly

Models: See RPI website for Serial # breakdown of the models
 this part fits: 100, 104 & 304



RPI Part #MIC296
OEM Part #053-0048-00

STEP COVER (BLACK)

- Material: Black slip-resistant tread with PSA backing
- Size: 12" x 16"

Fits: Step Assembly

Category: Exam Tables & Chairs
OEM: MIDMARK® - RITTER

Model: See RPI website for Serial # breakdown of the model this part fits: 100