



**New Parts to fit**

**Midmark • Ritter**

**Evolution 75, 75L, 200, 300, 400 & 600 Series Tables!**



**RPI Part #MIC248  
OEM Part #002-0604-00**

**HAND CONTROL (PROGRAMMABLE)**


- Includes: Coiled Cord (RPI Part #MIC251), PC Board (RPI Part #MIB250) and Membrane Panel (Programmable) (RPI Part #MIP249)
- Also available: Hand Control Holster Kit (RPI Part #MIK252)

**Models and Serial Numbers this part fits:**

- **411 (75L)-011 thru -018:** Serial #HZ1000 to Present, JY1000 to Present, LS1000 to Present.
- **Evolution 75-014:** Serial # HZ: All Serial #'s.
- **Evolution 75-017:** Serial # JY: All Serial #'s.

**All Parts also sold separately! - see Page 2 for details**

- Membrane Panel (Programmable) (RPI Part #MIP249)
- PC Board (RPI Part #MIB250)
- Coiled Cord (RPI Part #MIC251)



**RPI Part #MIC262  
OEM Part #002-0603-00**

**HAND CONTROL (NON-PROGRAMMABLE)**


- Includes: Coiled Cord (RPI Part #MIC251), PC Board (RPI Part #MIB264) and Membrane Panel (Non-Programmable) (RPI Part #MIP263)
- Also available: Hand Control Holster Kit (RPI Part #MIK252)

**Models and Serial Numbers this part fits:**

- **411 (75L)-011 thru -018:** Serial #HY1000 to Present, JX1000 to Present.
- **Evolution 75-013:** Serial # HY: All Serial #'s.
- **Evolution 75-016:** Serial # JX: All Serial #'s.

**All Parts also sold separately! - see Page 2 for details**

- Membrane Panel (Non-Programmable) (RPI Part #MIP263)
- PC Board (RPI Part #MIB264)
- Coiled Cord (RPI Part #MIC251)



**replacement parts industries, inc.**

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Midmark®•Ritter® equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Flyer Code 1605MI SHOW 07/17

# Tables & Chairs (Exam/Surgical)

Models: Evolution 75, 75L, 200, 300, 400, and 600 Series

## MEMBRANE PANEL (PROGRAMMABLE), PC BOARD, & COILED CORD

Serial Numbers these parts fit:

- Models 411 (75L)-011 thru -018: Serial #HZ1000 to Present, JY1000 to Present, LS1000 to Present.
- Model Evolution 75-014: Serial #HZ: All Serial #'s.
- Model Evolution 75-017: Serial #JY: All Serial #'s.



RPI Part #MIP249  
OEM Part #015-1116-00


### MEMBRANE PANEL (PROGRAMMABLE)

- Size: 2.47" x 5.05"
- Adhesive backing
- Included in: Hand Control (Programmable) (RPI Part #MIP248)
- Also available: Coiled Cord (RPI Part #MIC251) and PC Board (RPI Part #MIB250)

Fits: Hand Control

Models: 411, 75L & Evolution 75

(See **Serial Numbers these parts fit above**)



RPI Part #MIB250  
OEM Part #015-1066-00


### PC BOARD

- Includes: Mounting hardware
- Included in: Hand Control (Programmable) (RPI Part #MIC248)
- Also available: Membrane Panel (Programmable) (RPI Part #MIP249) and Coiled Cord (RPI Part #MIC251)

Fits: Hand Control Case

Models: 411, 75L & Evolution 75

(See **Serial Numbers these parts fit above**)



RPI Part #MIC251  
OEM Part #015-0839-00


### COILED CORD

- Material: 8 conductor cable, black jacket, RJ45 connectors with protective boots
- Length: 25 ft. extended
- Included in: Hand Control (Programmable) (RPI Part #MIC248) and Hand Control (Non-Programmable) (RPI Part #MIC262)

Fits: Hand Control and Foot Control

Models and Serial Numbers this part fits:

- Models Evolution 75 -009 thru -014, 016 & 017: All Serial #'s.
- Models 75L (411) -004 thru -018: All Serial #'s.
- Models 319-001 and -002: All Serial #'s.
- Models 405-001 thru -018: Serial #LE1000 to present, Serial #LF1000 to present, Serial #LG1000 to present, Serial #V2200 to present.
- Models 411-009 thru -018: All Serial #'s.
- Models 411 (75L) -009 thru -018: All Serial #'s.
- Models 419-003 thru -006: All Serial #'s.



## MEMBRANE PANEL (NON-PROGRAMMABLE), PC BOARD, & COILED CORD

Serial Numbers these parts fit:

- Models 411 (75L)-011 thru -018: Serial #HY1000 to Present, JX1000 to Present.
- Model Evolution 75-013: Serial #HY: All Serial #'s.
- Model Evolution 75-016: Serial #JX: All Serial #'s.



RPI Part #MIP263  
OEM Part #015-1134-00,  
015-1288-00, 015-1438-00


### MEMBRANE PANEL (NON-PROGRAMMABLE)

- Size: 2.47" x 5.05"
- Adhesive backing
- Included in: Hand Control (Non-Programmable) (RPI Part #MIC262)
- Also available: Coiled Cord (RPI Part #MIC251) and PC Board (RPI Part #MIB264)

Fits: Hand Control

Models: 411, 75L & Evolution 75

(See **Serial Numbers these parts fit above**)



RPI Part #MIB264  
OEM Part #015-1291-00


### PC BOARD

- Includes: Mounting hardware
- Included in: Hand Control (Non-Programmable) (RPI Part #MIC262)
- Also available: Membrane Panel (Non-Programmable) (RPI Part #MIP263) and Coiled Cord (RPI Part #MIC251)

Fits: Hand Control Case

Models: 411, 75L & Evolution 75

(See **Serial Numbers these parts fit above**)



RPI Part #MIC251  
OEM Part #015-0839-00


### COILED CORD

- Material: 8 conductor cable, black jacket, RJ45 connectors with protective boots
- Length: 25 ft. extended
- Included in: Hand Control (Programmable) (RPI Part #MIC248) and Hand Control (Non-Programmable) (RPI Part #MIC262)

Fits: Hand Control and Foot Control

Models and Serial Numbers this part fits:


- Models Evolution 75 -009 thru -014, 016 & 017: All Serial #'s.
- Models 75L (411) -004 thru -018: All Serial #'s.
- Models 319-001 and -002: All Serial #'s.
- Models 405-001 thru -018: Serial #LE1000 to present, Serial #LF1000 to present, Serial #LG1000 to present, Serial #V2200 to present.
- Models 411-009 thru -018: All Serial #'s.
- Models 411 (75L) -009 thru -018: All Serial #'s.
- Models 419-003 thru -006: All Serial #'s.



## HAND CONTROL HOSTER

### Serial Numbers these parts fit:

- Models 230-001 thru -004: All serial #'s.
- Models 411-009 thru -018: All serial #'s.
- Models 419-003 thru -006: All serial #'s.
- Models 630-001 thru -009: All serial #'s.
- Models 641-002 thru -005: All serial #'s.
- Models Evolution 75-009 thru -014: All serial #'s.
- Models Evolution 75-016 thru -017: All serial #'s.
- Models Evolution 75-019 thru -021: All serial #'s.




RPI Part #MIK252  
OEM Part #002-1474-00

### HAND CONTROL HOLSTER KIT

- Upgrade kit for all tables using teardrop style hand control
- Accommodates side rail or wall mount attachment
- Includes: Hand Control Holster (RPI Part #MIH253) and Side Rail Hand Control Bracket (RPI Part #MIB254)

**Fits:** Side Rail and Wall Mount

**Models:** 230, 411, 419, 630, 641 & Evolution 75  
(See *Serial Numbers these parts fit above*)




RPI Part #MIH253  
OEM Part #002-1472-00

### HAND CONTROL HOLSTER

- Material: Black ABS
- Includes: Mounting hardware
- Included in: Hand Control Holster Kit (RPI Part #MIK252)
- Also available: Side Rail Hand Control Bracket (RPI Part #MIB254)

**Fits:** Side Rail and Wall Mount

**Models:** 230, 411, 419, 630, 641 & Evolution 75  
(See *Serial Numbers these parts fit above*)



RPI Part #MIB254  
OEM Part #053-1388-00 & 053-2268-00-312

### SIDE RAIL HAND CONTROL BRACKET

- Material: Black ABS
- Included in: Hand Control Holster Kit (RPI Part #MIK252)
- Also available: Hand Control Holster (RPI Part #MIH253)

**Fits:** Side Rail

**Models:** 230, 411, 419, 630, 641 & Evolution 75  
(See *Serial Numbers these parts fit above*)

## CONNECTION PORT HARNESES

### Serial Numbers these parts fit:

- Models 230-001 thru -004: All serial #'s.
- Models 630-001 thru -009: All serial #'s.
- Models Evolution 75-019 thru -021: All serial #'s.

RPI Part #MIH257  
OEM Part #002-1110-00




### CONNECTION PORT HARNESS (SHORT)

- Material: 18 conductor cable; black braided sleeve; (2) 18 Pin connectors
- Length: 13"

**Fits:** Hand Control and Foot Control  
**Models:** 230, 630 & Evolution 75

(See *Serial Numbers these parts fit above*)



RPI Part #MIH258  
OEM Part #002-1110-01




### CONNECTION PORT HARNESS (LONG)

- Material: 18 conductor cable; black braided sleeve; (2) 18 Pin connectors
- Length: 26"

**Fits:** Hand Control and Foot Control  
**Models:** 230, 630 & Evolution 75

(See *Serial Numbers these parts fit above*)



## LIMIT SWITCHES

Serial Numbers these parts fit: See Models and Serial Numbers below

RPI Part #MIS255  
OEM Part #015-1055-00 & 62977400




### LIMIT SWITCH

- Momentary (on - off)
- 0.1A @ 125VAC, SPDT; 3/16" tabs
- Panel cutout: 1.125" x .550"

**Models and Serial Numbers this part fits:**

- Models 222-001 thru -017: All serial #'s.
- Models 223-001 thru -016: All serial #'s.
- Models 411 (75L)-009 thru -018: All serial #'s.
- Models 622-001 thru -008: All serial #'s.
- Models 623-001 thru -010: All serial #'s.
- Models Evolution 75-009 thru -010: All serial #'s.



RPI Part #MIS259  
OEM Part #015-1024-00




### LIMIT SWITCH

- 15A @ 277 VAC; SPDT; 3/16" Tabs

**Fits:** Pan Switch Bracket

**Models and Serial Numbers this part fits:**

- Models 411 (75L)-009 thru -018: All serial #'s.
- Models Evolution 75-009 thru -014: All serial #'s.
- Models Evolution 75-016 thru -017: All serial #'s.




New Parts to fit

Midmark® • Ritter®

Evolution 75, 75L, 200, 300, 400 & 600 Series Tables!

RPI is pleased to offer:

- Hand Controls
- PC Boards
- Membrane Panels
- Hand Control Holster and Kit
- Side Rail Hand Control Bracket
- Limit Switches
- Connection Port Harnesses




ADDRESS SERVICES REQUESTED

PERMITTED NO. 2399  
VAN NUYS, CA  
U.S. POSTAGE PAID  
PRESORTED STANDARD


replacement parts industries, inc.   
20338 Corisco Street • Chatsworth, California • 91311

Don't forget, RPI offers Actuators to fit Midmark® 411 (75L)\* Exam Tables!


RPI Part #MIA215  
OEM Part #002-0495-00  
ACTUATOR (FOOT)  
• 2.7A @ 115VAC  
Models: 411




RPI Part #MIA217  
OEM Part #002-0497-00  
ACTUATOR (BASE)  
• 2.7A @ 115VAC  
Models: 411



RPI Part #MIA216  
OEM Part #002-0496-00  
ACTUATOR (BACK)  
• 2.7A @ 115VAC  
Models: 411



RPI Part #MIA218  
OEM Part #002-0498-00  
ACTUATOR (TILT)  
• 2.7A @ 115VAC  
Models: 411



Capacitors also available to fit these Models and more! See our website [www.rpiparts.com](http://www.rpiparts.com) for more information.

\* Fits Models 411 (75L)-009 and -010: Serial # prefixes FA & FB, and 411 (75L)-011 thru -018: Serial # prefixes GT, GV, JX, JY & LS