

New Parts to fit

Midmark® • Ritter® UltraClave® M9 & M11 Sterilizers!




RPI Part #MIH327
OEM Part #015-1639-00

HEATER WIRE HARNESS

- Size: #14 AWG, 14.00" lg. wire, 1/4" QC terminals
- Also available: Jumper Wire (RPI Part #MIJ328) and Main Harness (RPI Part #MIH238)

Fits: PC Board to Heater (Newer Style Units Only)

Models: M9-020 thru -022, M9D-020 thru -022, M11-020 thru -022 & M11D-020 thru -022




RPI Part #MIJ328
OEM Part #015-1640-00

JUMPER WIRE

- Size: #14 AWG, 2.85" lg. wire, 1/4" QC flag terminals
- Included in: Main Harness (RPI Part #MIH238)
- Also available: Heater Wire Harness (RPI Part #MIH327)

Fits: Hi Limit Thermostat to Thermostat (Newer Style Units Only)

Models: M9-020 thru -022, M9D-020 thru -022, M11-020 thru -022 & M11D-020 thru -022



Drain Hose Kits to fit M9 & M11 Sterilizers with two different Hose Clip options!
Check Serial Numbers for which Kit you need!



Fits
Newer Models

RPI Part #MIK334
OEM Part # (No OEM Part # Available)

DRAIN HOSE KIT

- 13" lg. drain hose can be cut to fit as needed.
- Includes: (2 pcs) Hose Clip (RPI Part #MIC333), (1 pc) Kwik Clamp (RPI Part #RPC286), (13" Long) Tubing (3/8" ID x 5/8" OD Clear) (RPI Part #RPH285)

• Also available: Drain Hose Kit (RPI Part #MIK061) for older style units

Fits: Reservoir to Front Panel with Hole Diameter Size: .252" thru .267" for the "Revised Hose Clip Replacement"

Models: M9/M9D (-020 thru -022 Serial #'s V1464210-Present) & M11/M11D (-020 thru -022 Serial #'s V1479341-Present)



Fits
Older Models

RPI Part #MIK061
OEM Part # (No OEM Part # Available)


DRAIN HOSE KIT

- 13" lg. drain hose can be cut to fit as needed.
- Includes: (2 pcs) Hose Clip (RPI Part #RPC283), (2 pcs) Pop Rivet (RPI Part #RPR284), (1 pc) Kwik Clamp (RPI Part #RPC286), (13" Long) Tubing (3/8" ID x 5/8" OD Clear) (RPI Part #RPH285)

• Also available: Drain Hose Kit (RPI Part #MIK334) for newer style units

Fits: Reservoir to Front Panel (Older Style Units Only)

Models: M9/M9D (-001 thru -019), M11/M11D (-001 thru -019), M9/M9D (-020 thru -022 Serial #'s V2200-V1464209) & M11/M11D (-020 thru -022 Serial #'s V2200-V1479340)



The Grommet is now included in the Steam Trap Kit!

RPI Part #RXG099
OEM Part #053-0394-02

GROMMET

- 6 per package
- Size: .250" ID x .500" OD x .250" Thk.; Groove Size: .375" dia. x .125" wide
- Material: Black SBR Rubber, 60 Duro

Fits: Water Reservoir to 1/4" Copper Tubing, Groove Size: .375" dia. x .125" wide (Older Style Units Only)

Models: M9-001 thru -019, M9D-001 thru -019, M11-001 thru -019 & M11D-001 thru -019



RPI Part #MIK078
OEM Part #002-0654-00

STEAM TRAP KIT (New Style/Upgrade Kit)

- Includes: Steam Trap Valve (RPI Part #MIV081), (1 pc) Grommet (RPI Part #RXG099), (3 pcs) Compression Nut (RPI Part #RCN059), and copper tubing

Fits: Chamber to Condensating Tank (Older Style Units Only)

Models: M9-001 thru -019, M9D-001 thru -019, M11-001 thru -019 & M11D-001 thru -019



replacement parts industries, inc. *A Healthcare Components Group Company*
CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. Flyer Code 1703MI Web 11/21

Make Your Midmark®·Ritter® UltraClave® Look as Good as New with RPI's Decal Labels!




RPI Part #MIK329
OEM Part #002-10231-00, 061-0311-00

DECAL KIT

- Ideal for bringing an old sterilizer back to life.
- Includes: (2 pcs) Decal (Hot Surface) (RPI Part #MID330), Decal (Distilled Water) (RPI Part #MID331), Decal (Water Level) (RPI Part #MID332), and Decal (Caution) (RPI Part #MID335)

Fits: Front and Rear Panels
Models: M9, M9D, M11, M11D



RPI Part #MID330
OEM Part #061-0726-00

DECAL (HOT SURFACE)

- 2 per package
- Size: 1.50" x 1.56"
- Included in: Decal Kit (RPI Part #MIK329)

Fits: Front Panel and Inside of Door Panel
Models: M9, M9D, M11 & M11D




RPI Part #MID331
OEM Part #061-0827-00

DECAL (DISTILLED WATER)

- Size: 2.00" x 2.75"
- Included in: Decal Kit (RPI Part #MIK329)

Fits: Front Panel
Models: M9, M9D, M11 & M11D




RPI Part #MID332
OEM Part #061-0311-00

DECAL (WATER LEVEL)

- Size: 1" x 4.25"
- Included in: Decal Kit (RPI Part #MIK329)

Fits: Front Panel Behind the Drain Hose
Models: M9, M9D, M11 & M11D




RPI Part #MID335
OEM Part #MPM60455

DECAL (CAUTION)

- Size: .75" x 1.50"
- Included in: Decal Kit (RPI Part #MIK329)

Fits: Rear Panel
Models: M9, M9D, M11 & M11D



Don't Forget, RPI offers other parts you need most to fit Midmark·Ritter M9 & M11 UltraClaves including:

- | | | | |
|----------------------|--------------------------|------------------------------|--------------------------------|
| • Condensation Coils | • Door Gaskets | • Heater Element Assemblies | • Sterilizer PM Kits |
| • Control PC Boards | • Door Handles | • Rack & tray Kits | • Top Cover Kits |
| • Dam Gaskets | • Door Panel Kits | • Solenoid Valves | • Water Level Sensor Harnesses |
| • Display Assemblies | • Fill/Vent Manifold Kit | • Solenoid Valve Repair Kits | • Many more! |