



... Proud to be a **Healthcare Components Group** company ...



**HEALTHCARE
Components Group**

STERRAD

New Oil Return Solenoid Valves to fit **STERRAD® NX® & 100NX®** Hydrogen Peroxide Sterilization Systems!



RPI Part #SDV105
OEM Part #28-51971-001

SOLENOID VALVE (OIL RETURN)

- 12VDC, 7W (80Ω); 24" lg. wires with 2 pin connector
- Material: Stainless steel body; Viton® (FKM) seals
- Includes: (2 pcs) Male Connector (1/4" Push-In x 1/8" MPT) (RPI Part #RXF117)
- Also available: Solenoid Valve Repair Kit (RPI Part #SDK107) and Coil (RPI Part #SDC108)

Fits: Vacuum Pump Chassis
Models: NX



RPI Part #SDV106
OEM Part #28-51498-001

SOLENOID VALVE (OIL RETURN)

- 12VDC, 7W (80Ω); 24" lg. wires with 2 pin connector
- Material: Stainless steel body; Viton® (FKM) seals
- Includes: Vacuum fitting, face seal O-Ring (RPI Part #RXO091) and (1 pc) Male Connector (1/4" Push-In x 1/8" MPT) (RPI Part #RXF117)
- Also available: Solenoid Valve Repair Kit (RPI Part #SDK107) and Coil (RPI Part #SDC108)

Fits: Vacuum Pump Chassis
Models: 100NX



RPI Part #SDC108
OEM Part # (No OEM Part # Available)

COIL

- 12VDC, 7W (80Ω); 24" lg. wires with (2) installed contacts and 2 pin connector (not installed)
- Also available: Solenoid Valve (Oil Return) (RPI Part #SDV105), Solenoid Valve (Oil Return) (RPI Part #SDV106), and Solenoid Valve Repair Kit (RPI Part #SDK107)

Fits: Solenoid Valve (Oil Return)
Models: NX & 100NX



RPI Part #SDK107
OEM Part # (No OEM Part # Available)

SOLENOID VALVE REPAIR KIT

- Viton® (FKM) seals
- Includes: Spanner wrench
- Also available: Solenoid Valve (Oil Return) (RPI Part #SDV105), Solenoid Valve (Oil Return) (RPI Part #SDV106), and Coil (RPI Part #SDC108)

Fits: Solenoid Valve (Oil Return)
Models: NX & 100NX



RPI Part #RXF117
OEM Part #26-00391-043
**MALE CONNECTOR
(1/4" PUSH-IN X 1/8" MPT)**

- 2 per package
- Material: Nickel plated brass
- Included in: Solenoid Valve (Oil Return) (RPI Part #SDV105) and Solenoid Valve (Oil Return) (RPI Part #SDV106)

Fits: Solenoid Valve (Oil Return)
Models: NX & 100NX

x2



RPI Part #RXO091
OEM Part # (No OEM Part # Available)

O-RING

- 6 per package
- Size: .239" ID x .379" OD x .070" C/S; Duro 75
- Material: Viton® (FKM)
- Included in: Solenoid Valve (Oil Return) (RPI Part #SDV106)

Fits: Solenoid Valve (Oil Return)
Models: 100NX

x6



replacement parts industries, inc. A Healthcare Components Group Company

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. 2109SD Web 01/23

**SECTION
10
PAGE
538
UPDATED
01/23**