


**New Parts to fit Tuttnauer Automatic Autoclaves!**

Models: EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EA/EK/EKA, 2540E/EA/EHS/EK/EKA & 3870EHS




**Introducing the Diagnostic Smart Kit®  
For Servicing Tuttnauer® Automatic Autoclaves  
RPI Part #TUK108 (OEM Part #TEST-9)**

- Calibrate/verify displayed pressure values.
- Calibrate/verify the Temperature Sensor values.
- Check and monitor all test point values from the AJUNC or ANL boards.
- Plus, Kit includes a durable plastic case and a troubleshooting information booklet.



Kit includes all parts below!

**RPI Part #TUT114  
OEM Part #TEST-7**




**SIMULATOR (PT100)**

- For Temperature Sensor calibration (simulates high and low temperatures)
- Use with Wire Harness No. 1, 2, or 3 (RPI Part #TUH111, TUH112, or TUH113) depending on the model number.

**RPI Part #TUH111  
OEM Part # (No OEM Part # Available)**

**WIRE HARNESS (No. 1)**


- Small-sized 2 Pin Connector without jumper
- Use with Simulator (PT100) (RPI Part #TUT114)



**RPI Part #TUH112  
OEM Part # (No OEM Part # Available)**

**WIRE HARNESS (No. 2)**


- Medium-sized 2 Pin Connector without jumper
- Use with Simulator (PT100) (RPI Part #TUT114)



**RPI Part #TUH113  
OEM Part # (No OEM Part # Available)**

**WIRE HARNESS (No. 3)**


- Large-sized 3 Pin Connector with jumper
- Use with Simulator (PT100) (RPI Part #TUT114)



**RPI Part #TUG110  
OEM Part #TEST-2**

**PRESSURE GAUGE (TEST)**


- See reverse side for more information



**RPI Part #TUB109  
OEM Part #TEST-1**

**TEST POINT BOARD**

- Use with Ribbon Cable (RPI Part #TUC117)



**RPI Part #TUC117  
OEM Part # (No OEM Part # Available)**

**RIBBON CABLE**


- Use with Test Point Board (RPI Part #TUB109)



**RPI Part #RPT460  
OEM Part # (No OEM Part # Available)**

**TRIM POT ADJUSTMENT TOOL**


- 5" Long; Plastic with stainless steel blades



**RPI Part #RPT113  
OEM Part # (No OEM Part # Available)**

**MAX REGISTER THERMOMETER**


- Use for temperature verification in sterilizers



**RPI Part #RPC476  
OEM Part # (No OEM Part # Available)**

**CARRYING CASE**

- Includes "Pick-n-Pluck" foam insert




**replacement parts industries, inc.**

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Tuttnauer equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Web 05/10

# Table Top Sterilizers

Models: EZ9, EZ10, EZ10k, 1730, 2340, 2540 & 3870 Series




RPI Part #TUG110  
OEM Part #TEST-2

### PRESSURE GAUGE (TEST)

- Includes all parts listed below:
  - Pressure Gauge (RPI Part #RPG693)**
    - Dual scale: 0-60 PSI/0-400 kpa
    - Maximum indicating pointer for highest PSI reading
  - Gauge Guard (RPI Part #RPG694)**
    - Soft rubber; Color: Blue
  - Male Connector (RPF688)**
    - 1/4 MPT x 3/16" Brass Compression Fitting
  - Tube Adapter (RPI Part #RPF689)**
    - 1/4 MPT x 3/16" Brass Compression Tube
    - Includes (1) piece each: Compression Nut (RPI Part #RPF219) and Compression Sleeve (RPI Part #RPF215)

- All parts also sold separately
- Also included in Diagnostic Smart Kit® (RPI Part #TUK108)


**Fits:** Chamber manifold  
**Models:** EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EA/EK/EKA, 2540E/EA/EK/EKA & 3870EHS



RPI Part #TUT092  
OEM Part #01630101(AJUNC 2 PCB)/  
01630101-1 (AJUNC 3 PCB)

### PRESSURE TRANSDUCER

- Fits units with AJUNC 3 PCB; and can be easily modified to fit units with AJUNC 2 PCB
- Includes Pressure Tube Kit (RPI Part #PRK282)
- Fits:** AJUNC 2 PCB & AJUNC 3 PCB
- Models:** EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EA/EK/EKA, 2540E/EA/EK/EKA & 3870E/EA



RPI Part #TUS093  
OEM Part #01610501

### TEMP SENSOR (PT100)

- Fits: Chamber Manifold AJUNC 3 PCB at the J11 connector
- Models:** EZ9, EZ10, EZ10k, 1730E/EK, 2340E/EA/EKA, & 2540E/EA/EKA

*Machine must be calibrated after replacing Pressure Transducer or Temp Sensor.  
For Diagnostic Smart Kit® see RPI Part #TUK108 on reverse side.*



RPI Part #TUP106  
OEM Part #04400299

### POWER SUPPLY


- Protek brand
- Fits:** Chassis
- Models:** EZ9, EZ10, 1730E, 2340E/EA/EK/EKA, 2540E/EA/EKA & 3870E/EA



RPI Part #TUF103  
OEM Part #FIL175-0006

### FILTER SCREEN (New Style)

- .70" dia. x 1.26" lg.; S.S. mesh
- Fits:** Filter Cap Assembly
- Models:** 1730E/EK, 2340E/EA/EK/EKA, 2540E/EA/EK/EKA & 3870E/EA



RPI PART #	DESCRIPTION
RPC476	CARRYING CASE
RPF688	MALE CONNECTOR
RPF689	TUBE ADAPTER
RPG693	PRESSURE GAUGE
RPG694	GAUGE GUARD
RPK282	PRESSURE TUBE KIT
RPT113	MAX REGISTER THERM.
RPT460	TRIM POT ADJUST. TOOL
TUB109	TEST POINT BOARD
TUC117	RIBBON CABLE

RPI PART #	DESCRIPTION
TUF103	FILTER SCREEN (NEW STYLE)
TUG110	PRESSURE GAUGE (TEST)
TUH111	WIRE HARNESS (NO. 1)
TUH112	WIRE HARNESS (NO. 2)
TUH113	WIRE HARNESS (NO. 3)
TUK108	DIAGNOSTIC SMART KIT®
TUP106	POWER SUPPLY
TUT092	PRESSURE TRANSDUCER
TUT093	TEMP SENSOR (PT100)
TUT114	SIMULATOR (PT100)

Look inside for new  
Diagnostic Smart Kit®  
to fit Tuttnauer  
Automatic Autoclaves

ADDRESS SERVICES REQUESTED

PERMIT NO. 2399  
VAN NUYS, CA  
U.S. POSTAGE PAID  
PRESORTED STANDARD

20338 Corisco Street • Chatsworth, California • 91311  
**replacement parts industries, inc.**

