


## CONDENSATION COILS ... RPI HAS THE ONE YOU NEED!

*Over the years, Tuttnauer has changed the configuration of their reservoirs by reversing the coil.*

*RPI has you covered – we offer all three styles!*

RPI Part #TUC040  
OEM Part #CU836101


Fits EZ9, EZ10,  
1730E/EK & 1730M/MK



RPI Part #TUC041  
OEM Part #CT836101 (Prior to 1993)

Fits 2340E/EK/M/MK  
& 2540E/EK/M/MK  
& Valueklave (1730MKV)

*Coil joins reservoir at the left rear and  
vents towards the front of the machine*



RPI Part #TUC063  
OEM Part #CT836101 (After 1993)

Fits 2340E/EK/M/MK  
& 2540E/EK/M/MK

*Coil joins reservoir at the left front and  
vents towards the rear of the machine*

