



replacement parts industries, inc.

"The Alternate Source"®

CALL 800-221-9723 or 818-882-8611 • FAX 818-882-7028
E-MAIL: techsupport@rpiparts.com • WEBSITE: www.rpiparts.com

RPI Part #SDK021
OEM Part #05-06899-2-200
PM-2 KIT
FITS STERRAD
Model: 100S

Includes:

- (1) Sterile Air HEPA Filter (RPI Part #SSF011)
- (1) Door Seal (RPI Part #SDS018)
- (1) Oil Mist Filter (RPI Part #SDF022)
- (1 pkg of 18) Spacer (Electrode) (RPI Part #SDS023)
- (1 pkg of 2) Shelf Stop Kit (RPI Part #SDK024)

All parts also sold separately

Please Note: This kit does not include the Catalytic Converter.



WARNING - Properly dispose of all collected oil as required by law, and according to local, state and federal regulations.