


ACCESSORIES

RPI PART #	DESCRIPTION	ACCESSORY	DETAILS
LAK001	LEADWIRE ADAPTER KIT	 <div style="display: flex; align-items: center; justify-content: center;"> <div style="display: flex; flex-direction: column; gap: 2px;"> RA V1 V2 V3 V4 V5 V6 LA RL LL </div> <p style="font-size: small; margin: 5px 0;">Includes 10 adapters individually labeled with leadwire identifiers</p> </div>	<ul style="list-style-type: none"> • Includes 10 adapters individually labeled with leadwire identifiers • Adapts banana plugs or pins to tab or snap electrodes • Widest contact area for superior trace quality • Flat, wide bodied design helps prevent twisting • Durable construction - built to last