



LIST OF RPI PARTS TO FIT A-DEC DENTAL CHAIRS			PRIORITY	DECADE	DECADE	DECADE	DECADE	DECADE	DECADE	POSITIONER	CASCADE	PERFORMER	PERFORMER	PERFORMER
RPI Part #	OEM Part #	Description	1005	1010	1011	1015	1020	1021	1030	1040	I	II	III	
ADA213	61.0625.00	Switch Assy (Paddles LT)	*											
ADA214	61.0626.00	Switch Assy (Paddles RT)	*											
ADA215	61.0661.00 & 61.0656.00	Switch Assy (Rockers RT)	*											
ADA216	61.0661.01 & 61.0656.01	Switch Assy (Rockers LT)	*											
ADA228	61.0453.00	Switch Assy (Tilt Limit)	*											
ADB011	023-001-01	Barb - 10/Pkg	•											
ADB203	61.0688.00, 61-0833-00 & 61-0616-00	Circuit Board ('A' Box)	•											
ADB219	61.0645.00	Circuit Board ('B' Box)	•											
ADC035	025-015-01	Sleeve Clamp - 20/Pkg	•	•	•	•	•	•	•	•	•	•	•	
ADC174	61-0860-00	Tilt Cylinder Kit	•											
ADC175	61-0850-00	Lift Cylinder Kit	•											
ADC176	61-1287-00	Lift Cylinder Kit		•	•	•	•	•		•				
ADC177	61-2050-00	Tilt Cylinder								•				
ADC178	61-1267-00	Tilt Cylinder		•	•	•	•	•						
ADC204	61.0612.00	Cable Assy (Main)	*											
ADC205	61.0624.00	Cable Assy (Ctrl Switches)	*											
ADC206	61.0613.00	Cable Assy (Lift Limit)	*											
ADC207	61.0599.00	Cable Assy (Main)	*											
ADC209	61.0693.00	Cable (Pre-Position Base)	*											
ADC210	61.0614.00	Cable Assy (Stop Plate)	*											
ADC226	61.0025.00 (cable) & 006.023.00 (nut)	Brake Cable	*											
ADH217	61.0217.01	Hose Assy (Tilt Cylinder)	•											
ADH218	61.0218.01	Hose Assy (Lift Cylinder)	•											
ADK208	61.0668.00	Cable Assy (Main Complete)	*											
ADK220	60.0601.00	Solenoid Valve Repair Kit	•											
ADK221	041.164.00 & 90.1031.00	Capacitor Kit	•	•	•	•	•	•	•	•			•	
ADK222	61.0030.00 (latch) & 013-017-00 (compression spring)	Latch Kit (Armrest)	•											
ADM225	61.0081.00 (mount) & 010.028.00 (retaining ring)	Mount (Brake Cable)	*											
ADS045	013-022-00	Conical Compression Spring - 10/Pkg	•											
ADS052	001-002-00	Truss Head Slotted Screw - 10/Pkg	•											
ADS062	19-0230-00	Helical Compression Spring - 10/Pkg	•											
ADS211	044.184.01 & 044.049.01	Limit Switch (Modified)	•	•	•		•	•	•	•			•	
ADS212	044.184.00 & 044.049.00	Micro-Switch		•	•		•	•			•	•	•	
ADS223	041.119.00	Switch (Miniature)	•											
ADS224	041.059.00	Switch (Push-Button)	*					•		•				
ADS229	013.017.00	Compression Spring - 10/Pkg	•											
ADW038	004-005-00	Flat Washer - 50/Pkg	•											
RPF449	61-0197-00	Hydraulic Fluid - 16oz. Bottle	•	•	•	•	•	•	•	•	•	•	•	
RPF453	61-0197-00	Hydraulic Fluid (Case - 12 bottles)	•	•	•	•	•	•	•	•	•	•	•	
RPF737	046-070-00	Fuse (0.6A, 250V) - 5/Pkg	•											
RPH726	010.028.00	Retaining Ring (External) - 10/Pkg	*											
RPH727	006.023.00	Nut, Hex (1/4-20) - 10/Pkg	•											
RPO303	030-012-00	O-Ring - 12/Pkg	•											
RPO305	030-010-00	O-Ring - 12/Pkg	•	•	•	•	•	•						
RPO312	030-004-00	O-Ring - 12/Pkg	•	•	•	•	•	•						
RPO316	030-015-00	O-Ring - 12/Pkg	*											
RPO317	030-016-00	O-Ring - 12/Pkg	•	•	•	•	•	•						
RPO738	35.041.02	O-Ring - 10/Pkg	*											
RPP032	No OEM # Available	Hospital Grade Male Plug (Right Angle)	•	•	•	•	•	•	•	•	•	•	•	
RPT017	No OEM # Available	Wire Nut - 15/Pkg	•	•	•	•	•	•	•	•	•	•	•	
RPT083	No OEM # Available	Cable Tie (4" White) - 25/Pkg	•	•	•	•	•	•	•	•	•	•	•	
RPT295	No OEM # Available	1/4" Sleeve Tool	•	•	•	•	•	•	•	•	•	•	•	
RPT296	No OEM # Available	Universal Sleeve Tool	•	•	•	•	•	•	•	•	•	•	•	
RPT398	No OEM # Available	Cable Tie (12" Gray) - 25/Pkg	•	•	•	•	•	•	•	•	•	•	•	

* This part does not fit all serial numbers within the stated model, please find the detailed part description inside this flyer or see our website www.rpiparts.com for serial number breakdowns.