




Parts to Fit A-dec Utility Module and Filter Regulators

UTILITY MODULE OEM Part #24-0352-00 Exploded View

LOOK!
Black Valve Body Repair Kit (RPI Part #ADK173)
also fits the delivery unit's optional accessories.




*It's clear!
The RPI Filter Housing
is clear so you can see
the filter inside!*

FILTER REGULATORS

GOOD NEWS!
These parts also fit the Regulators:

- Filter Kit (RPI Part #ADK201)
- Black Valve Body Repair Kit (RPI Part #ADK173)
- White Valve Body Repair Kit (RPI Part #ADK153)



* The Black Valve Body (RPI Part #ADB198) and the White Valve Body (RPI Part #ADB197) have been replaced by the Clear Valve Body (RPI Part #ADB202). RPI's valve body is clear so you can troubleshoot without taking apart the entire unit! It's the RPI Advantage!