

# LIST OF RPI PARTS TO FIT AIR TECHNIQUES PERI-PRO® , PERI-PRO II® & PERI-PRO III® AS OF FEBRUARY, 2009.

RPI PART #	OEM PART #	DESCRIPTION	PERI-PRO		
			I	II	III
<a href="#">ATA617</a>	90520	"L" TRANSFER ARM ASSY	•	•	•
<a href="#">ATA618</a>	90564	"U" TRANSFER ARM	•	•	
<a href="#">ATA625</a>	90530	"U" TRANSFER ARM ASSY	•	•	
<a href="#">ATA626</a>	90550	"U" TRANSFER DRIVE ASSY	•	•	•
<a href="#">ATA655</a>	94530	"L" TRANSFER ARM ASSY			•
<a href="#">ATB600</a>	90548	SHUTTER RELEASE BAR	•	•	•
<a href="#">ATB609</a>	90546	RUBBER BUMPER - 6/PKG	•	•	•
<a href="#">ATB616</a>	90522	TRANSFER ARM BEARING - 3/PKG	•	•	•
<a href="#">ATB621</a>	90523	"L" TRANSFER ARM BUSHING - 6/PKG	•	•	•
<a href="#">ATB622</a>	90524	TRANSFER ARM BUSHING - 6/PKG	•	•	•
<a href="#">ATF633</a>	90102	DRYER FAN BLADE	•	•	•
<a href="#">ATF649</a>	94278	THERMAL FUSE (98° C)			•
<a href="#">ATF651</a>	90278	THERMAL FUSE (109° C)	•	•	
<a href="#">ATG605</a>	90521	TIMER GEAR	•	•	•
<a href="#">ATG611</a>	90995	MAIN DRIVE GEAR	•		
<a href="#">ATG612</a>	91993	DRYER DRIVE GEAR	•	•	
<a href="#">ATG613</a>	90525	"L" TRANSFER ARM GEAR - 3/PKG	•	•	•
<a href="#">ATG614</a>	90519	IDLER GEAR - 5/PKG	•	•	•
<a href="#">ATG615</a>	90518	TRANSFER ARM GEAR - 3/PKG	•	•	•
<a href="#">ATG619</a>	90517	TRANSFER ARM GEAR	•	•	•
<a href="#">ATG620</a>	90516	TRANSFER GEAR - 3/PKG	•	•	•
<a href="#">ATG627</a>	90996	MAIN DRIVE GEAR		•	
<a href="#">ATG650</a>	94995	MAIN DRIVE GEAR			•
<a href="#">ATG653</a>	94501	SPEED REDUCER GEAR			•
<a href="#">ATG659</a>	90555	FILM INLET GRILLE	•	•	•
<a href="#">ATH604</a>	90540	HUB & PIN ASSY	•	•	
<a href="#">ATH634</a>	90610 <small>(Old OEM Part #90228)</small>	HEATING ELEMENT (DRYER)	•	•	
<a href="#">ATH648</a>	94630	HEATING ELEMENT (DRYER)			•
<a href="#">ATH660</a>	94970	HEATER PLATE			•
<a href="#">ATH665</a>	94540	HUB & PIN ASSY			•
<a href="#">ATK628</a>	<small>(No OEM Part # Available)</small>	SHUTTER RELEASE REPAIR KIT	•	•	

RPI PART #	OEM PART #	DESCRIPTION	PERI-PRO		
			I	II	III
<a href="#">ATK629</a>	90988	TRANSPORT REPAIR KIT	•	•	
<a href="#">ATK663</a>	94988	TRANSPORT REPAIR KIT			•
<a href="#">ATK664</a>	<small>(No OEM Part # Available)</small>	SHUTTER RELEASE REPAIR KIT			•
<a href="#">ATM630</a>	90147	MAIN DRIVE MOTOR	•	•	•
<a href="#">ATM631</a>	90223	DRYER DRIVE MOTOR	•	•	
<a href="#">ATM632</a>	90880	DRYER FAN & MOTOR ASSY	•	•	•
<a href="#">ATP602</a>	30391	SPRING PIN - 6/PKG	•	•	•
<a href="#">ATP610</a>	90505	TOP PLATE PIN	•	•	•
<a href="#">ATR606</a>	90543	PLASTIC SNAP RING - 10/PKG	•	•	•
<a href="#">ATR658</a>	90810	FILM RECEPTACLE	•	•	•
<a href="#">ATS601</a>	90545	COMPRESSION SPRING - 6/PKG	•	•	•
<a href="#">ATS603</a>	90512	SHUTTER RETURN SPRING - 3/PKG	•	•	•
<a href="#">ATS607</a>	90502	IDLER GEAR SHAFT	•	•	•
<a href="#">ATS623</a>	90504	TRANSFER ARM SHAFT	•	•	•
<a href="#">ATS624</a>	90506	TRANSFER DRIVE SHAFT	•	•	•
<a href="#">ATS636</a>	90221	POWER SWITCH ASSY	1		
<a href="#">ATS652</a>	<small>(No OEM Part # Available)</small>	SPACER - 3/PKG			•
<a href="#">ATS671</a>	116521	PROCESS SWITCH		2	
<a href="#">ATS672</a>	116520	POWER SWITCH		2	
<a href="#">ATS673</a>	94269 (BLACK) & 90262 (WHITE)	POWER SWITCH	3	4	•
<a href="#">ATS674</a>	94977 <small>(Old OEM Part #94266)</small>	PROCESS SWITCH		4	•
<a href="#">AIT635</a>	90210	THERMOSTAT/BACKET ASSY	5		
<a href="#">ATT662</a>	94975	THERMISTOR PROBE			•
<a href="#">ATV661</a>	94994	DRAIN VALVE			•
<a href="#">PCB134</a>	117327	ELECTRICAL TERMINAL BUSHING - 3/PKG	6	•	•
<a href="#">RPF051</a>	90265	FUSE (5A, 250VAC) - 5/PKG	•	•	•
<a href="#">RPH638</a>	90219	FUSE HOLDER (MAIN)	•	•	•
<a href="#">RPO304</a>	60604	O-RING - 12/PKG			•
<a href="#">RPT604</a>	94124	DRAIN TUBING			•

NEW PARTS ARE HIGHLIGHTED AND IN **BOLD TYPE**

1=Serial #14999 and below  
5=Serial #20479 and below

2= Serial # 15142 and below  
6=Serial #20480 and above

3=Serial #15000 and above

4=Serial #15143 and above