



Castle/Getinge Bulk Sterilizers

Models: 3522, 3523, 3525, 3533, 3622, 3623 & 3633

PM PACKAGES & SPARE PARTS PACKAGES

| RPI PART # | OEM PART # | DESCRIPTION | (QTY) & RPI PART #'S INCLUDED IN KIT |
|------------|-------------------------|---|--|
| CSP036 | 61301602085 | PM PKG (SOLENOID VALVES) Models: 3522, 3523 & 3533 | • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183 & (2) CSK048 |
| CSP037 | 61301602090 | PM PKG (SOLENOID VALVES) Models: 3525 | • Diaphragm Replacement: (1) CSD051 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183, (1) CSK048 & (1) CSK050 |
| CSP038 | 61301602091 | PM PKG (FILTER & SCREEN) Models: 3522, 3523, 3525 & 3533 | • Filter: (1) CSF052 • Screen Filter: (2) CSS053 |
| CSP039 | 61301602087 | PM PKG (STM TRAP & CHK VLV) Models: 3525 | • Check Valve Rebuild Kit: (3) CSK066 • Steam Trap: (2) CST059 |
| CSP040 | 61301602092 | PM PKG (STM TRAP & CHK VLV) Models: 3525 | • Check Valves: (1) CSV056, (1) CSV057 & (1) CSV058 • Check Valve Rebuild Kit: (1) CSK066 |
| CSP068 | 61301602079 | PM PKG (STM TRAP & CHK VLV) Models: 3622 & 3633 | • Check Valve Rebuild Kits: (1) CSK065 & (3) CSK066 • Steam Trap: (1) CST059 • Steam Trap Rebuild Kit: (1) CSK060 |
| CSP069 | No OEM Part # Available | PM PKG (STM TRAP & CHK VLV) Models: 3522, 3523, 3525 & 3533 | • Check Valve Rebuild Kit: (3) CSK066 • Steam Trap Rebuild Kit: (1) CSK067 |
| CSP070 | No OEM Part # Available | PM PKG (STM TRAP & CHK VLV) Models: 3622 & 3633 | • Check Valve Rebuild Kits: (1) CSK065 & (3) CSK066 • Steam Trap Rebuild Kit: (1) CSK060 & (1) CSK067 |
| CSP041 | 61301602510 | SPARE PARTS PACKAGE Models: 3525 | • Check Valves: (1) CSV056, (1) CSV057 & (1) CSV058 • Check Valve Rebuild Kit: (3) CSK066 • Diaphragm Replacement: (1) CSD051 • Door Gasket: (1) CSG002 • Filter: (1) CSF052 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183, (3) CSK048 & (1) CSK050 • Screen, Detergent: (1) CSS054 • Solid State Relay: (1) CSR028 & (1) CSR029 |
| CSP042 | 61301602498 | SPARE PARTS PACKAGE Models: 3522 | • Check Valve Rebuild Kit: (3) CSK066 • Door Gasket: (1) CSG002 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177 (1) AMK183 & (2) CSK048 • Solid State Relay: (1) CSR028 & (1) CSR029 • Steam Trap: (2) CST059 |
| CSP043 | 61301602499 | SPARE PARTS PACKAGE Models: 3523 | • Check Valve Rebuild Kit: (3) CSK066 • Door Gasket: (1) CSG002 • Filter: (1) CSF052 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183 & (3) CSK048 • Solid State Relay: (1) CSR028 & (1) CSR029 • Steam Trap: (3) CST059 |
| CSP044 | 61301602500 | SPARE PARTS PACKAGE Models: 3533 | • Check Valve Rebuild Kit: (3) CSK066 • Door Gasket: (1) CSG002 • Filter: (1) CSF052 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183, (1) CSK018 & (2) CSK048 • Solid State Relay: (1) CSR028 & (1) CSR029 • Steam Trap: (2) CST059 |
| CSP045 | 61301602501 | SPARE PARTS PACKAGE Models: 3522 | • Check Valve Rebuild Kit: (3) CSK066 • Door Gasket: (1) CSG001 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183 & (2) CSK048 • Solid State Relay: (1) CSR028 & (1) CSR029 • Steam Trap: (2) CST059 |
| CSP046 | 61301602502 | SPARE PARTS PACKAGE Models: 3523 | • Check Valve Rebuild Kit: (3) CSK066 • Door Gasket: (1) CSG001 • Filter: (1) CSF052 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183 & (3) CSK048 • Solid State Relay: (1) CSR028 & (1) CSR029 • Steam Trap: (3) CST059 |
| CSP047 | 61301602503 | SPARE PARTS PACKAGE Models: 3533 | • Check Valve Rebuild Kit: (3) CSK066 • Door Gasket: (1) CSG001 • Filter: (1) CSF052 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183, (1) CSK018 & (2) CSK048 • Solid State Relay: (1) CSR028 & (1) CSR029 • Steam Trap: (2) CST059 |