

## MOTORS & BRUSHES TO FIT CLAY ADAMS CENTRIFUGES

MODEL	MOTOR REPLACEMENT	RPI MOTOR BRUSH & CAP REPLACEMENT		OEM MOTOR BRUSH & CAP REPLACEMENT	
		Brush	Cap	Brush	Cap
Autocrit	CAM115	CAB018	RPC036	CAB018	RPC036
Autocrit II	CAM115	CAB018	RPC036	CAB018	RPC036
Autocrit Ultra III	CAM115	CAB018	RPC036	CAB018	RPC036
Compact II	CAM102				
Dynac	CAM026	CAB018	RPC036*	CAB019	
Dynac II		CAB018			
Micro-Hematocrit	CAM115	CAB018	RPC036	CAB018	RPC036
Micro-Hematocrit II	CAM115	CAB018	RPC036	CAB018	RPC036
Micro-Hematocrit III	CAM115	CAB018	RPC036	CAB018	RPC036
Readacrit	CAM024	CAB022	RPC039	CAB022	RPC039
Readacrit MP	CAM024	CAB022	RPC039	CAB022	RPC039
Safeguard		CAB018			
Serofuge					
Serofuge II	CAM101	CAB083			
Triac	CAM025	CAB018	RPC036	CAB018	RPC036
Triac II	CAM027	CAB018	RPC036	CAB018	RPC036

*\* This Cap (RPI Part #RPC036) only fits the RPI Motor (RPI Part #CAM026). This Cap does not fit the OEM Motor.*