Air Techniques, Apollo/Midmark®, DentalEZ®/Custom Air/Jun-Air/RamVac, Matrx/Midmark®, Midmark® & Tech West



OIL-LESS COMPRESSOR PM KITS (continued on page 2)

<i>OEM</i> Models	RPI Part #	OEM Part #	Parts Included In Compressor PM Kit	
AIR TECHNIQUES AirStar 1 AirStar 10 AirStar 2-1 AirStar 2-2 AirStar 3 AirStar 21 AirStar 22	CMK141	87351 87352	Air Filter Element Kit (RPI Part #CMK065), Filter Element (RPI Part #CME102) and Filter Element (RPI Part #CME103) (2) Air Filter Element Kit (RPI Part #CMK065), Filter Element (RPI Part #CME102) and Filter Element (RPI Part #CME103)	RP0829 RP0830 CME102 CME103
AirStar 30 AirStar 40 AirStar 40 NEO	CMK199	87355	(3) Air Filter Element Kit (RPI Part #CMK065), Filter Element (RPI Part #CME102) and Filter Element (RPI Part #CME103)	
AirStar 5 AirStar 50	CMK143	87353	(4) Air Filter Element Kit (RPI Part #CMK065), Filter Element (RPI Part #CME102) and Filter Element (RPI Part #CME103)	
AirStar 7 AirStar 70	CMK144	87354	(6) Air Filter Element Kit (RPI Part #CMK065), (2) Filter Element (RPI Part #CME102) and (2) Filter Element (RPI Part #CME103)	
AirStar 100M	CMK278	87356	(8) Air Filter Element Kit (RPI Part #CMK065), (2) Filter Element (RPI Part #CME102) and (2) Filter Element (RPI Part #CME103)	CMK065



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

OIL-LESS COMPRESSOR PM KITS (continued from page 1)

<i>OEM</i> Models	RPI Part #	OEM Part #	Parts Included In Compressor PM Kit	
APOLLO/MIDMARK®				
AOCOS11/D/DQ/Q AOCOS12/D/DQ/Q AOCOS21/D/DQ/Q AOCOS22/D/DQ/Q AOCRT21DQ/Q	CMK170	ACA85403 & ACA85417	(2) Replacement Filter Element (RPI Part #CME004) and Coalescing Filter Element (RPI Part #CMK024)	O-RING CME004
AOCRT22DQ/Q R22 (with 10 gal tank)				RP0685 RP0685
AOCOS31/D/DO/O AOCOS32/D/DO/O AOCOS42/D/DO/O AOCOT32/D/DO/O AOCOT42/D/DO/O	CMK171	ACA85404	(4) Replacement Filter Element (RPI Part #CME004) and Coalescing Filter Element (RPI Part #CMK025)	
AOCOT62/D/DQ/Q AOCOT82/D/DQ/Q	CMK172	ACA85406	(8) Replacement Filter Element (RPI Part #CME004) and Coalescing Filter Element (RPI Part #CMK025)	CME043 CME044
AOCRS11 (Manufactured before Apollo was acquired by Midmark) AOCRS12 (Manufactured before Apollo was acquired by Midmark)	CMK270	No OEM Part # Available	Air Filter Element (RPI Part #CME288) and Particle Filter Element (RPI Part #CME085)	CME067
AOCRS11 (Manufactured after Apollo was acquired by Midmark) AOCRS12 (Manufactured after Apollo was acquired by Midmark)	CMK271	No OEM Part # Available	Replacement Filter Element (RPI Part #CME004) and Particle Filter Element (RPI Part #CME085)	RP0685
AOCRR42DQ R22 (with 20, 30, or 40 gal tank) R42	CMK273	ACA85418	(3) Replacement Filter Element (RPI Part #CME004) and Coalescing Filter Element (RPI Part #CMK025)	CME074 RP0831
AOCOS52 AOCOS52D	CMK274 ¹	ACA85415	Air Intake Filter Element (RPI Part #CME281) and Coalescing Filter Element (RPI Part #CMK025)	CMERCE
AOCOS52 AOCOS52D	CMK2861	ACA85416	Air Intake Filter Element (RPI Part #CME281), Coalescing Filter Element (RPI Part #CMK025) and Compressor Tip Seal Kit (RPI Part #CMK287)	CME085 RP0837
DENTALEZ®/CUSTOM AIR/JU				
1013D 1023D 1025D 1027D	CMK076 ²	004016	(2) Air Filter Element (RPI Part #CME067), Coalescing Filter Element (RPI Part #CME043) and Coalescing Filter Element (RPI Part #CME044)	CME105
1022	CMK092 ²	004015	Air Filter Element (RPI Part #CME067), (Coalescing Filter Element (RPI Part #CME043) and Coalescing Filter Element (RPI Part #CME044)	CME107 CME281
1029T	CMK101 ²	004017	(3) Air Filter Element (RPI Part #CME067), Coalescing Filter Element (RPI Part #CME043) and Coalescing Filter Element (RPI Part #CME044)	
CA-812 CA-822 CA-912 CA-922	CMK145	65468179	Intake Filter Element (RPI Part #CME074) and (2) Coalescing Filter Element (RPI Part #CME105)	CME288
CA-922 CA-813-D (SN 440537 and below) CA-823-D (SN 440537 and below)	CMK146	65468180	(2) Air Intake Filter (RPI Part #CMF073) and (2) Coalescing Filter Element (RPI Part #CME105)	CMF073
CA-813/D (SN 440538 - 631602) CA-823/D (SN 440538 - 631602) CA-913 (SN 440538 - 631602) CA-923 (SN 440538 - 631602)	CMK147	65468186	(4) Intake Filter (RPI Part #CMF075) and (2) Coalescing Filter Element (RPI Part #CME105)	RP0676 RP0677
CA-813/D (SN 631603 and above) CA-823/D (SN 631603 and above) CA-913 (SN 631603 and above) CA-923 (SN 631603 and above)	CMK191	004444	(4) Intake Filter (RPI Part #CMF075) and (2) Filter Element (RPI Part #CME107)	
CA-825/D/Q (SN 631602 and below) CA-925/D (SN 631602 and below)	CMK148	65468181	(2) Intake Filter Element (RPI Part #CME074) and (2) Coalescing Filter Element (RPI Part #CME105)	CMK024 CMK025
CA-825/D/Q (SN 631603 and above) CA-925/D (SN 631603 and above)	CMK192	004445	(2) Intake Filter Element (RPI Part #CME074) and (2) Filter Element (RPI Part #CME107)	
CA-827-D (SN 631602 and below) CA-927-D (SN 631602 and below)	CMK149	65468182	(3) Intake Filter Element (RPI Part #CME074) and (2) Coalescing Filter Element (RPI Part #CME105)	
CA-827-D (SN 631603 and above) CA-927-D (SN 631603 and above)	CMK193	004446	(3) Intake Filter Element (RPI Part #CME074) and (2) Filter Element (RPI Part #CME107)	RXT002 RXG070
CA-829-D (SN 631602 and below) CA-929-D (SN 631602 and below)	CMK150	65468183	(4) Intake Filter Element (RPI Part #CME074) and (2) Coalescing Filter Element (RPI Part #CME105) (continued on page 3)	000
Footnotes: ¹ This Kit only fits Compres ² Fits only Dental EZ/Custor			l Air Intake Upgrade Kit (Round Style Housing) (RPI Part #CMK280)	CMK287

OIL-LESS COMPRESSOR PM KITS (continued from page 2)

OFF CEN	RPI	OEM	(continued from page 2)	
Models	Part #	Part #	Parts Included In Compressor PM Kit	
DENTALEZ®/CUSTOM AIR/JU				
CA-829-D (SN 631603 and above)	CMK194	004449	(4) Intake Filter Element (RPI Part #CME074)	
CA-929-D (SN 631603 and above) CA-8212-T (SN 631602 and below)	CMK151	65468184	and (2) Filter Element (RPI Part #CME107) (6) Intake Filter Element (RPI Part #CME074) and	— O-RING
CA-0212-1 (SN 631602 and below)	CIVIN 13 I	03400104	(2) Coalescing Filter Element (RPI Part #CME105)	CME004
CA-8212-T (SN 631603 and above)	CMK195	004447	(6) Intake Filter Element (RPI Part #CME074) and	CIVIEU04
CA OZ 1Z-1 (SN 031003 and above)	OWIN 175	001111	(2) Filter Element (RPI Part #CME107)	
Osprey 13	CMK223	005147SP	(2) Air Filter Element (RPI Part #CME067),	
Osprey 22			Coalescing Filter Element (RPI Part #CME227) and	CME067 CME069
Osprey 23			Particulate Filter Element (RPI Part #CME228)	- OMEGO?
Osprey 24	CMK226	005350SP	(3) Air Filter Element (RPI Part #CME067),	
			Coalescing Filter Element RPI Part #CME227) and	
			Particulate Filter Element (RPI Part #CME228)	
Osprey 25	CMK224	005148SP	(4) Air Filter Element (RPI Part #CME067),	
			Coalescing Filter Element (RPI Part #CME227) and	
	0111/005	00544000	Particulate Filter Element (RPI Part #CME228)	
Osprey 28	CMK225	005149SP	(6) Air Filter Element (RPI Part #CME067),	
			Coalescing Filter Element (RPI Part #CME227) and	
MATDY/MDMADI/®			Particulate Filter Element (RPI Part #CME228)	CME071 CME074
<i>MATRX/MIDMARK</i> ® 0L-100	CMK161	No OEM Part # Available	(3) Foam Element (RPI Part #CMF013) and	RP0831 RP0837
OL-100 OL-101	CIVIN 101	NU UEIVI PARI # AVAIIADIE	Coalescing Filter Element (RPI Part #CMK024)	
0L-100-2	CMK162	No OEM Part # Available	(5) Foam Element (RPI Part #CMF013) and	
OL 100 Z	CIVILLION	NO OLIVI FAIL # AVAIIADIC	Coalescing Filter Element (RPI Part #CMK024)	
0L-100-3	CMK163	No OEM Part # Available	(6) Foam Element (RPI Part #CMF013) and	
02 100 0	0	no ozmi art ii manabio	Coalescing Filter Element (RPI Part #CMK024)	
0L-3	CMK164	No OEM Part # Available	(2) Air Filter Element (RPI Part #CME071) and	
0L-2000			Coalescing Filter Element (RPI Part #CMK024)	OMETOE OMETOZ
0L-2001			,	CME105 CME107
0L-7	CMK165	No OEM Part # Available	(4) Air Filter Element (RPI Part #CME071) and	DD0023 CMC2E4
OL-4000			Coalescing Filter Element (RPI Part #CMK024)	RP0832 CMG254
0L-10	CMK166	No OEM Part # Available	(6) Air Filter Element (RPI Part #CME071) and	
OL-6000			Coalescing Filter Element (RPI Part #CMK024)	
MIDMARK®				
P21	CMK155	77001636	(2) Replacement Filter Element (RPI Part #CME004)	
P22			and Coalescing Filter Element (RPI Part #CMK025)	
P32	CMK198	77001636	(3) Replacement Filter Element (RPI Part #CME004)	
DEG	014/45/	77004/40	and Coalescing Filter Element (RPI Part #CMK025)	CME110 CME227
P52	CMK156	77001610	(4) Air Filter Element (RPI Part #CME071) and	
P72	CMK157	77001611	Coalescing Filter Element (RPI Part #CMK025) (6) Air Filter Element (RPI Part #CME071) and	CMG254
P12	CIVIN 137	77001011	Coalescing Filter Element (RPI Part #CMK025)	
S52	CMK2741	ACA85415	Air Intake Filter Element (RPI Part #CME281) and	
332	GIVINZ/4	A0A03413	Coalescing Filter Element (RPI Part #CMK025)	CME281
S52	CMK2861	ACA85416	Air Intake Filter Element (RPI Part #CME281),	
002	OWNZOO	710/100 110	Coalescing Filter Element (RPI Part #CMK025)	
			and Compressor Tip Seal Kit (RPI Part #CMK287)	
TECH WEST			· · · · · · · · · · · · · · · · · · ·	CMESSE
ACO2S1/ND (2005 and earlier)	CMK174	CSKO-1	(2) Air Filter Element (RPI Part #CME067) and	CME228 CMF013
ACO2S2/ND (2005 and earlier)			Coalescing Filter Element (RPI Part #CMK025)	RP0676 RP0677
ACO2S1/ND (2006 and later)	CMK175	CSKO-1P	(2) Air Filter Element (RPI Part #CME067) and	
ACO2S2/ND (2006 and later)			Coalescing Filter Element (RPI Part #CME110)	
ACO4S2/ND (2005 and earlier)	CMK176	CSKO-2	(4) Air Filter Element (RPI Part #CME067) and	
ACO4D2/ND (2005 and earlier)		001/0	Coalescing Filter Element (RPI Part #CMK025)	
ACO4D2/ND (2006 and later)	CMK177	CSKO-2P	(4) Air Filter Element (RPI Part #CME067) and	
ACO/TO/ND	CM//170	0010 0	Coalescing Filter Element (RPI Part #CME110)	
ACO6T2/ND (2005 and earlier)	CMK178	CSKO-3	(6) Air Filter Element (RPI Part #CME067) and	CMK024 CMK025
ACO6T2/ND (2006 and later)	CMK179	CSKO-3P	Coalescing Filter Element (RPI Part #CMK025) (6) Air Filter Element (RPI Part #CME067) and	
ACOUT STAID (2009 and later)	CIVIN 1/9	OSKU-SF	Coalescing Filter Element (RPI Part #CME1007) and	
ACOR2D1	CMK275	No OEM Part # Available	Air Filter Element (RPI Part #CME067) and	
ACOR2D1Q	OIVII\Z I J	IND OLIVI FAIL# AVAIIADIC	Coalescing Filter Element (RPI Part #CME110)	
ACOR2D12			Sound on Element (the Fruit # ON/LTTO)	DVIOO2 DVOOZO
ACOR2D2Q				DVT000 DV0076
ACOR3T2				RXT002 RXG070
ACOR3T2Q				
ACOR4Q2	CMK277	No OEM Part # Available	Air Filter Element (RPI Part #CME069) and	
			Coalescing Filter Element (RPI Part #CME110)	
Footnotes: 1 This Kit only fits Compr	essors that have	been upgraded with the F	RPI Air Intake Upgrade Kit (Round Style Housing) (RPI Part #CMK280)	CMK287
,				4

LUBRICATED COMPRESSOR PM KITS

<i>OEM</i> Models	RPI Part #	OEM Part #	Parts Included In Compressor PM Kit	
APOLLO/MIDMARK® ALCQL82/D	CMK173	ACA85407	(4) Replacement Filter Element (RPI Part #CME004), Coalescing Filter Element (RPI Part #CMK025), (4) Compressor Oil (RPI Part #CML001) and	O-RING
ALCRL62/D	CMK160	ACA85402	(1 pc) Disposal Bag (RPI Part #RPB838) (3) Replacement Filter Element (RPI Part #CME004), Coalescing Filter Element (RPI Part #CMK025), (3) Compressor Oil (RPI Part #CML001) and (1 pc) Disposal Bag (RPI Part #RPB838)	CME004 RP0831 RP0831
ALCSL11/D ALCSL12/D ALCSL21/D ALCSL22/D	CMK167	ACA85400	Replacement Filter Element (RPI Part #CME004), Coalescing Filter Element (RPI Part #CMK024), Compressor Oil (RPI Part #CML001) and (1 pc) Disposal Bag (RPI Part #RPB838)	O-RING
ALCTL31/D ALCTL32/D ALCTL41/D ALCTL42/D	CMK159	ACA85401	(2) Replacement Filter Element (RPI Part #CME004), Coalescing Filter Element (RPI Part #CMK025), (2) Compressor Oil (RPI Part #CML001) and (1 pc) Disposal Bag (RPI Part #RPB838)	CME104 CME105
DENTALEZ®/CUSTOM AI CA-711 CA-721	R/JUN-AIR/R/ CMK152	4MVAC 65468176	Air Intake Filter (RPI Part #CMF073), Carbon Filter Element (RPI Part #CME104), Coalescing Filter Element (RPI Part #CME105), (2) Synthetic Compressor Oil (RPI Part #CML119), (1 pc) Disposal Bag (RPI Part #RPB838), O-Ring (Rib Cover Gasket) (RPI Part #RP0848) and (1 pc) Drain Plug Gasket (RPI Part #CMG197)	CMF013 CMF073
CA-713 CA-723	CMK153	65468177	(2) Air Intake Filter (RPI Part #CMF073), Carbon Filter Element (RPI Part #CME104), Coalescing Filter Element (RPI Part #CME105), (4) Synthetic Compressor Oil (RPI Part #CML119), (1 pc) Disposal Bag (RPI Part #RPB838), (2) O-Ring (Rib Cover Gasket) (RPI Part #RP0848) and (2 pcs) Drain Plug Gasket (RPI Part #CMG197)	CMG197
CA-725	CMK154	65468178	(4) Air Intake Filter (RPI Part #CMF073), Carbon Filter Element (RPI Part #CME104), Coalescing Filter Element (RPI Part #CME105), (6) Synthetic Compressor Oil (RPI Part #CML119), (1 pc) Disposal Bag (RPI Part #RPB838), (4) O-Ring (Rib Cover Gasket) (RPI Part #RP0848) and (4 pcs) Drain Plug Gasket (RPI Part #CMG197)	CMK024 CMK025
MATRX/MIDMARK®				
AMD-2 AMD-2-115	CMK185	No OEM Part # Available	(3) Foam Element (RPI Part #CMF013), Coalescing Filter Element (RPI Part #CMK024), Compressor Oil (RPI Part #CML001) and (1 pc) Disposal Bag (RPI Part #RPB838)	
AMD-4	CMK186	No OEM Part # Available	(3) Foam Element (RPI Part #CMF013), Coalescing Filter Element (RPI Part #CMK024), (2) Compressor Oil (RPI Part #CML001) and (1 pc) Disposal Bag (RPI Part #RPB838)	CML001 CML119
AMD-6	CMK187	No OEM Part # Available	(4) Foam Element (RPI Part #CMF013), Coalescing Filter Element (RPI Part #CMK024), (3) Compressor Oil (RPI Part #CML001) and (1 pc) Disposal Bag (RPI Part #RPB838)	GINESO! GINE!!
AMD-100 AMD-101	CMK167	77000958	Replacement Filter Element (RPI Part #CME004), Coalescing Filter Element (RPI Part #CMK024), Compressor Oil (RPI Part #CML001) and (1 pc) Disposal Bag (RPI Part #RPB838)	
AMD-100-2	CMK168	77000959	(2) Replacement Filter Element (RPI Part #CME004), Coalescing Filter Element (RPI Part #CMK024), (2) Compressor Oil (RPI Part #CML001) and (1 pc) Disposal Bag (RPI Part #RPB838)	RPB838
AMD-100-3	CMK169	77000960	(3) Replacement Filter Element (RPI Part #CME004), Coalescing Filter Element (RPI Part #CMK024), (3) Compressor Oil (RPI Part #CML001) and (1 pc) Disposal Bag (RPI Part #RPB838)	Not Shown Actual Size
AMD-100 (Manufactured before Midmark acquired Matrx) AMD-102 (Manufactured before Midmark acquired Matrx)	CMK263	No OEM Part # Available	(3) Foam Element (RPI Part #CMF013), Coalescing Filter Element (RPI Part #CMK024), Compressor Oil (RPI Part #CML001) and (1 pc) Disposal Bag (RPI Part #RPB838) (continued on page 5)	RP0848

LUBRICATED COMPRESSOR PM KITS (continued from page 4)


<i>OEM</i> Models	RPI Part #	OEM Part #	Parts Included In Compressor PM Kit	
MATRX/MIDMARK® (contin	ued)			
AMD-100-2 (Manufactured	CMK264	No OEM Part # Available	(5) Foam Element (RPI Part #CMF013),	
before Midmark acquired Matrx)			Coalescing Filter Element (RPI Part #CMK024),	0-RING
			(2) Compressor Oil (RPI Part #CML001) and (1 pc) Disposal Bag (RPI Part #RPB838)	0115004
AMD-100-3 (Manufactured	CMK265	No OFM Port # Available	(6) Foam Element (RPI Part #CMF013),	CME004
before Midmark acquired Matrx)	CIVINZOS	No OEM Part # Available	Coalescing Filter Element (RPI Part #CMK024),	
before witdiffank acquired watrk)			(3) Compressor Oil (RPI Part #CML001) and	
			(1 pc) Disposal Bag (RPI Part #RPB838)	
AC-075	CMK266	No OEM Part # Available	(3) Foam Element (RPI Part #CMF013),	
AC-100			Compressor Oil (RPI Part #CML001) and	CME067 CME069
AC-101			(1 pc) Disposal Bag (RPI Part #RPB838)	
AC-100-2	CMK267	No OEM Part # Available	(5) Foam Element (RPI Part #CMF013),	
			(2) Compressor Oil (RPI Part #CML001) and	
10D 07F	0141/0/0		(1 pc) Disposal Bag (RPI Part #RPB838)	
ACD-075	CMK268	No OEM Part # Available	(2) Foam Element (RPI Part #CMF014),	
ACD-100 ACD-101			Coalescing Filter Element (RPI Part #CMK024), Compressor Oil (RPI Part #CML001) and	CMF013 CMF014
ACD-101			(1 pc) Disposal Bag (RPI Part #RPB838)	CIVII 013 CIVII 014
ACD-100-2	CMK269	No OEM Part # Available	(4) Foam Element (RPI Part #RFB036)	
7100 100 2	OMITEO	NO CENT LIT # MUNICIPAL	Coalescing Filter Element (RPI Part #CMK024),	RP0832
			(2) Compressor Oil (RPI Part #CML001) and	KF0632
			(1 pc) Disposal Bag (RPI Part #RPB838)	
MIDMARK®				
CL21	CMK158	77001633	Replacement Filter Element (RPI Part #CME004),	
CL22			Coalescing Filter Element (RPI Part #CMK025),	
			Compressor Oil (RPI Part #CML001) and	
			(1 pc) Disposal Bag (RPI Part #RPB838)	
CL32	CMK159	77001634	(2) Replacement Filter Element (RPI Part #CME004),	CME110
			Coalescing Filter Element (RPI Part #CMK025),	
			(2) Compressor Oil (RPI Part #CML001) and (1 pc) Disposal Bag (RPI Part #RPB838)	
CL52	CMK160	77001635	(3) Replacement Filter Element (RPI Part #CME004),	— RP0676 RP0677
0102	OWNCTOO	77001033	Coalescing Filter Element (RPI Part #CMK025),	
			(3) Compressor Oil (RPI Part #CML001) and	
			(1 pc) Disposal Bag (RPI Part #RPB838)	
TECH WEST				
ACL2S1/ND (2005 and earlier)	CMK188	CSKL-1	Air Filter Element (RPI Part #CME067),	
ACL2S2/ND (2005 and earlier)			Coalescing Filter Element (RPI Part #CMK025),	CMK024 CMK025
			Compressor Oil (RPI Part #CML001) and	
			(1 pc) Disposal Bag (RPI Part #RPB838)	
ACL2S1 (2006 and later)	CMK182	CSKL-1P	Air Filter Element (RPI Part #CME067),	
ACL2S2 (2006 and later)			Coalescing Filter Element (RPI Part #CME110),	
			Compressor Oil (RPI Part #CML001) and (1 pc) Disposal Bag (RPI Part #RPB838)	
ACL4D2/ND (2005 and earlier)	CMK189	CSKL-2	Air Filter Element (RPI Part #CME067),	
ACCADZ/ND (2003 and earlier)	CIVIL 107	OSINE-Z	Coalescing Filter Element (RPI Part #CMK025),	To continue the second of the continue the c
			(2) Compressor Oil (RPI Part #CML001) and	
			(1 pc) Disposal Bag (RPI Part #RPB838)	
ACL4D2 (2006 and later)	CMK183	CSKL-2P	Air Filter Element (RPI Part #CME067),	CML001
			Coalescing Filter Element (RPI Part #CME110),	
			(2) Compressor Oil (RPI Part #CML001) and	
101 (70):-		-	(1 pc) Disposal Bag (RPI Part #RPB838)	
ACL6T2/ND (2005 and earlier)	CMK190	CSKL-3	Air Filter Element (RPI Part #CME069),	and
			Coalescing Filter Element (RPI Part #CMK025),	7
			(3) Compressor Oil (RPI Part #CML001) and	
ACL6T2 (2006 and later)	CMK184	CSKL-3P	(1 pc) Disposal Bag (RPI Part #RPB838) Air Filter Element (RPI Part #CME069),	
MOLUTZ (ZUU6 AND TATET)	CIVIN 104	CONT-OL	Coalescing Filter Element (RPI Part #CME110),	
			(3) Compressor Oil (RPI Part #CML001) and	RPB838
			(1 pc) Disposal Bag (RPI Part #RPB838)	
<u> </u>			() () (Free () () () () () () () () () (

Don't Forget RPI's Free Dental Compressor PM Guide!

RPI's easy-to-use Dental Compressor PM Guide walks you through the steps of what to look for and what to replace at monthly and annual intervals. What's more, this PM Guide is easy to personalize with your company's name and phone number which makes a perfect leave behind for your customers as a handy reminder! Download your free copy of RPI's PM Guide from our website, *www.rpiparts.com* today!



ESSENTIAL PM TOOLS FOR SCROLL STYLE COMPRESSORS!




RPI Part #RXK069 OEM Part #HPA30530

GREASE GUN KIT

- OEM recommends regreasing of compressor head and motor bearings every 5,000 hours
- Includes: Grease Gun (RPI Part #RXG071), Grease Gun Nozzle (RPI Part #RXN072), Grease Cartridge (RPI Part #RXG070)
- Used to grease the bearings in the scroll compressor head.

Fits: Scroll Compressor Head

Models: Apollo/Midmark: AOCOS52 & AOCOS52D; Midmark: S52



4 ESSENTIAL TOOLS FOR AN EASIER OIL CHANGE!

SYRINGE PUMP/EXTRACTOR (RPI Part #RPT842)

This one-of-a-kind tool allows complete removal of oil from a compressor without needing to tilt or disassemble the unit!

2 DRAIN VALVE (RPI Part # RPV841)

Upgrade to RPI's one-of-a-kind Drain Valve which has a special Tee-handle that allows you to easily open and close the valve to control the flow of oil as it is draining into the Disposal Baq!


DISPOSAL BAG (RPI Part #RPB838)

The perfect solution for carrying, storing and disposing of used oil. This one gallon bag is constructed of durable clear plastic with an expandable base, a built in handle, and tight sealing cap.

COLLECTION CONTAINER KIT (RPI Part #'s CMK180 & CMK181)

These one-of-kind collection bottles can replace the purge drain lines and be directly connected with tubing to the compressor so that the purge waste water and oil can be collected and disposed of properly according to local, state and federal regulations.

- For Pressurized Purge Systems use RPI Part #CMK180
- For Gravity Drain Systems use RPI Part #CMK181



Don't Forget ...

ABSORBENT PAD (RPI Part #RPP958)

This 9-1/4" x 14-1/2" pad is perfect for oil drip collection and spills and is sold in a package of 3!

