



RPI Parts to fit Pelton & Crane Magnaclave Sterilizers

RPI PART #	OEM PART #	DESCRIPTION
PCA118	004231	Locator Arm - 3/pkg
PCB001	004048	Bellows
PCB088	004491	Buzzer
PCB237	004800 (Bellows Only)	Bellows (Pressure Switch)
PCC005	004228	Bellows Cap
PCC073	004030	Main Valve Cap
PCC241	004598 (Cap) & 004864 (O-Ring)	Drain Cap With O-Ring
PCF009	1881031 & 004326	Filter (Fill-Line)
PCG019	004497	Door Gasket
PCG040	3336356 & 015316	Pressure Gauge
PCG079	4207796 & 015363	Temperature Gauge
PCH083	004502	Main Heater Strip (240VAC)
PCH089	004516	Wall Heater Strip (208VAC)
PCH091	004515	Main Heater Strip (208VAC)
PCH092	004503	Wall Heater Strip (240VAC)
PCK026	017706 (OEM has only 1 set)	Main Valve Kit
PCK216	004497 (Door Gasket Only)	Sterilizer PM Kit
PCK224	011047	Bellows Kit
PCK226	<i>OEM Part # Not Available</i>	Filter with Cotter Pin
PCK236	5473041	Solenoid Kit (Door Latch)
PCK238	004630 & A65-004630	Condensation Coil Kit
PCK240	004598, 004864, 004599, 005929 & 091007	Drain Kit
PCK242	<i>OEM Part # Not Available</i>	Bushing Kit (Door Hinge)
PCL029	004501	Lamp (Open Door)
PCL030	004499	Lamp (Heat On)
PCL031	004500	Lamp (Sterile)
PCL032	004498	Lamp (Power On)

RPI PART #	OEM PART #	DESCRIPTION
PCN072	004027	Main Valve Nut
PCN136	004076	Speednut - 10/pkg
PCR037	004039	Valve Button Retainer
PCR196	004463	Latching Relay
PCR234	041002 & A62-011081	Relay (Elect. Contactor)
PCS043	004079	Function Switch (3 Lead)
PCS044	004073	Power Switch (2 Lead)
PCS074	004028	Main Valve Stem
PCS086	004468	Door Open Switch
PCS087	004461	1 LB. Pressure Switch
PCS090	004545	Door Interlock Switch
PCS094	004478	Pressure Switch
PCS235	004837 & B02-004837	Relay Socket
PCS239	004525 & A63-004525	Spring (Door Interlock)
PCT013	004640	Fill Vent Tube
PCT081	004531 & 041001	Electric Timer (60 Minute)
PCT082	004527	Wall Thermostat
PCT197	004505	Overheat Thermostat Switch
PCV075	004310	Main Valve Body
PCV109	004529	Pressure Relief Valve
RPC442	<i>OEM Part # Not Available</i>	Uniclamp - 6/pkg
RPF217	<i>OEM Part # Not Available</i>	Compression Sleeve - 5/pkg
RPF221	<i>OEM Part # Not Available</i>	Compression Nut - 5/pkg
RPH285	<i>OEM Part # Not Available</i>	Drain Hose - sold by the foot
RPO354	004004	Air Valve O-Ring - 6/pkg
RPO355	004000	Main Valve O-Ring - 6/pkg
RPO386	004864	O-Ring - 6/pkg
RPT113	<i>OEM Part # Not Available</i>	Max Register Thermometer

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Pelton & Crane equipment. All product names used in this document are trademarks or registered trademarks of their respective holders.