



More New Parts to fit Midmark • Ritter 7 & M7 SpeedClave® Sterilizers




INDEPENDENT BUZZER-STYLE TIMER


- Use this Kit to replace or upgrade a Timer with an independent buzzer.
(Use *Time Switch* (RPI Part #RCT048) if you need a timer with a built-in bell.)
- Timer Knob sold separately – see RPI Part #'s RCK111 and RCK092.
- The Wire Harness in this Kit may be needed to extend the length of the original Wire Harness when upgrading to this style timer.

RPI Part #RCK116
OEM Part #002-0373-00

TIMER KIT (30 MIN)
Fits: Front Panel
Model: M7



Includes all parts as shown





RPI Part #RCK111
OEM Part #053-1008-01

TIMER KNOB
Fits: Timer
Model: M7

RPI Part #RCK092
OEM Part #H98062

TIMER KNOB
Fits: Timer
Models: 7 & M7




RPI Part #RCT107
OEM Part # (No OEM Part # Available)

TIMER (30 MIN)
• See "Independent Buzzer-Style" note above with Timer Kit (RPI Part #RCK116)
Fits: Front Panel
Model: M7



RPI Part #RCB096
New OEM Part #015-1813-00
Old OEM Part #015-0745-02

BUZZER
Fits: Machine Base
Model: M7



RPI Part #RCH117
OEM Part # (No OEM Part # Available)


WIRE HARNESS
• See "Independent Buzzer-Style" note above with Timer Kit (RPI Part #RCK116)
Fits: Original Wire Harness
Model: M7

BUILT-IN BELL-STYLE TIMER

- Use this Timer to replace a Timer with a built-in bell. This timer will not work in conjunction with machines that utilize an independent buzzer. (Use *Timer Kit* (RPI Part #RCK116) if you need a timer with an independent buzzer.)

RPI Part #RCT048
OEM Part #H284667-1

TIME SWITCH
Fits: Front Panel
Models: 7 & M7-without buzzer




Includes all parts as shown




RPI Part #RPK213
OEM Part # (No OEM Part # Available)

KNOB
Fits: Time Switch (RPI Part #RCT048) only
Models: 7 & M7




RPI Part #RCA099
OEM Part #H98064
FLEX CABLE ASSEMBLY
Fits: Front Panel and Relay Assembly
Models: 7 & M7

Includes all parts as shown



RPI Part #RCK110
OEM Part #053-1008-00
TEMPERATURE KNOB
Fits: Flex Cable Assembly
Model: M7



RPI Part #RCK093
OEM Part #H98063
CONTROL KNOB
Fits: Flex Cable Assembly
Models: 7 & M7



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Midmark-Ritter equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. WEB flyer with newsletter 04/09




More New Parts to fit Midmark•Ritter 7 & M7 SpeedClave® Sterilizers

RPI Part #RCH118
Old RPI Part #RCH024
OEM Part #002-0839-00, H284687-1
& 002-0245-00

HEATER ELEMENT ASSEMBLY


- Includes all parts as shown
- See **Service Tip**
- Fits:** Chamber
- Models:** 7 & M7

Attaching Hardware (RPI Part #RCH119)
 for the Heater Element Assembly (RPI Part #RCH024) is also sold separately. Includes (2) each Fiber Washer, Flat Washer, Internal Tooth Washer and Hex Nut



Service Tip

To help prevent clogging of the 7 or M7, install RPI Filter (RPI Part #MIF062) into the chamber end of the Water Fill/Exhaust line.




RPI Part #MIF062
New OEM Part #002-0360-00 (Mesh Style)
Old OEM Part #016-0790-00

FILL/VENT MESH CHAMBER FILTER

- See **Service Tip**
- Fits:** Chamber
- Models:** 7 & M7


RPI Part #RCK106
NEW OEM Part #002-0248-00
OLD OEM Part #053-0513-00

DOOR HANDLE KIT
Fits: Door Assembly
Models: 7 & M7



RPI Part #RCP120
OEM Part #042-0005-16

CAM PIN
Fits: Door Handle Assembly
Models: 7 & M7




RPI Part #RCN103
OEM Part #K906034

PUSH CAP NUT (5/PKG)
Fits: Cam Pin (RPI Part #RCP120)
Models: 7 & M7




RPI Part #MIT047
New OEM Part #002-0370-00
Old OEM Part #015-0647-00

OVERTEMP THERMOSTAT
Fits: Bracket that mounts to under side of chamber
Model: M7




RPI Part #RCT050
OEM Part #H225639

MANUAL THERMOSTAT
Fits: Bracket that mounts to under side of chamber
Models: 7 & M7



RPI Part #RCL121
OEM Part #015-1304-00

PILOT LIGHT (GREEN)
Fits: Front Panel
Model: M7



Don't forget, RPI also carries these parts to fit Midmark•Ritter 7 & M7 sterilizers and more!

RPI PART #	OEM PART #	DESCRIPTION	RPI PART #	OEM PART #	DESCRIPTION
PCC003	H98159	INDUSTRIAL GRADE POWER CORD (7 & M7)	RCS012	040-0010-103, H96673 & 3046	HINGE SCREW (7 & M7)
RCB002	H284709	BELLOWS ASSEMBLY (7 Only)	RCS042	002-0263-00 & 33053	RELAY SWITCH ASSEMBLY (7 & M7)
RCB089	002-0375-00	AIR VENT BELLOWS ASSEMBLY (7 & M7) (Old OEM Part #014-0193-00)	RCS086	H225730	SLEEVE (7 & M7)
RCB100	002-0375-00	AIR VENT BELLOWS ASSEMBLY (7 & M7) (Old OEM Part #014-0193-00)	RCS098	002-0247-00	VALVE STEM KIT (7 & M7)
RCC010	34646	CAM (7 Only)	RCS108	H106330	DOOR SPRING (7 & M7)
RCC095	H98009	FILL COVER (7 & M7)	RCS115	015-1791-00	FILL/VENT SWITCH (M7 Only)
RCC114	014-0236-10	COIL (FILL/VENT) (M7 Only)	RCV051	H98089 & 11683	SAFETY VALVE (31 PSI) (7 & M7)
RCF016	H98021	SUCTION CUP FOOT (7 & M7)	RCV112	014-0236-00	SOLENOID VALVE (FILL/VENT) (M7 Only)
RCG019	H98134 & 33166	DOOR GASKET (7 & M7)	RCW009	H96771 & 3059	DOOR WASHER (7 & M7)
RCG085	002-0242-00	TEMPERATURE GAUGE (DIAL THERMOMETER) (7 & M7)	RPB167	H98160	STRAIN RELIEF BUSHING (7 & M7)
RCG087	H98135	RESERVOIR GASKET (7 & M7)	RPF426	053-0990-00	FEMALE QUICK CONNECT FITTING (GREY) (M7 Only)
RCG104	053-1579-00	FIBER GASKET (7 & M7)	RPF427	053-0990-01	MALE QUICK CONNECT FITTING (GREY) (M7 Only)
RCK113	014-0236-11	REPAIR KIT (FILL/VENT) (M7 Only)	RPF679	015-0346-35	FUSE (12A, 250VAC) - FAST ACTING (M7)
RCL090	H98036	PILOT LIGHT (RED) (7 & M7)	RPH674	015-1259-04	FUSE HOLDER (M7 Only)
RCN013	041-0010-08 & 3047	HINGE SCREW NUT (7 & M7)	RPK431	002-0653-00	DRAIN TUBE ASSEMBLY KIT (M7)
RCP011	042-0007-07 & H97989	CAM PIN (7 & M7)	RPO306	H98136 & 36906	O-RING (7 & M7)

See the RPI Catalog or website for more parts to fit the 7 & M7 sterilizers!

Note: Parts fit model(s) noted to the right of the descriptions above.



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Midmark•Ritter equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. WEB flyer with newsletter 04/09