




Models: 421105 Nutator, 421106 Nutator, H5500 Mini LabRoller™, S0500 GyroMini™ & S0600 Enduro™ MiniMix™

# RPI Now Offers a New Line Of Parts To Fit LABORATORY NUTATORS AND ROTATORS!




**RPI Part #CAM118**  
OEM Part #42110601/421507/NE110

- MOTOR (20 RPM, 230V)**
- Bi-directional operation
  - Offset shaft



**Fits:** Motor Housing  
**Model:** 421106 Nutator




**RPI Part #CAM116**  
OEM Part #42110501/4210504/NE100

- MOTOR (24 RPM, 115V)**
- Bi-directional operation
  - Offset shaft



**Fits:** Motor Housing  
**Models:** 421105 Nutator & S0500 GyroMini




**RPI Part #CAM117**  
OEM Part #NE111

- MOTOR (18 RPM, 115V)**
- Bi-directional operation
  - Offset shaft



**Fits:** Motor Housing  
**Model:** S0600 Enduro MiniMix




**RPI Part #CAS122**  
OEM Part #NE210

**POWER SWITCH**

- Rocker switch
- 16A @ 250VAC; SPST; .187" Tabs
- Color: Black Body with white I/O marking

**Fits:** Chassis  
**Models:** S0500 GyroMini, S0600 Enduro MiniMix & H5500 Mini LabRoller




**RPI Part #CAS119**  
OEM Part #42110502/421505

**POWER SWITCH**

- Rotary style: On-Off-On-Off
- 3A@125VAC; SPST; Wire length: ~6" lg.

**Fits:** Chassis  
**Model:** 421105 Nutator




**RPI Part #RPB866**  
OEM Part # (No OEM Part # Available)

**BUMPER/FOOT**

- 4 per package
- Size: .81" Square x .52" Ht.
- Material: Grey Polyurethane
- Self adhesive backing

**Fits:** Base  
**Model:** 421105 Nutator




**RPI Part #CAK121**  
OEM Part # (No OEM Part # Available)

**POWER CORD KIT**

- 220VAC; Length: ~8 ft.
- Also available: Strain Relief Bushing Pliers (RPI Part #RPT172)

**Fits:** Chassis  
**Model:** 421106 Nutator



**RPI Part #CAK120**  
OEM Part #42110503/421503

**POWER CORD KIT**

- 115VAC; Length: ~8 ft.
- Also available: Strain Relief Bushing Pliers (RPI Part #RPT172)

**Fits:** Chassis  
**Model:** 421105 Nutator



**replacement parts industries, inc.**

20338 Corisco Street • Chatsworth, California • 91311

PRESORTED STANDARD  
U.S. POSTAGE PAID  
VAN NUYS, CA  
PERMIT NO. 2399

ADDRESS SERVICES REQUESTED

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

# RPI Mixes Things Up With Our New Line Of Parts To Fit Laboratory Nutators and Rotators!

## NUTATORS

New Parts to fit Clay Adams 421105 & 421106 Nutators including:

- Motor (20 RPM, 230V)
- Motor (24 RPM, 115V)
- Power Cord Kit
- Power Switch
- Bumper/Foot

New Parts to fit Labnet International, Inc. S0500 GyroMini™ & S0600 Enduro™ MiniMix™ Nutators including:

- Motor (18 RPM, 115V)
- Motor (24 RPM, 115V)
- Power Switch

## ROTATORS

New Part to fit Labnet International, Inc. H5500 Mini LabRoller™ Rotators including:

- Power Switch

