



**replacement parts industries, inc.**

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

# PM Kits for Table Top Sterilizers

RPI Part #	Fits Models	Parts Included in PM Kits
<b>AMSCO</b>		
AMK234	Eagle® 10 & 10 Plus	Door Gasket (RPI Part #AMG161), Water Guard (RPI Part #AMG172), Filter (Fill Line) (RPI Part #PCF009)
<b>MIDMARK®-RITTER®</b>		
MIK072	UltraClave® M9 & M9D (door without metal flap)	Dam Gasket (RPI Part #MIG035), Door Gasket (RPI Part #MIG027), (2 pcs) Door Spring (RPI Part #MIS076), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Filter (Sintered Style) (RPI Part #MIF038)
MIK080	UltraClave® M11 & M11D	Dam Gasket (RPI Part #MIG036), Door Gasket (RPI Part #MIG028), (2 pcs) Door Spring (RPI Part #MIS076), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Filter (Sintered Style) (RPI Part #MIF038)
MIK134	SpeedClave® M7 & 7	Door Gasket (RPI Part #RCG019), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062)
MIK158	UltraClave® M9 & M9D (door with metal flap)	Dam Gasket (RPI Part #MIG157), Door Gasket (RPI Part #MIG027), (2 pcs) Door Spring (RPI Part #MIS076), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Filter (Sintered Style) (RPI Part #MIF038)
<b>PELTON &amp; CRANE</b>		
PCK216	Magnaclave	Door Gasket (RPI Part #PCG019), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), Filter with Cotter Pin (RPI Part #PCK226), Filter (Fill Line) (RPI Part #PCF009)
PCK217	Sentry	Door Gasket (RPI Part #PCG020), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), Filter (Fill Line) (RPI Part #PCF009)
PCK218	OCM	Door Gasket (RPI Part #PCG014), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), Filter (Fill Line) (RPI Part #PCF009)
PCK219	OCR and OCR+	Door Gasket (RPI Part #PCG015), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), Filter (Fill Line) (RPI Part #PCF009)
PCK220	Validator 8	Door Gasket (RPI Part #PCG065), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), Filter (Fill Line) (RPI Part #PCF009), Spacer (Door) (RPI Part #PCS067)
PCK221	Validator 10	Door Gasket (RPI Part #PCG066), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), Filter (Fill Line) (RPI Part #PCF009), Spacer (Door) (RPI Part #PCS067)
PCK222	Validator Plus 8	Door Gasket (RPI Part #PCG065), Filter (Fill Line) (RPI Part #PCF009), Reservoir Filter Kit (RPI Part #PCK167), Spacer (Door) (RPI Part #PCS067)
PCK223	Validator Plus 10	Door Gasket (RPI Part #PCG066), Filter (Fill Line) (RPI Part #PCF009), Reservoir Filter Kit (RPI Part #PCK167), Spacer (Door) (RPI Part #PCS067)
PCK231	Delta AE8, XL8 and Q8	Door Gasket (RPI Part #PCG065), Compressor Filter (RPI Part #SCF002), Filter (Fill Line) (RPI Part #PCF009), Filter In-Line Assembly (RPI Part #PCF233), HEPA Filter (RPI Part #RPF644), Reservoir Filter Kit (RPI Part #PCK167), Spacer (Door) (RPI Part #PCS067)
PCK232	Delta AF10, XL10 and Q10	Door Gasket (RPI Part #PCG066), Compressor Filter (RPI Part #SCF002), Filter (Fill Line) (RPI Part #PCF009), Filter In-Line Assembly (RPI Part #PCF233), HEPA Filter (RPI Part #RPF644), Reservoir Filter Kit (RPI Part #PCK167), Spacer (Door) (RPI Part #PCS067)
<b>SCICAN®</b>		
SCK035	Statim® 2000	Cassette Seal (RPI Part #SCS001), Biological Filter (RPI Part #SCF032), Compressor Filter (RPI Part #SCF002), Liquid Soap (RPI Part #RPS287), Water Filter (In-Line) (RPI Part #SCF034)
SCK036	Statim® 5000	Cassette Seal (RPI Part #SCS029), Air Compressor Filter (RPI Part #SCF031), Biological Filter (RPI Part #SCF032), Liquid Soap (RPI Part #RPS287), Water Filter (In-Line) (RPI Part #SCF034)
<b>TUTTNAUER®</b>		
TUK121	1730E/M	Door Gasket (RPI Part #TUG001), Air Jet Valve (Black Top) (RPI Part #TUJ034), Door Bellows (RPI Part #TUB019), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Rubber Boot for Door Switch (RPI Part #TUB072)
TUK122	1730EK/MK/MKV	Door Gasket (RPI Part #TUG001), Air Jet Valve (Red Top) (RPI Part #TUJ033), Door Bellows (RPI Part #TUB019), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Rubber Boot for Door Switch (RPI Part #TUB072)
TUK123	2340M (Serial #8805 and below)	Door Gasket (RPI Part #TUG002), Air Jet Valve (Black Top) (RPI Part #TUJ034), Door Bellows (RPI Part #TUB019), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Rubber Boot for Door Switch (RPI Part #TUB072)
TUK124	2340E & 2340M (Serial #8806 and above)	Door Gasket (RPI Part #TUG021), Air Jet Valve (Black Top) (RPI Part #TUJ034), Door Bellows (RPI Part #TUB019), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Rubber Boot for Door Switch (RPI Part #TUB072)
TUK125	EZ9 & 2340EA	Door Gasket (RPI Part #TUG021), Air Jet Valve (Black Top) (RPI Part #TUJ034), Door Bellows (RPI Part #TUB019), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), HEPA Filter (RPI Part #RPF644), Rubber Boot for Door Switch (RPI Part #TUB072)
TUK126	2340EK/MK	Door Gasket (RPI Part #TUG021), Air Jet Valve (Red Top) (RPI Part #TUJ033), Door Bellows (RPI Part #TUB019), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Rubber Boot for Door Switch (RPI Part #TUB072)
TUK127	2340EKA	Door Gasket (RPI Part #TUG021), Air Jet Valve (Red Top) (RPI Part #TUJ033), Door Bellows (RPI Part #TUB019), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), HEPA Filter (RPI Part #RPF644), Rubber Boot for Door Switch (RPI Part #TUB072)
TUK128	2540E/M	Door Gasket (RPI Part #TUG003), Air Jet Valve (Black Top) (RPI Part #TUJ034), Door Bellows (RPI Part #TUB019), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Rubber Boot for Door Switch (RPI Part #TUB072)
TUK129	EZ10k & 2540EKA	Door Gasket (RPI Part #TUG003), Air Jet Valve (Red Top) (RPI Part #TUJ033), Door Bellows (RPI Part #TUB019), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), HEPA Filter (RPI Part #RPF644), Rubber Boot for Door Switch (RPI Part #TUB072)
TUK130	2540EK/MK	Door Gasket (RPI Part #TUG003), Air Jet Valve (Red Top) (RPI Part #TUJ033), Door Bellows (RPI Part #TUB019), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Rubber Boot for Door Switch (RPI Part #TUB072)
TUK131	3850E/M & 3870E/M	Door Gasket (RPI Part #TUG074), Air Jet Valve (Red Top) (RPI Part #TUJ033), Door Bellows (RPI Part #TUB019), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Rubber Boot for Door Switch (RPI Part #TUB072)
TUK132	3850EA & 3870EA	Door Gasket (RPI Part #TUG074), Air Jet Valve (Red Top) (RPI Part #TUJ033), Door Bellows (RPI Part #TUB019), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), HEPA Filter (RPI Part #RPF644), Rubber Boot for Door Switch (RPI Part #TUB072)
TUK133	EZ10 & 2540EA/EHS	Door Gasket (RPI Part #TUG003), Air Jet Valve (Black Top) (RPI Part #TUJ034), Door Bellows (RPI Part #TUB019), (2 pcs) Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), HEPA Filter (RPI Part #RPF644), Rubber Boot for Door Switch (RPI Part #TUB072)