

# Listing of parts to fit Pelton & Crane Spirit® Delivery Units

RPI PART #	OEM PART #	DESCRIPTION	SPIRIT	SPIRIT II	SPIRIT IIE	SPIRIT AC	SPIRIT EC
PCA752	33 03 042/51 74 052	MICRO-SWITCH & DIAPHRAGM ASSEMBLY	•	•	•		•
PCC756	94 33 525	IN-LINE COUPLING (1/8")	•	•	•		•
PCC757	94 33 517	MALE QUICK CONNECT (1/8")	•	•	•		•
PCC772	94 33 541	FEMALE QUICK CONNECT (1/16")				•	
PCC773	94 33 533	MALE QUICK CONNECT (1/16")				•	
PCC777	(No OEM Part # Available)	WATER SOLENOID COIL	•	•	•		•
PCC790	33 02 762	HANDPIECE BLOCK CONNECTOR (8/pkg)	•	•	•		•
PCD748	33 02 978	REGULATOR DIAPHRAGM	•	•	•		•
PCD749	33 02 853	FOOT SWITCH DIAPHRAGM	•	•	•		•
PCG760	30 04 484	PRESSURE GAUGE				•	
PCI753	020892	INSERT (SYRINGE)	•	•	•	•	•
PCK732	021814 (Assy.) <sup>1</sup>	HANDPIECE HANGER ASSEMBLY (BLACK)	•	•			
PCK741	021813 (Assy.) <sup>2</sup>	SYRINGE HANGER ASSEMBLY (BLACK)	•	•			
PCK754	(No OEM Part # Available)	WATER SOLENOID REPAIR KIT	•	•	•		•
PCK759	15 32 535	HANDPIECE HANGER ASSEMBLY					7
PCK775	51 51 779	ANTI-SIPHON VALVE KIT				5	
PCK776	33 22 315 (Assy.) <sup>3</sup>	HANDPIECE HANGER ASSEMBLY (WHITE)			•		6
PCK778	33 39 889 (Assy.) <sup>4</sup>	SYRINGE HANGER ASSEMBLY (WHITE)			•		6
PCK780	(No OEM Part # Available)	FOOT CONTROL VALVE KIT	•	•	•		•
PCK781	(No OEM Part # Available)	FOOT CONTROLLER SPACER KIT	•	•	•		•
PCK783	33 39 970	HANDPIECE HANGER ASSEMBLY				•	
PCK784	(No OEM Part # Available)	FOOT CONTROLLER REPAIR KIT	•	•	•		•
PCL768	30 03 402	AIR INDICATOR (GREEN)				•	
PCS750	33 02 879	CUSHION SPRING	•	•	•		•
PCS779	33 03 067	STEM	•	•	•		•
PCS788	021833	AIR VALVE STEM	•	•	•		•
PCS789	020177	WATER VALVE STEM	•	•	•		•
PCS791	94 47 897	WATER VALVE STEM			•		•
PCT785	023725	FOOT CONTROLLER TUBING ASSEMBLY (8')	•	•	•		•
PCT786	(No OEM Part # Available)	FOOT CONTROLLER TUBING (25')	•	•	•		•
PCV755	33 03 000	CHECK VALVE ASSEMBLY	•	•	•		•
PCV765	(No OEM Part # Available)	VALVE STEM				•	
PCV766	15 33 764	VALVE TOGGLE SWITCH				•	
PCV767	18 98 902	CHECK VALVE				•	
PCV787	33 07 431	AUTOHOLDER VALVE				•	
RPK628	(No OEM Part # Available)	COMPRESSION NUT & SLEEVE-1/8" (10/pkg)	•	•	•		•
RPO302	(No OEM Part # Available)	O-RING (12/pkg)	•	•	•		•
RPO305	18 97 151	O-RING (12/pkg)	•	•	•		•
RPO310	33 03 075	O-RING (12/pkg)	•	•	•	•	•
RPO312	(No OEM Part # Available)	O-RING (12/pkg)	•	•	•		•
RPO314	51 74 052	O-RING (12/pkg)	•	•	•		•
RPO617	70 32 741	O-RING (12/pkg)				•	

1 = 020282 (Top); 020185 (Bottom)

2 = 020547 (Top); 020185 (Bottom)

3 = 15 35 546 (Top); 15 35 827 (Bottom)

4 = 15 24 813 (Top); 15 35 827 (Bottom)

5 = Fits Spirit AC Serial #4839 and below

6 = Spirit EC prior to October, 1993

7 = Spirit EC after October, 1993

**replacement parts industries, inc.**  
 20338 Corisco Street • Chatsworth, California 91311

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)