


## Safety Valves and Kits to fit Tuttnauer Autoclaves



**SERVICE TIP:** When a Safety Valve needs replacement, replace it with the same rated PSI Valve – in other words, replace a 37 PSI Valve with a 37 PSI Valve, and a 40 with a 40. The PSI cracking pressure is actually etched onto the body of the valve for your reference.

### THE RPI ADVANTAGE! Safety Valve & Kits

- Safety Valve now available for 37 PSI and 40 PSI !
- Updated Safety Valve Kits include both the Elbow Fitting and Threaded Adapter for easy plumbing !

DESCRIPTION & PART #	EZ9/EZ10/EZ10k 1730E/M, 2340E/M, 2540E/M & 3870E/M		1730EK/MK 2340EK/EKA/MK 2540 EK/EKA/MK		VALUEKLAVE 1730 MKV
	37 PSI	40 PSI	37 PSI	40 PSI	40 PSI
<b>AIR JET VALVE</b> (Stainless steel; 10 mm hex) RPI Part # TUJ033 (OEM PART #CB842010) RPI Part # TUJ034 (OEM PART #CT842010)	TUJ034 <i>BLACK TOP</i>	TUJ034 <i>BLACK TOP</i>	TUJ033 <i>RED TOP</i>	TUJ033 <i>RED TOP</i>	TUJ033 <i>RED TOP</i>
<b>SAFETY VALVE HOLDER KIT</b> RPI Part #'s TUK053/TUK054/TUK077/TUK078 OEM PART #'s (No OEM Part #'s Available) • Includes Air Jet Valve, Safety Valve & Safety Valve Holder (see below)	TUK054	TUK078	TUK053	TUK077	TUK077
<b>SAFETY VALVE HOLDER</b> RPI Part #TUH031 (OEM PART #CT841010) RPI Part #TUH032 (OEM PART #CT841020) • Includes Safety Valve Holder, Threaded Adapter, Elbow Fitting & Hardware	TUH032	TUH032	TUH031	TUH031	TUH031
<b>SAFETY VALVE</b> (Brass; 1/4" MPT) RPI Part # TUV011 (OEM PART #03110003) RPI Part # TUV065 (OEM PART #03110002)	TUV011	TUV065	TUV011	TUV065	TUV065
<b>THREADED ADAPTER</b> (Brass) RPI Part #TUA060 OEM PART # (No OEM Part #'s Available)	TUA060	TUA060	TUA060	TUA060	TUA060
<b>ELBOW FITTING</b> (Brass) RPI Part #TUF079 OEM PART # (No OEM Part #'s Available)	TUF079	TUF079	TUF079	TUF079	TUF079
<b>MOUNTING HARDWARE</b> RPI Part #TUK055 OEM PART # (No OEM Part #'s Available)	TUK055	TUK055	TUK055	TUK055	TUK055



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Tuttnauer equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. WEB 02/08