



replacement parts Industries, Inc.

CALL: (800) 221-9723 or (818) 882-8611 • FAX: (818) 882-7028
 techsupport@rpiparts.com • www.rpiparts.com

LIST OF RPI PARTS AVAILABLE TO FIT
 TUTTNAUER MANUAL UNITS AS OF MAY, 2008.

RPI PART #	OEM PART #	DESCRIPTION	SEE FOOTNOTE	VALUE-KLAVE	1730	2340	2540	3870
MDL019	01960007	SIGNAL LIGHT - RED (125VAC)			M	M	M	
MDL021	01910257	SIGNAL LIGHT - AMBER (125VAC)		MKV	M	M	M	
MDL022	01910258	SIGNAL LIGHT - GREEN (125VAC)		MKV	M	M	M	
MIF062	(No OEM Part # Available)	FILL/VENT MESH CHAMBER FILTER		MKV	M MK	M MK	M MK	M
RPC582	02819996	POWER CORD (230VAC)			MK	MK	MK	
RPH285	02620016A	DRAIN HOSE		MKV	M MK	M MK	M MK	
RPO386	02610030/Inner	O-RING (Drain Valve)			M MK	M MK	M MK	
RPO387	02610027/Outer	O-RING (Drain Valve)			M MK	M MK	M MK	
RPR583	02819993	AC INLET RECEPTACLE		MKV	M	M MK	M MK	M
TUA060	(No OEM Part # Available)	THREADED ADAPTER (Fits Safety Valve Holder)		MKV	M MK	M MK	M MK	M
TUA066	CT312036	DOOR SWITCH ACTIVATOR		MKV	M MK	M MK	M MK	M
TUB019	(No OEM Part # Available)	DOOR BELLOWS		MKV	M MK	M MK	M MK	M
TUB023	01910098	CIRCUIT BREAKER (15A)		MKV	M	M	M	
TUB024	01910097	CIRCUIT BREAKER (10A)			MK	MK	MK	
TUB047	01910100	CIRCUIT BREAKER (15A)			M MK	M MK	M MK	M
TUB048	01910099	CIRCUIT BREAKER (7A)			MK			
TUB064	CT245010	DOOR BELLOW HOUSING BOLT		MKV	M MK	M MK	M MK	M
TUB072	(No OEM Part # Available)	RUBBER BOOT (Fits Door Switch)		MKV	M MK	M MK	M MK	M
TUC028	02819995	POWER CORD (125VAC)		MKV	M	M	M	M
TUC040	CJ836101	CONDENSATION COIL			M MK			
TUC041	CT836101	CONDENSATION COIL	1	MKV		M MK	M MK	
TUC063	CT836101	CONDENSATION COIL	1			M MK	M MK	
TUC067	02550019	WATER RESERVOIR/FILTER COVER		MKV	M MK	M MK	M MK	M
TUF079	(No OEM Part # Available)	ELBOW FITTING (Fits Safety Valve Holder)		MKV	M MK	M MK	M MK	M
TUG001	02610020	DOOR GASKET		MKV	M MK			
TUG002	02610005	DOOR GASKET	2			M		
TUG003	02610023	DOOR GASKET					M MK	
TUG012	02300011	PRESSURE GAUGE W/ INDICATOR (2-1/2")			M MK	M MK	M MK	M
TUG020	02300012	PRESSURE GAUGE (1-1/2")		MKV				
TUG021	02610118	DOOR GASKET	3			M MK		
TUG022	02610029	WATER RESERVOIR GASKET		MKV	M MK	M MK	M MK	M
TUG074	02610019	DOOR GASKET						M
TUH004	01720011	HEATER ELEMENT (120VAC, 450W)		MKV	M MK			
TUH005	01720002	HEATER ELEMENT (120VAC, 350W)				M		
TUH006	01720003	HEATER ELEMENT (120VAC, 350W)					M	
TUH015	(No OEM Part # Available)	ATTACHING HARDWARE (Fits Heaters)		MKV	M MK	M MK	M MK	M
TUH016	01720012	HEATER ELEMENT (230VAC, 450W)			MK			
TUH017	01720013	HEATER ELEMENT (230VAC, 550W)				MK		
TUH018	01720014	HEATER ELEMENT (230VAC, 550W)					MK	
TUH027	01720001	HEATER ELEMENT (120VAC, 350W)			M			
TUH031	CT841010	SAFETY VALVE HOLDER	4	MKV	M MK	MK	MK	
TUH032	CT841020	SAFETY VALVE HOLDER	4		M	M	M	M
TUH043	CJ900012	WIRE HARNESS	5		M MK			
TUH044	CT900012	WIRE HARNESS	6			M MK	M MK	
TUJ033	CB842010	AIR JET VALVE (Red Top)	4	MKV	MK	MK	MK	
TUJ034	CT842010	AIR JET VALVE (Black Top)	4		M	M	M	M
TUK030	CT241111	DOOR BELLOWS ASSEMBLY KIT		MKV	M MK	M MK	M MK	M
TUK037	(No OEM Part # Available)	MULTI-PURPOSE VALVE REPAIR KIT			M MK	M MK	M MK	M
TUK049	02450002	TIMER & MULTI-PURPOSE VALVE KNOB		MKV	M MK	M MK	M MK	M
TUK050	02450003	THERMOSTAT KNOB		MKV	M MK	M MK	M MK	M
TUK053	(No OEM Part # Available)	SAFETY VALVE HOLDER KIT (37 PSI)	4		MK	MK	MK	
TUK054	(No OEM Part # Available)	SAFETY VALVE HOLDER KIT (37 PSI)	4		M	M	M	M
TUK055	(No OEM Part # Available)	MOUNTING HARDWARE (Fits Safety Valve Holder)		MKV	M MK	M MK	M MK	M
TUK061	01910197/Switch Only	MICRO-SWITCH KIT		MKV	M MK	M MK	M MK	M
TUK075	(No OEM Part # Available)	RESERVOIR COVER WITH DIPSTICK		MKV	M MK	M MK	M MK	M
TUK077	(No OEM Part # Available)	SAFETY VALVE HOLDER KIT (40 PSI)	4	MKV	MK	MK		
TUK078	(No OEM Part # Available)	SAFETY VALVE HOLDER KIT (40 PSI)	4		M	M	M	M
TUS009	01910172	POWER SWITCH		MKV	M MK	M MK	M MK	M
TUS013	01910191	MICRO-SWITCH (CHERRY E13/BLACK)			M MK	M MK	M MK	M
TUS014	01910190	DOOR SWITCH		MKV	M MK	M MK	M MK	M
TUS035	01620301	OVER TEMPERATURE SAFETY SWITCH	7		MK	MK	MK	
TUS057	01910197	MICRO-SWITCH (CHERRY D48X/GREY)		MKV	M MK	M MK	M MK	M
TUS068	2550043	WATER RESERVOIR DIP STICK		MKV	M MK	M MK	M MK	M
TUT007	01620101	CONTROL THERMOSTAT		MKV	M MK	M MK	M MK	M
TUT008	01910011	TIMER (Old OEM Part #01910005)			M MK	M MK	M MK	M
TUT038	01620004	SAFETY THERMOSTAT (MANUAL RESET)			M MK	M MK	M MK	M
TUT039	01620103	SAFETY THERMOSTAT (AUTOMATIC RESET)			M MK	M MK	M MK	M
TUT073	TEST-3	TAP (3/8-19 BSPP)		MKV	M MK	M MK	M MK	M
TUV011	03110003	SAFETY VALVE (37 PSI)	4		M MK	M MK	M MK	M
TUV025	CT810013	MULTI-PURPOSE VALVE ASSEMBLY			M MK	M MK	M MK	M
TUV042	CT844180	DRAIN VALVE ASSEMBLY		MKV	M MK	M MK	M MK	M
TUV065	03110002	SAFETY VALVE (40 PSI)	4	MKV	M MK	M MK	M MK	M

IMPORTANT NOTE: Over the years, Tuttnauer has substituted parts from what has been noted in their manuals. As a precaution, please verify parts before replacing or servicing them.

FOOTNOTES 1) TUC041 & TUC063: Both parts fit Models 2340M/MK & 2540M. However, see page 7, Condensation Coils, for specifics. 2) TUG002: Fits Models 2340M S/N 8805 and below. 3) TUG021: Fits Models 2340M/MK S/N 8806 and above. 4) TUJ033/TUJ034, TUH031/TUH032, TUK053/TUK054, TUK077/TUK078 and TUV011/TUV065: See page 7, Safety Valves & Safety Valve Holders. 5) TUH043: For Models 1730M/MK, existing connections should accommodate most units manufactured after January 1993. However, some modifications may be necessary. 6) TUH044: For Models 2340M/MK and 2540 M/MK, existing connections should accommodate most units manufactured after January 1993. However, some modifications may be necessary. 7) TUS035: Fits units prior to 1993.