




replacement parts industries, inc.

Call (800) 221-9723 • Fax (818) 882-7028 • E-mail: order@rpiparts.com • Website: www.rpiparts.com




RPI parts to fit AIR SHIELDS Infant Incubators

This newly redesigned gasket is totally symmetrical, which makes installation easier and more precise!




RPI Part #AIG027
OEM Part #68 120 01
ACCESS DOOR GASKET
Fits: Access Panel



RPI Part #AIT086
OEM Part #26 681 00

RIGHT ANGLE THERMOMETER


- Graduated in Celsius 18° to 40°
 - Graduated in Fahrenheit 60° to 105°
- Fits: Isolette Hood



RPI Part #AIA055
OEM Part #78 931 82

HEATER SOCKET ASSEMBLY

Fits: Heater/Reflector Assembly



RPI Part #AIH040
New OEM Part #78 263 00
Old OEM Part #78 263 10

HEATING ELEMENT

- 120V, 600W, ~24Ω
- Fits: Heater Connector Assembly


RPI Part#	OEM Part#	DESCRIPTION	C-86	C-100 & C-200	C-100QT & C-200QT	C-300	C-400QT & C450QT	C-500QT & C-500QT Model XL	C-550QT & C-550QT Model XL	ICS	IICS
AIG015	12 216 00	HOOD GASKET	•	•	•		•				
AIK068	17 061 89	INSTRUMENT KNOB		•	•						
AIP063	17 575 14	SKIN TEMPERATURE PROBE	•								
AIS051	17 681 33	SLIDE SWITCH								•	
AIS049	17 681 61	SLIDE SWITCH								•	
AIS052	17 681 87	SLIDE SWITCH								•	
AIS050	17 681 88	SLIDE SWITCH								•	
AIS047	17 682 63	ROCKER SWITCH				•					•
AIS042	17 724 32	BULB SOCKET								•	
AIM021	17 776 59	PANEL METER	•							•	
AMB028	17 800 30	LIGHT BULB	•								
AIC037	17 802 30	LENS CAP (AMBER)	•								
AIC035	17 802 60	LENS CAP (RED)	•								
AIC036	17 804 00	LENS CAP (WHITE)	•								
AIS034	17 806 27	LIGHT SOCKET	•								
AIC043	17 806 44	LENS CAP (WHITE)								•	
AIB041	17 806 45	BULB								•	
AIC044	17 806 46	LENS CAP (RED)								•	
LMP014	17 807 41	INCANDESCENT LAMP									•
AIG017	26 501 00	THERMOSTAT PLUG GASKET	•	•	•		•	•	•		
AIG018	26 503 00	POWER PACK GASKET	•	•	•		•	•	•		
AIG019	26 630 11	ACCESS PANEL GASKET	•								
AID001	26 650 00	ACCESS DOOR	•	•		•					
AIG025	26 654 00	ACCESS DOOR GASKET	•								
AIT086	26 681 00	RIGHT ANGLE THERMOMETER	•								
AIT011	26 928 71	101.5°F THERMOSTAT	•								
AIT012	26 928 72	103°F THERMOSTAT	•								
AIF024*	26 945 70	MICRO FILTER	•	•	•		•	•	•		
AIG016	68 120 00	ACCESS DOOR GASKET		•		•					
AIG027	68 120 01	ACCESS DOOR GASKET			•		•	•			
AIK062	68 156 62	SMALL RETAINER KNOB			•		•				
AIK061	68 156 63	LARGE RETAINER KNOB			•		•				
AIS060	68 157 56	"QT" HEAT SHIELD			•		•				
AIH008	68 205 45	IMPELLER		•	•						
AIH006	68 205 46	IMPELLER	•								
AIH069	68 208 25	HEATER		•	•		•				
AIP064	68 209 70	SKIN TEMPERATURE PROBE		•	•		•				•
AIP070	68 209 80	AIR AUXILIARY PROBE		•	•		•				
AIS067	68 210 75	RESET SWITCH ASSEMBLY		•	•						
AIG028	68 512 10	SEALING STRIP		•	•		•	•			
AIR013	68 903 84	RELAY (Old OEM Part #17 652 78)		•							•
AIK007	68 905 70	MOTOR KIT (SERVO) - RPI "Grey" New Style		•	•						
AIK010	68 905 71	MOTOR KIT - 220/240V		•	•						
AIK005	68 905 72	MOTOR KIT (SERVO) - RPI "Grey" New Style	•								
AIK009	68 925 73	MOTOR KIT - 220/240V	•								
AIH040	78 263 00	HEATING ELEMENT									•
AIH039	78 902 70	HEATING ELEMENT								•	
AIA055	78 931 82	HEATER SOCKET ASSEMBLY									•
AIS048	78 932 00	PUSH-BUTTON SWITCH									•
AIC058	N/A	MOTOR CAPACITOR	•	•	•						
AIC059	N/A	MOTOR CAPACITOR	•	•	•						

*Fits models with Dew-Ette®

RPI Part #AIK007
OEM Part #68 905 70

MOTOR KIT (SERVO)
RPI "GREY" NEW STYLE


- 120V/60 Hz., 1800 RPM
 - This new style motor has a "grey" colored base
- Fits: Chassis Controller



RPI Part #AIK005
OEM Part #68 905 72

MOTOR KIT (SERVO)
RPI "GREY" NEW STYLE

- 120V/60 Hz., 1800 RPM
 - This new style motor has a "grey" colored base
- Fits: Chassis Controller



RPI Part #AIP070
OEM Part #68 209 80


AIR AUXILIARY PROBE
Fits: Controller



RPI Part #AIF024
OEM Part #26 945 70

MICRO FILTER

- 4 per package
 - Totally synthetic media; contains *no* glass fibers
- Fits: Shell and Deck Assembly



RPI Part #AIG025
OEM Part #26 654 00


ACCESS DOOR GASKET
Fits: Access Panel



RPI Part #AIT011
OEM PART #26 928 71
(101.5°F)

RPI Part #AIT012
OEM PART #26 928 72
(103°F)

THERMOSTAT
Fits: Thermostat Capsule



RPI Part #AIP064
OEM Part #68 209 70

SKIN TEMPERATURE PROBE
Fits: Controller

