



RPI Parts to fit Midmark•Ritter 7 & M7 Speedclave™ and M9 & M11 UltraClave™

RPI Part #	OEM Part #	DESCRIPTION	7 SPEEDCLAVE	M7 SPEEDCLAVE	M9 ULTRACLAVE	M11 ULTRACLAVE
RCT050	61854	MANUAL THERMOSTAT, O.T.	•	•		
RCG085	002-0242-00	TEMPERATURE GAUGE (DIAL THERM)	•	•		
RCS098	002-0247-00	VALVE STEM WITH O-RINGS (Fits OEM Fill-Vent Valve)	•	•		
RCS042	002-0263-00, 33053	RELAY SWITCH ASSEMBLY	•	•		
MIP050	002-0357-00	TEMPERATURE PROBE ASSEMBLY			•	•
MIS075	002-0358-00	WATER LEVEL SENSOR ASSY			•	•
PCV052	002-0359-00	PRESSURE RELIEF VALVE			•	•
MIF038	002-0360-00	FILL FILTER			•	•
MIS039	002-0362-00	DOOR SWITCH ASSEMBLY			•	•
MIK072	002-0361-00	STERILIZER PM KIT			•	
MIS044	002-0363-00	PULSE SOLENOID			•	
MIS079	002-0363-03	PULSE SOLENOID				•
MIS042	002-0365-00	SOLENOID VALVE (FILL)			•	•
MIS043	002-0366-00	SOLENOID VALVE (VENT)			•	•
MIH048	002-0367-00	HEATER ELEMENT			•	
MIT047	002-0370-00	OVERTEMP THERMOSTAT			•	•
MIK080	002-0504-00	STERILIZER PM KIT - M11				•
MIH049	002-0505-00	HEATER ELEMENT ASSEMBLY				•
MIV073	002-0519-00	SOLENOID VALVE (VENT)			•	•
MIK078	002-0654-00	STEAM TRAP KIT (New Style)			•	•
MIV081	002-0654-00	STEAM TRAP (New Style)			•	•
MIF062	002-0360-00	FILL/VENT MESH CHAMBER			•	•
RCB089	014-0193-00	AIR VENT DIAPHRAGM ASSEMBLY			•	•
RCB100	014-0193-00	AIR VENT BELLOWS ASSEMBLY			•	•
MIS045	014-0199-00	SOLENOID VALVE (VENT)			•	
MIS046	014-0200-00	SOLENOID VALVE (FILL)			•	
RPF049	015-0346-08	FUSE-PRINTER (3A, 250V)			•	•
RPF281	015-0346-07	FUSE (15 A, 250V)			•	•
RPR290	015-0639-00	AC RECEPTACLE, SNAP-IN			•	•
RPC291	015-0640-00	POWER CORD W/RIGHT ANGLE PLUG			•	•
RPB289	015-1137-00	STRAIN RELIEF BUSHING			•	•
RPC288	015-1139-00	POWER CORD W/CONNECTORS			•	•
MIP041	016-0395-02	DOOR INSULATOR PAD				•
MIF051	016-0523-00	LEVELER FEET			•	•
RCN013	041-0010-08	HINGE SCREW NUT	•	•		
MIS076	050-2347-00	DOOR SPRING			•	•
MIA086	050-2362-00	SPRING ARM			•	
MIA077	050-3970-00	SPRING ARM				•
MIC037	052-0173-00	CONDENSATION COIL			•	•
MIG027	053-0366-00	DOOR GASKET				•
RPI268	053-0385-00	CHAMBER INSULATION			•	
MIG035	053-0508-00	DAM GASKET			•	
MIG028	053-0527-00	DOOR GASKET				•
RPI269	053-0585-00	CHAMBER INSULATION				•
MIG036	053-0784-00	DAM GASKET				•
RCB096	15-0745-00	BUZZER		•		
MIP040	16-0395-01	DOOR INSULATOR PAD				•
RCS086	H225730	RUBBER SLEEVE	•	•		
RCT048	H284667-1	MECHANICAL TIME SWITCH	•	•		
RCH024	H284687-1, 002-0245-00	HEATING ELEMENT ASSEMBLY	•	•		
RCB002	H284709	BELLOWS ASSEMBLY	•	•		
RCS012	H96673, 3046	HINGE SCREW	•	•		
RCW009	H96771, 3059	DOOR WASHER	•	•		
RCB088	H96772	T-BOLT	•	•		
RCB003	H97948, H97076, H97075 H97966/GASKET	THERMAL DIAPHRAGM BELLOWS ASSY	•	•		
RCP011	H97989, 34501	CAM PIN & E-RING	•	•		
RCC010	34646	CAM	•			
RCC095	H98009	FILL COVER	•	•		
RCF016	H98021	SUCTION CUP FOOT	•	•		
RCL090	H98036	PILOT LIGHT	•	•		
RCK092	H98062	TIMER KNOB	•	•		
RCK093	H98063	CONTROL KNOB	•	•		
RCV051	H98089, 11683	PRESSURE RELIEF SAFETY VALVE	•	•		
RCG019	H98134, 33166	DOOR GASKET	•	•		
RCG087	H98135	RESERVIOR GASKET	•	•		
RPO306	H98136, 36906	O-RING (Old RPI Part #RCO039)	•	•		
RPO307	H98138, 146180	O-RING (Old RPI Part #RCO040)	•	•		
PCC004	H98159, 33071	POWER CORD W/O CONNECTORS	•	•		
MIC057	N/A	COIL (FILL VALVE)			•	•
MIC058	N/A	COIL (VENT VALVE)			•	•
MIC059	N/A	COIL (VENT VALVE)			•	
MIC060	N/A	COIL (FILL VALVE)			•	
MIC083	N/A	COIL (VENT VALVE)			•	•
MIH052	N/A	ATTACHING HARDWARE			•	•
MIK053	N/A	REPAIR KIT (FILL VALVE)			•	•
MIK054	N/A	REPAIR KIT (VENT VALVE)			•	•
MIK055	N/A	REPAIR KIT (VENT VALVE)			•	
MIK056	N/A	REPAIR KIT (FILL VALVE)			•	
MIK061	N/A	DRAIN HOSE KIT			•	•
MIK082	N/A	REPAIR KIT (VENT VALVE)			•	•
PCN136	N/A	SPEEDNUT	•	•	•	•
RCH025	N/A	ATTACHING HARDWARE	•	•		
RCS097	N/A	VALVE STEM WITH O-RINGS (Fits RPI Fill-Vent Valve)	•	•		
RPK282	N/A	PRESSURE TUBE KIT			•	•
RPT113	N/A	MAX REG THERMOMETER	•	•	•	•