

RPI Part #	OEM Part #	Description	Section 2 Page #
------------	------------	-------------	------------------

GOMCO

DIAPHRAGM PUMPS

Models: 280, 2801, 3020, 3040 & 6030

● CAS063	01-90-2304	ON/OFF SWITCH	200
----------	------------	---------------	-----

Model 300 Series

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
● GOG046	01-90-3059	VACUUM GAUGE (300 SERIES PORT)	200
● GOK043	01-90-2295	DIAPHRAGM PARTS KIT (300 SERIES)	200
GOS001	01-90-0550	NO. 8 STOPPER	200

All Diaphragm Pump Models

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOS001	01-90-0550	NO. 8 STOPPER	200

ROTARY PUMPS

Models: 400, 402, 4010, 4030, 4040, 4050 & 4070

● CAS063	01-90-2304	ON/OFF SWITCH	200
----------	------------	---------------	-----

All Models with Rotary Pumps with 5/8" Rotor and 1/20 hp Motor

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOL022	01-90-3106	ROTARY PUMP OIL	201
GOS001	01-90-0550	NO. 8 STOPPER	200
● GOS036	01-90-2322	5/8" ROTOR SLIDE	201
● GOW037	01-90-3067	5/8" OIL WICK	201

All Models with Rotary Pumps with 3/4" Rotor and 1/6 hp Motor

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOK017	01-90-2416	3/4" ROTOR SLIDE ASSEMBLY KIT	201
GOL022	01-90-3106	ROTARY PUMP OIL	201
GOR028	01-90-2417	3/4" ROTOR	201
GOS001	01-90-0550	NO. 8 STOPPER	200
GOS029	01-90-2418	3/4" ROTOR SLIDE	201
● GOW035	01-90-2875	3/4" OIL WICK	201

All Models with Rotary Pumps with 1" Rotor and 1/4 or 1/6 hp Motor

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOK017	01-90-2416	3/4" ROTOR SLIDE ASSEMBLY KIT	201
GOL022	01-90-3106	ROTARY PUMP OIL	201
GOR028	01-90-2417	3/4" ROTOR	201
GOS001	01-90-0550	NO. 8 STOPPER	200
● GOS038	01-90-0853	1" ROTOR SLIDE	201
● GOW039	01-90-2876	1" OIL WICK	201

All Models with Rotary Pumps and Overflow Valve

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOD018	01-90-3108	FILTER DISC	201

● NEW PART

RPI Part #	OEM Part #	Description	Section 2 Page #
------------	------------	-------------	------------------

GOMCO (continued)

All Models with Rotary Pumps and Overflow Valve (continued)

GOF002	01-90-3107	FELT FILTER DISC	201
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOG014	01-90-2330	VALVE GASKET	201
GOK027	01-90-2291	ROTARY UNIT PARTS KIT	201
GOS001	01-90-0550	NO. 8 STOPPER	200
GOS007	01-90-2902	VALVE SPRING	201

All Models with Rotary Pumps requiring 1/8 MPT with Center Back Mount

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOG006	01-90-2338	VACUUM GAUGE	201
GOL022	01-90-3106	ROTARY PUMP OIL	201
GOS001	01-90-0550	NO. 8 STOPPER	200

All Models with Rotary Pumps using 600ml Bottle Top

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOK027	01-90-2291	ROTARY UNIT PARTS KIT	201
GOL022	01-90-3106	ROTARY PUMP OIL	201
GOS001	01-90-0550	NO. 8 STOPPER	200

All Models with Rotary Pumps with Trap Bottle

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOL022	01-90-3106	ROTARY PUMP OIL	201
GOS001	01-90-0550	NO. 8 STOPPER	200

All Models with Rotary Pumps incorporating Bacteria Filter as Standard Equipment

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOL022	01-90-3106	ROTARY PUMP OIL	201
GOS001	01-90-0550	NO. 8 STOPPER	200

THERMOTIC PUMPS

Model 270

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOL032	01-90-2310	PUMP CYCLE LIGHT	202
GOS001	01-90-0550	NO. 8 STOPPER	200

Model 765A

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOB016**	01-90-3045	PC BOARD	202
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOG014	01-90-2330	VALVE GASKET	201
GOH004*	01-90-2713	HEATER 115 VAC 149 Ω (OLD STYLE)	202
GOH005**	01-90-2713	HEATER 115 VAC 160 Ω	202
GOL026	01-90-3570	PILOT LAMP	202

(continued)

Aspirators & Pumps

RPI Part #	OEM Part #	Description	Section 2 Page #
------------	------------	-------------	------------------

GOMCO

Model 765A (continued)

GOS001	01-90-0550	NO. 8 STOPPER	200
GOS010*	01-90-3050	MICRO-SWITCH	202
GOS030	(None)	ON/OFF SWITCH	202
GOS031	(None)	HI/LOW SWITCH	202
GOV011	01-90-3433	SUCTION VALVE (INLET)	202
GOV012	01-90-3432	PRESSURE VALVE (OUTLET)	202
GOV013	01-90-8140	AEROVENT RUBBER VALVE	202

* Indicates unit with electro-mechanical control
 ** Indicates unit with electronic control

Model 764

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOB016**	01-90-3045	PC BOARD	202
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOG014	01-90-2330	VALVE GASKET	201
GOH004*	01-90-2713	HEATER 115 VAC 149 Ω (OLD STYLE)	202
GOH005**	01-90-2713	HEATER 115 VAC 160 Ω	202
GOL026	01-90-3570	PILOT LAMP	202
GOS001	01-90-0550	NO. 8 STOPPER	200
GOS010*	01-90-3050	MICRO-SWITCH	202
GOS030	(None)	ON/OFF SWITCH	202
GOS031	(None)	HI/LOW SWITCH	202
GOV011	01-90-3433	SUCTION VALVE (INLET)	202
GOV012	01-90-3432	PRESSURE VALVE (OUTLET)	202
GOV013	01-90-8140	AEROVENT RUBBER VALVE	202

* Indicates unit with electro-mechanical control
 ** Indicates unit with electronic control

Model 965A

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOG014	01-90-2330	VALVE GASKET	201
GOH004*	01-90-2713	HEATER 115 VAC 149 Ω (OLD STYLE)	202
GOH005**	01-90-2713	HEATER 115 VAC 160 Ω	202
GOS001	01-90-0550	NO. 8 STOPPER	200
GOS010*	01-90-3050	MICRO-SWITCH	202
GOV011	01-90-3433	SUCTION VALVE (INLET)	202
GOV012	01-90-3432	PRESSURE VALVE (OUTLET)	202
GOV013	01-90-8140	AEROVENT RUBBER VALVE	202

* Indicates unit with electro-mechanical control
 ** Indicates unit with electronic control

Model 6000 (600)

● CAS063	01-90-2304	ON/OFF SWITCH	200
GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOB023	01-90-2311	PC BOARD	202
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOH005**	01-90-2713	HEATER 115 VAC 160 Ω	202
GOL032	01-90-2310	PUMP CYCLE LIGHT	202
GOS001	01-90-0550	NO. 8 STOPPER	200
GOV011	01-90-3433	SUCTION VALVE (INLET)	202
GOV012	01-90-3432	PRESSURE VALVE (OUTLET)	202

** Indicates unit with electronic control

Model 6010 (601)

● CAS063	01-90-2304	ON/OFF SWITCH	200
GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOB023	01-90-2311	PC BOARD	202
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOH005**	01-90-2713	HEATER 115 VAC 160 Ω	202
GOL032	01-90-2310	PUMP CYCLE LIGHT	202
GOS001	01-90-0550	NO. 8 STOPPER	200
GOV011	01-90-3433	SUCTION VALVE (INLET)	202
GOV012	01-90-3432	PRESSURE VALVE (OUTLET)	202

** Indicates unit with electronic control

RPI Part #	OEM Part #	Description	Section 2 Page #
------------	------------	-------------	------------------

GOMCO (continued)

Model 6020

● CAS063	01-90-2304	ON/OFF SWITCH	200
GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOL032	01-90-2310	PUMP CYCLE LIGHT	202
GOS001	01-90-0550	NO. 8 STOPPER	200

Model 6650

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOL032	01-90-2310	PUMP CYCLE LIGHT	202
GOS001	01-90-0550	NO. 8 STOPPER	200
GOV011	01-90-3433	SUCTION VALVE (INLET)	202
GOV012	01-90-3432	PRESSURE VALVE (OUTLET)	202

All Electro-Mechanical Thermotic Pumps with Bimetallic Arm

Prior to PC Board

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOH004*	01-90-2713	HEATER 115 VAC 149 Ω	202
GOS001	01-90-0550	NO. 8 STOPPER	200
GOS010	01-90-3050	MICRO-SWITCH	202

* Indicates unit with electro-mechanical control

All Models with Aerovent Valve

GOA019	01-90-2000	TUBING KIT (DISPOSABLE)	200
GOC020	01-90-2674	600ml RUBBER BOTTLE CAP	200
● GOC044	01-90-3695	DISPOSABLE COLLECTION CANISTER (1100ML)	200
● GOC045	01-90-3712	DISPOSABLE COLLECTION CANISTER (2100ML)	200
GOF024	01-90-3928	BACTERIA FILTER	200
GOF025	01-90-3100	BACTERIA FILTER	200
GOS001	01-90-0550	NO. 8 STOPPER	200
GOV013	01-90-8140	AEROVENT RUBBER VALVE	202

SKLAR

Models: All J.R. Tompkins (100-65)

SKA009	02-6521	COMPRESSOR OUTLET ADAPTOR	300
SKC011	02-4513	FILTER CAP	300
SKF005	04-5001	FELT FILTER	300
SKF008	(None)	RUBBER FOOT	300
SKF012	02-4007	SUCTION AND PRESSURE FILTER ASSEMBLY	300
SKG004	02-4506	FILTER CAP GASKET	300
SKG010	02-5101	VACUUM GAUGE	300
SKN003	02-4512	OPEN FILTER NUT	300
SKR007	2-6519	ROTOR BODY	300
SKS001	02-9031	TOGGLE SWITCH	300
SKS006	02-6502	ROTOR SLIDE	300
SKS013	02-6506	ROTOR HOUSING PLATE SCREW	300
SKS014	02-6505	ROTOR HOUSING SCREW	300
SKS021	04-3501	NEEDLE VALVE REGULATOR SPRING	300
SKW015	02-4510	FILTER WINDOW	300

RPI PART #GOC044
OEM PART #01-90-3695

DISPOSABLE COLLECTION CANISTER (1100 ml)

- 1100 ml
- Clear plastic graduated canister with blue lid
- Hydro-filter shut-off
- Large pour spout

Fits: Stand or mounting bracket ring
Diaphragm Pump Models: 300 Series

Rotary Pump Models: Rotary pumps with 5/8" rotor and 1/20 HP motor; with 3/4" rotor and 1/6 HP motor; and, with 1" rotor and 1/4 or 1/6 HP motor

Thermotic Pump Models: 6000, 6020 and 6650



RPI Part #GOF024
OEM Part #01-90-3928

BACTERIA FILTER

- 3 per package
- 0.3 micron
- 1/4" Barb x 1/4" hose barb
- Disk size 2.100"

Models:

All models incorporating bacteria filter as standard equipment



RPI PART #GOC045
OEM PART #01-90-3712

DISPOSABLE COLLECTION CANISTER (2100 ml)


- 2100 ml
- Clear plastic graduated canister with blue lid
- Hydro-filter shut-off
- Large pour spout; Tandem ports

Fits: Stand or mounting bracket ring

Diaphragm Pump Models: 300 Series and all diaphragm pump models

Rotary Pump Models: Rotary pumps with 5/8" rotor and 1/20 HP motor; with 3/4" rotor and 1/6 HP motor; and, with 1" rotor and 1/4 or 1/6 HP motor

Thermotic Pump Models: 6000, 6020 and 6650



RPI Part #GOF025
OEM Part #01-90-3100

BACTERIA FILTER

- 3 per package
- 0.3 micron
- 1/4" Barb x 1/8" MPT
- Disk size 2.100"

Models: All models incorporating bacteria filter as standard equipment



RPI Part #CAS063
OEM Part #01-90-2304

ON/OFF SWITCH


- SPDT
- 15 A @ 125 VAC; 10 A @ 250 VAC;
- 1/2 HP @ 250 VAC

Fits: Front panel of base

Diaphragm Pump Models: 280, 2801, 3020, 3040 & 6030

Rotary Pump Models: 400, 402, 4010, 4030, 4040, 4050 & 4070

Thermotic Pump Models: 6000, 6010 & 6020
Motor: 6050




New RPI Part #GOA019
Old RPI Part #GOK019
OEM Part #01-90-2000

TUBING KIT (DISPOSABLE)

- Includes (1) each of 15" and 6' long tubing

Models: All models



RPI Part #GOS001
OEM Part #01-90-0550

NO. 8 STOPPER


Fits: Bottle Top
Models: All models with trap bottle



RPI Part #GOC020
OEM Part #01-90-2674

600ml RUBBER BOTTLE CAP

Fits: All models using 600ml bottle top



Also available to fit Diaphragm Pumps


- Disp. Collection Canisters
- Bacteria Filters
- Rubber Bottle Cap
- No. 8 Stopper
- On/Off Switch

See parts above.

RPI PART #GOK043
OEM PART #01-90-2295

DIAPHRAGM PARTS KIT (300 SERIES)

Fits: Pump head
Models: 300 Series




Kit includes: Head gasket, (2) valves, (2) valve retainers, (2) valve screws, and (3) filters

RPI PART #GOG046
OEM PART #01-90-3059

VACUUM GAUGE (300 SERIES PORT)

Fits: Muffler bar
Models: 300 Series




Also available to fit Rotary Pumps

- Disp. Collection Canisters
 - Bacteria Filters
 - Rubber Bottle Cap
 - No. 8 Stopper
 - On/Off Switch
- See Section 2, page 200.

RPI Part #GOG006
OEM Part #01-90-2338

VACUUM GAUGE
• 1.500" diameter dial
• 0-30" Hg. (0-76 mm Hg.)
• Black ABS case
Fits: Muffler bar
Models: All models requiring 1/8" MPT with center back mount



RPI Part #GOG014
OEM Part #01-90-2330

VALVE GASKET
Fits: Overflow Valve Cover (back)
Models: All models with overflow valve



RPI Part #GOS007
OEM Part #01-90-2902

VALVE SPRING
Fits: Overflow Valve Back
Models: All models with overflow valve




RPI PART #GOW037
OEM PART #01-90-3067

5/8" OIL WICK
• 3 per package; White F-1 felt
Models: Fits rotary pumps with 5/8" rotor and 1/20 HP motor



RPI PART #GOS036
OEM PART #01-90-2322


5/8" ROTOR SLIDE
• 3 per package; brass
Models: Fits rotary pumps with 5/8" rotor and 1/20 HP motor



RPI Part #GOK027
OEM Part #01-90-2291

ROTARY UNIT PARTS KIT
• Includes 15" long tubing, 6' long tubing, rotary pump oil (4 fl. oz.), 600 ml rubber bottle cap, (4) felt filter disc, and (4) filter disc
Models: All rotary pumps with overflow valve and 600ml bottle

All parts in Kit also sold separately




RPI Part #GOD018
OEM Part #01-90-3108

FILTER DISC
• 12 per package
• Replaces old style chamois disc
• When replacing this part, it is also a good time to replace the Felt Filter Disc (RPI Part# GOF002)
Fits: Between Valve Spring and Overflow Valve Cover Gasket inside Overflow Valve
Models: All models with overflow valve



RPI Part #GOF002
OEM Part #01-90-3107

FELT FILTER DISC
• 12 per package
Fits: Overflow Valve Cover (front)
Models: All models with overflow valve




RPI Part #GOL022
OEM Part #01-90-3106

ROTARY PUMP OIL
• 4 fl. oz.
• ISO68
• Safety Data Sheet included
For: All Rotary Pumps (dispenses in oil reservoir)




RPI PART #GOW039
OEM PART #01-90-2876

1" OIL WICK
• 3 per package; White F-1 felt
Models: Fits rotary pumps with 1" rotor and 1/4 or 1/6 HP motor




RPI PART #GOS038
OEM PART #01-90-0853

1" ROTOR SLIDE
• 3 per package; brass
Models: Fits rotary pumps with 1" rotor and 1/4 or 1/6 HP motor




RPI Part #GOK017
OEM Part #01-90-2416

3/4" ROTOR SLIDE ASSEMBLY KIT
• Includes 3/4" rotor, and (3) 3/4" rotor slide
Fits: Pump Housing
Models: All 1/6 HP rotary pumps




RPI Part #GOR028
OEM Part #01-90-2417

3/4" ROTOR
Fits: Pump Housing
Models: All 1/6 HP rotary pumps




RPI PART #GOW035
OEM PART #01-90-2875

3/4" OIL WICK
• 3 per package; White F-1 felt
Models: Fits rotary pumps with 3/4" rotor and 1/6 HP motor



RPI Part #GOS029
OEM Part #01-90-2418

3/4" ROTOR SLIDE
• 3 per package; brass
Models: Fits rotary pumps with 3/4" rotor and 1/6 HP motor



Also available to fit Thermotic Pumps

- Disp. Collection Canisters
 - Bacteria Filters
 - Rubber Bottle Cap
 - No. 8 Stopper
 - On/Off Switch
- See Section 2, page 200.*

RPI Part #GOS010
OEM Part #01-90-3050

MICRO-SWITCH

- 125/230/480 VAC
- 15 A
- 1/8 HP @ 125 VAC,
1/4 HP @ 250 VAC

Fits: Thermostat Heater Arm


Models: All electromechanical thermotic pumps with Bimetallic arm prior to PC board



RPI Part #GOV012
OEM Part #01-90-3432

PRESSURE VALVE (OUTLET)


Fits: Cylinder Assembly
Models: 764, 765A, 965A, 6000 (600), 6010 (601) & 6650



RPI Part #GOV011
OEM Part #01-90-3433

SUCTION VALVE (INLET)

Fits: Cylinder Assembly
Models: 764, 765A, 965A, 6000 (600), 6010 (601) & 6650



RPI Part #GOV013
OEM Part #01-90-8140

AEROVENT RUBBER VALVE


Fits: Aerovent Assembly
Models: All models with Aerovent valve



RPI Part #GOS031
OEM Part #(None)

HI/LOW SWITCH


- SPST
 - 3 A @ 250 VDC; 6 A @ 125 VAC
- Fits:** PC Board (RPI Part #GOB016) only
Models: 764 & 765A



RPI Part #GOS030
OEM Part #(None)

ON/OFF SWITCH


- DPST
 - 3 A @ 250 VDC; 6 A @ 125 VAC
- Fits:** PC Board (RPI Part #GOB016) only
Models: 764 & 765A



RPI Part #GOL026
OEM Part #01-90-3570

PILOT LAMP


- Lens: Red
 - Type: Neon
 - 125 VAC, 1/3 W
- Fits:** Fror
Models:



RPI Part #GOL032
OEM Part #01-90-2310

PUMP CYCLE LIGHT


- Lens: Amber
 - Type: Neon
 - 125 VAC
- Fits:** Front Panel
Models: 270, 6000 (600), 6010 (610), 6020 & 6650



RPI Part #GOB023
OEM Part #01-90-2311

PC BOARD


- New style
 - 115 VAC
- Fits:** Case
Models: 6000 (600) & 6010 (601)



RPI Part #GOB016
OEM Part #01-90-3045

PC BOARD


- Old style
 - 115 VAC
- Fits:** Front Panel
Models: 764 & 765A



RPI Part #GOH005
OEM Part #01-90-2713

HEATER 115 VAC (160 Ω)


Fits: Cylinder Assembly
Models: 765A, 764, 6000 (600) & 6010 (601) all models using PC board



RPI Part #GOH004
OEM Part #01-90-2713

HEATER 115 VAC (149 Ω)

Fits: Cylinder Assembly
Models: All electro-mechanical thermotic pumps with Bimetallic arm prior to PC board



Aspirators & Pumps


Models: All J.R. Tompkins (100-65)

RPI Part #SKS001
OEM Part #02-9031

TOGGLE SWITCH

- DPST
- 6 A @ 125 VAC
- Includes push-on terminals, and wire nuts
- Mounts in base of machine


Model: J.R. Tompkins (100-65)



RPI Part #SKF012
OEM Part #02-4007

SUCTION AND PRESSURE FILTER ASSEMBLY


- Mounts off compressor side
- Model:** J.R. Tompkins (100-65)



RPI Part #SKA009
OEM Part #02-6521

COMPRESSOR OUTLET ADAPTOR

- 1/8" MPT x 7/16-27
- Fits:** Compressor Section
Model: J.R. Tompkins (100-65)



RPI Part #SKN003
OEM Part #02-4512

OPEN FILTER NUT

- 5/16-24 thread
 - Chrome plated
- Fits:** Suction and Pressure Filter Assembly

Model: J.R. Tompkins (100-65)




RPI Part #SKW015
OEM Part #02-4510

FILTER WINDOW

- Assembly requires two
- Fits:** Suction and Pressure Filter Assembly

Model: J.R. Tompkins (100-65)



RPI Part #SKF008
OEM Part #None

RUBBER FOOT

- 4 per package
 - Mounts to bottom of base
- Model:** J.R. Tompkins (100-65)



RPI Part #SKG004
OEM Part #02-4506

FILTER CAP GASKET


- 2 per package
 - Assembly requires two
- Fits:** Filter Cap
Model: J.R. Tompkins (100-65)



RPI Part #SKG010
OEM Part #02-5101

VACUUM GAUGE


- 1.500" diameter dial
 - 1/8" MPT, bottom mount
 - Reads: 0-30 in. Hg.
 - Mounts off suction and pressure filter assembly
- Model:** J.R. Tompkins (100-65)



RPI Part #SKS021
OEM Part #04-3501

NEEDLE VALVE REGULATOR SPRING

- **Fits:** Regulator
- Model:** J.R. Tompkins (100-65)




RPI Part #SKF005
OEM Part #04-5001

FELT FILTER

- 12 per package
- Assembly requires four

Fits: Window portion of Suction and Pressure Filter Assembly
Model: J.R. Tompkins (100-65)



RPI Part #SKR007
OEM Part #2-6519

ROTOR BODY

- **Fits:** Compressor Section
- Model:** J.R. Tompkins (100-65)



RPI Part #SKS013
OEM Part #02-6506

ROTOR HOUSING PLATE SCREW

- 3 per package
 - Slotted fillister head
 - #8-32 x 1.230" long
 - Stainless steel
- Fits:** Rotor Housing Plate
Model: J.R. Tompkins (100-65)



RPI Part #SKC011
OEM Part #02-4513

FILTER CAP


- Assembly requires two
- Fits:** Suction and Pressure Filter Assembly
Model: J.R. Tompkins (100-65)



RPI Part #SKS006
OEM Part #02-6502

ROTOR SLIDE

- 3 per package
- Fits:** Rotor Body
Model: J.R. Tompkins (100-65)



RPI Part #SKS014
OEM Part #02-6505

ROTOR HOUSING SCREW

- 3 per package
 - Slotted fillister head
 - #8-32 x .975" long
 - Stainless steel
- Fits:** Rotor Housing
Model: J.R. Tompkins (100-65)

