

RPI Part #	Industry Part #	Description	Section 8 Page #
LAMPS & BULBS			
AMB028	NE51	LIGHT BULB	201
AIB041	GE345	BULB, NEON	201
LMP001	GE252	LAMP	200
LMP002	GE778	LAMP	200
LMP003	GE787	LAMP	200
LMP004	GE1460	LAMP	200
LMP005	15S11/102	LAMP	200
LMP006	FCR7023	LAMP	200
LMP007	ESA/FHD	LAMP	200
LMP008	ESB7388	LAMP	200
LMP009	GE1493	LAMP	200
LMP010	15S11/CL13	LAMP	200
LMP011	BRL7027	LAMP	200
LMP012	13347W	LAMP	200
LMP013	GE1634	LAMP	200
LMP014	EXN	LAMP	200
LMP015	25T10	LAMP	200
LMP101	GE1004	SPECTROPHOTOMETER LAMP	200
LMP102	D2A	SPECTROPHOTOMETER LAMP	200
LMP103	GE82	SPECTROPHOTOMETER LAMP	201
LMP104	BWD	SPECTROPHOTOMETER LAMP	201
LMP600	FCS	LAMP	201
LMP601	H-3	LAMP	201
LMP602	125T10/P	LAMP	201
LMP603	EJL	LAMP	201
LMP604	ENZ	LAMP	201
LMP605	ETJ	LAMP	201
LMP606	ELZ	LAMP	201
LMP607	FDS/DZE	LAMP	201
LMP608	75T3Q/CL	LAMP	201
PCL606	150T4Q/CL	LAMP	201
PCT685	N/A	SOCKET RETRACTOR TOOL	201

● NEW PART

Call Us Toll-Free Monday–Friday, 8:00 AM–4:30 PM (Pacific Time), (800) 221-9723

SECTION

8

PAGE

100

UPDATED

6/98

RPI Part #LMP001
Industry Part #GE252

LAMP

- 2.5V, .35A
- Base: Single contact, midget flanged



RPI Part #LMP007
Industry Part #ESA/FHD

LAMP


- 6V, 10W
- Base: 2-pin



RPI Part #LMP013
Industry Part #GE1634

LAMP


- 20V, 1.0A
- Base: Double contact, prefocus
- Halogen



RPI Part #LMP002
Industry Part #GE778

LAMP

- 6.0V, 20W
- Base: 2-pin
- Halogen



RPI Part #LMP008
Industry Part #ESB 7388

LAMP


- 6V, 20W
- Base: 2-pin



RPI Part #LMP014
Industry Part #EXN

LAMP

- 12V, 50W
- Base: 2-pin
- Incandescent



RPI Part #LMP003
Industry Part #GE787

LAMP


- 6.0V, 10W
- Base: 2-pin
- Halogen



RPI Part #LMP009
Industry Part #GE1493

LAMP


- 6.5V, 2.75A
- Base: Double contact, bayonet
- Halogen



New RPI Part #LMP015
Old RPI Part #HAL007
Industry Part #25T10

LAMP


- 120V, 25W
- Base: Medium screw, prefocus



RPI Part #LMP004
Industry Part #GE1460

LAMP

- 6.5V, 2.75A
- Base: Double contact, prefocus
- Halogen



RPI Part #LMP010
Industry Part #15S11/CL13

LAMP


- 115V, 15W
- Base: Candelabra



RPI Part #LMP101
Industry Part #GE1004

SPECTROPHOTOMETER LAMP


- 12.8V, .94A
- Base: Double contact, bayonet with mounting bracket
- Halogen
- 18 GA black and red wire with connectors



RPI Part #LMP005
Industry Part #15S11/102

LAMP

- 120V, 15W
- Base: Medium screw, prefocus



RPI Part #LMP011
Industry Part #BRL 7027

LAMP


- 12.V, 50W
- Base: Double contact, miniature



RPI Part #LMP012
Industry Part #13347W

LAMP


- 6V, 15W
- Base: Double contact, bayonet candelabra
- Halogen



RPI Part #LMP102
Industry Part #D2A

SPECTROPHOTOMETER LAMP

- 10.8V, 30W
- Base: Mini, 2-pin



Lamps & Bulbs

RPI Part #LMP103
Industry Part #GE82

SPECTROPHOTOMETER LAMP


- 6.5V, 1.02A
- Base: Double contact, bayonet
- Halogen



RPI Part #LMP104
Industry Part #BWD

SPECTROPHOTOMETER LAMP


- 8.5V, 4.2A
- Base: Double contact, bayonet with mounting bracket
- 18 GA black and red wire 8" with connectors



RPI Part #LMP600
Industry Part #FCS

LAMP


- 24V, 150W
- Base: 2-pin



RPI Part #LMP601
Industry Part #H-3

LAMP


- 12V, 55W
- Base: Wire
- Halogen



RPI Part #LMP602
Industry Part #125T10P

LAMP


- 120V, 125W
- Base: Medium, prefocus



RPI Part #LMP603
Industry Part #EJL

LAMP


- 24V, 200W
- Base: 2-pin



RPI Part #LMP604
Industry Part #ENZ

LAMP


- 30V, 50W
- Base: 2-pin



RPI Part #LMP605
Industry Part #ETJ

LAMP


- 120V, 250W
- Base: 2-pin, special



RPI Part #LMP606
Industry Part #ELZ

LAMP


- 21V, 150W
- Base: 2-pin



RPI Part #LMP607
Industry Part #FDS/DZE

LAMP


- 24V, 150W
- Base: 2-pin



RPI Part #LMP608
Industry Part #75T3Q/CL

LAMP


- 28V, 75W
- Base: Double contact, bayonet
- Halogen



RPI Part #AMB028
Industry Part #NE51

LIGHT BULB

- Standard NE51 bulb



RPI Part #AIB041
Industry Part #GE345

BULB, NEON


- T3, 1/4 E10 base
- 2 watts, 24/30 VDC



RPI Part #PCL606
Industry Part #150T4Q/CL

LAMP


- 24V, 150W
- Base: Double ended
- Halogen



RPI Part #PCT685

SOCKET RETRACTOR TOOL

- For installation of lamp (RPI Part #PCL606)



Use the handy
fax order form
located at the front
of the catalog.