



**replacement parts industries, inc.**

*"The Alternate Source™"*

CALL 800-221-9723 or 818-882-8611 • FAX 818-882-7028  
E-MAIL: techsupport@rplparts.com • WEBSITE: www.rplparts.com

## **GOS029, GOS036 & GOS038**

### **5/8", 3/4" & 1" ROTOR SLIDES**

#### **CUSTOMIZING INSTRUCTIONS**

To allow for universal fit, the slides have been manufactured slightly over-sized. To ensure optimal fit, the flat side of the slides can be polished on a surface plate using Very Fine to Ultra Fine wet/dry sand paper. Care should be taken not to over sand. Because of close tolerance clearances of the rotor and slides, pump performance may be compromised with excessive sanding of slides. All traces of loose metal as a result of sanding must be removed by thoroughly washing slides in Kerosene. NOTE: Refer to service manual for proper maintenance of pump.

***CAUTION: Do not sand radial sides of slides.***

