



replacement parts industries, inc.

"The Alternate Source"™

CALL 800-221-9723 or 818-882-8611 • FAX 818-882-7028

E-MAIL: techsupport@rplparts.com • WEBSITE: www.rplparts.com

**PCG038 & PCG040
GAUGES
IMPORTANT NOTICE**

The Gauge (RPI Part #PCG040) is now a center-back-mount gauge and can be used in place of the PCG038.

The manufacturer of both gauges no longer makes a lower-back-mount gauge so as a replacement, RPI offers the PCG040 as a center-back-mount that fits the **Pelton & Crane Magna-clave®**, **OCM**, **OCR**, and **OCR+**.