



"The Alternate Source™"

**replacement parts industries, inc.**

CALL 800-221-9723 or 818-882-8611 • FAX 818-882-7028

E-MAIL: techsupport@rpiparts.com • WEBSITE: www.rpiparts.com

**RPI Part #SDK003  
OEM Part #100798-02  
PM-2 KIT (ADIXEN/PFEIFFER PUMP)  
FITS STERRAD  
Model: NX**

Includes:

- (1) UV Lamp Assembly (RPI Part #SDA006)
- (5) Thermistor Assembly (RPI Part #SDA007)
- (1) Catalytic Converter (RPI Part #SDC005)
- (1) Air Filter (Base) (RPI Part #SDF001)
- (1) Oil Mist Filter (Red) (RPI Part #SDF002)
- (1) Conical Screen Filter (RPI Part #SDF017)
- (1 pc) Air Filter (25mm) (RPI Part #SDF025)
- (1 pkg of 2) Handle (Shelf Extension) (RPI Part #SDH011)
- (1) Feedthru Screw Washer Kit (RPI Part #SDK012)
- (1) Optical Window Kit (RPI Part #SDK014)
- (1 pkg of 2) Shelf Support (Upper) (RPI Part #SDS008)
- (1 pkg of 2) Shelf Support (Lower) (RPI Part #SDS009)
- (1 pkg of 6) Spacer (Electrode) (RPI Part #SDS013)
- (1 pkg of 2) Sleeve Gasket (RPI Part #SDS015)
- (1) Chamber/Door Seal (RPI Part #SDS018)
- (1) Delivery Valve (RPI Part #SDV010)
- (1) Sterile Air (HEPA) Filter (RPI Part #SSF011)
- (1) Bottle (1 liter/32 fl. oz.) (RPI Part #RPB939)
- (1) Thermal Grease (Thermalcote™) (RPI Part #RPG938)  
*SDS sheet available on RPI website, www.rpiparts.com*
- (1) Vacuum Pump Oil (Adixen/Pfeiffer Pump) (RPI Part #RPL926)  
*SDS sheet available on RPI website, www.rpiparts.com*
- (1 pc) Absorbent Pad (RPI Part #RPP958)
- (18" lg.) Tubing (3/8" ID x 1/2" OD Clear) (RPI Part #RPT646)

All parts also sold separately.



- Personal protective equipment, including eye protection and gloves, should be worn.
- Properly dispose of all collected oil and any contaminated items as required by law, and according to local, state and federal regulations.