



"The Alternate Source"™

replacement parts industries, inc.

CALL 800-221-9723 or 818-882-8611 • FAX 818-882-7028

E-MAIL: techsupport@rpiparts.com • WEBSITE: www.rpiparts.com

RPI Part #SDK037
OEM Part # (No OEM Part # Available)
CHAMBER PLASTICS KIT
FITS STERRAD
Model: 100NX

Includes:

- (1 pkg of 4) Shelf Handle (Extrusion) (RPI Part #SDH065)
- (1 pkg of 2) Shelf Support (Lower) (RPI Part #SDS038)
- (1 pkg of 2) Shelf Support (Upper) (RPI Part #SDS039)
- (1 pkg of 4) Electrode Spacer (Top) (RPI Part #SDS040)
- (1 pkg of 2) Electrode Spacer (Bottom) (RPI Part #SDS055)
- (1 pkg of 6) Electrode Spacer (Side) (RPI Part #SDS056)

All parts also sold separately.



WARNING

- Personal protective equipment, including eye protection and gloves, should be worn.