



"The Alternate Source"

replacement parts industries, inc.

CALL 800-221-9723 or 818-882-8611 • FAX 818-882-7028

E-MAIL: techsupport@rpiparts.com • WEBSITE: www.rpiparts.com

RPI Part #SDK058
OEM Part #05-53441-1-001A
PM-1 KIT
FITS STERRAD
Model: 100NX

Includes:

- (1) Oil Mist Filter Kit (RPI Part #SDF061)
- (1) Air Filter (Rear Panel) (RPI Part #SDF066)
- (2) Inlet Valve (RPI Part #SDV047)
- (1 pc) Disposal Bag (RPI Part #RPB838)
- (3) Bottle (1 liter/32 fl. oz.) (RPI Part #RPB939)
- (1 pkg of 3) Absorbent Pad (RPI Part #RPP958)
- (18" lg.) Tubing (3/8" ID x 1/2" OD Clear) (RPI Part #RPT646)
- (3) Vacuum Pump Oil (Leybold Pump) (RPI Part #RXL034)

SDS sheet available on RPI website, www.rpiparts.com

All parts also sold separately.



WARNING

- Personal protective equipment, including eye protection and gloves, should be worn.
- Properly dispose of all collected oil as required by law, and according to local, state and federal regulations.