



"The Alternate Source™"

replacement parts industries, inc.

CALL 800-221-9723 or 818-882-8611 • FAX 818-882-7028

E-MAIL: techsupport@rpiparts.com • WEBSITE: www.rpiparts.com

RPI Part #SDK059 OEM Part #05-53441-1-002A

PM-2 KIT

FITS STERRAD

Model: 100NX

Includes:

- (1) UV Lamp Assembly (RPI Part #SDA006)
- (7) Thermistor Assembly (RPI Part #SDA007)
- (1) Tubing Assembly (H₂O₂ Delivery) (RPI Part #SDA068)
- (1) Catalytic Converter (RPI Part #SDC005)
- (1 pc) Air Filter (25mm) (RPI Part #SDF025)
- (1) Oil Mist Filter Kit (RPI Part #SDF061)
- (1) Air Filter (Rear Panel) (RPI Part #SDF066)
- (1 pkg of 4) Shelf Handle (Extrusion) (RPI Part #SDH065)
- (1) Optical Window Kit (RPI Part #SDK014)
- (1 pkg of 2) Needle H₂O₂ Transfer (RPI Part #SDN062)
- (1 pkg of 2) Shelf Support (Lower) (RPI Part #SDS038)
- (1 pkg of 2) Shelf Support (Upper) (RPI Part #SDS039)
- (1 pkg of 4) Electrode Spacer (Top) (RPI Part #SDS040)
- (1 pkg of 2) Electrode Spacer (Bottom) (RPI Part #SDS055)
- (1 pkg of 6) Electrode Spacer (Side) (RPI Part #SDS056)
- (1) Chamber/Door Seal (RPI Part #SDS063)
- (1) Delivery Valve (RPI Part #SDV010)
- (2) Inlet Valve (RPI Part #SDV047)
- (1) Sterile Air (HEPA) Filter (RPI Part #SSF011)
- (1 pc) Disposal Bag (RPI Part #RPB838)
- (3) Bottle (1 liter/32 fl. oz.) (RPI Part #RPB939)
- (1) Thermal Grease (Thermalcote™) (RPI Part #RPG938)

SDS sheet available on RPI website, www.rpiparts.com

- (1 pkg of 3) Absorbent Pad (RPI Part #RPP958)
- (18" lg.) Tubing (3/8" ID x 1/2" OD Clear) (RPI Part #RPT646)
- (3) Vacuum Pump Oil (Leybold Pump) (RPI Part #RXL034)

SDS sheet available on RPI website, www.rpiparts.com

All parts also sold separately.



WARNING

- Personal protective equipment, including eye protection and gloves, should be worn.
- Properly dispose of all collected oil and any contaminated items as required by law, and according to local, state and federal regulations.