



"The Alternate Source"SM

replacement parts industries, inc.

CALL 800-221-9723 or 818-882-8611 • FAX 818-882-7028
E-MAIL: techsupport@rpiparts.com • WEBSITE: www.rpiparts.com

RPI Part #SDK060
OEM Part # 05-06899-2-200
PM-2 KIT
FITS STERRAD
Model: 100S

Includes:

- (1) Catalytic Converter (RPI Part #SDC005)
- (1) Extender (Electrode) (RPI Part #SDE036)
- (1) Oil Mist Filter (Red) (RPI Part #SDF002)
- (1) Mist Filter (RPI Part #SDF022)
- (1 pkg of 2) Handle (Shelf Extension) (RPI Part #SDH011)
- (1 pkg of 2) Vaporizer Plate (RPI Part #SDP064)
- (1 pkg of 2) Guide Rails (RPI Parts #SDR035)
- (1 pkg of 2) Shelf Support (Lower) (RPI Part #SDS033)
- (1 pkg of 2) Shelf Support (Upper) (RPI Part #SDS034)
- (1) Chamber/Door Seal (RPI Part #SDS057)
- (1) Sterile Air (HEPA) Filter (RPI Part #SSF011)
- (1 pc) Disposal Bag (RPI Part #RPB838)
- (1) Bottle (1 liter/32 fl. oz.) (RPI Part #RPB939)
- (1 pc) Absorbent Pad (RPI Part #RPP958)
- (18" lg.) Tubing (3/8" ID x 1/2" OD Clear) (RPI Part #RPT646)
- (2) Vacuum Pump Oil (Leybold Pump) (RPI Part #RXL034)

SDS sheet available on RPI website, www.rpiparts.com

All parts also sold separately.



WARNING

- **Personal protective equipment, including eye protection and gloves, should be worn.**
- **Properly dispose of all collected oil and any contaminated items as required by law, and according to local, state and federal regulations.**