



replacement parts industries, inc.

A Healthcare Components Group Company

CALL 800-221-9723 or 818-882-8611 • FAX 818-882-7028

E-MAIL: techsupport@rpiparts.com • WEBSITE: www.rpiparts.com

**RPI Part #SSP058
OEM Part #P764333-354
PM PACK (GENERAL)**

FITS STERIS®

Models: V-Pro 1, V-Pro 1+ & V-Pro Max

Includes:

- (1) Solenoid Valve Repair Kit For 1/4" Port Valves (RPI Part #AMK175)
- (1) Air Intake Filter (RPI Part #CMF282)
- (1) Filter Element (RPI Part #SSE061)
- (1) HEPA Filter (RPI Part #SSF059)
- (4) Exhaust Filter (RPI Part #SSF060)
- (2) Isolation Valve (3-Way) (RPI Part #SSV062)
- (1) Battery (Lithium) (RPI Part #RXB076)
SDS sheet available on RPI website, www.rpiparts.com
- (1 pc) Male Nut (1/8" OD Tubing) (RPI Part #RXF073)
- (1 pc) Ferrule (1/8" OD Tubing) (RPI Part #RXF074)
- (1) Tee (1/8" OD Tubing) (RPI Part #RXF075)

All parts also sold separately.