



replacement parts industries, inc.

A Healthcare Components Group Company

CALL 800-221-9723 or 818-882-8611 • FAX 818-882-7028

E-MAIL: techsupport@rpiparts.com • WEBSITE: www.rpiparts.com

**RPI Part #SSP065
OEM Part #P764334-481
PM PACK (VACUUM)**

FITS STERIS®

Models: V-Pro 1, V-Pro 1+ & V-Pro Max

Includes:

- (1) Air Intake Filter (Metal Style) (RPI Part #CMF002)
- (1) Filter Element (RPI Part #SSE066)
- (1) Exhaust Filter (RPI Part #SSF060)
- (2) Cartridge Filter (RPI Part #SSF067)
- (1) Drain Plug Gasket (RPI Part #SSG068)
- (1 pc) Disposal Bag (RPI Part #RPB838)
- (5 pcs) Absorbent Pad (RPI Part #RPP958)
- (3) Vacuum Pump Oil (Leybold Pump) (RPI Part #RXL034)

SDS sheet available on RPI website, www.rpiparts.com

All parts also sold separately.



WARNING

- Personal protective equipment, including eye protection and gloves, should be worn.
- Properly dispose of all collected oil as required by law, and according to local, state and federal regulations.