



replacement parts industries, inc.

A Healthcare Components Group Company

CALL 800-221-9723 or 818-882-8611 • FAX 818-882-7028

E-MAIL: techsupport@rpiparts.com • WEBSITE: www.rpiparts.com

**RPI Part #SSP069
OEM Part #P764335-557
PM PACK (VACUUM)
FITS STERIS®
Models: V-Pro 60**

Includes:

- (1) Diaphragm (Valve Umbrella) (RPI Part #SSD071)
- (1) Filter Element (RPI Part #SSE070)
- (1) Cartridge Filter (RPI Part #SSF067)
- (1 pkg) Sealing Washer (RPI Part #SSW072)
- (1 pc) Disposal Bag (RPI Part #RPB838)
- (1) Vacuum Pump Oil (Leybold Pump) (RPI Part #RXLO34)
SDS sheet available on RPI website, www.rpiparts.com
- (3 ft) Tubing (5/16" ID x 9/16" OD Natural) (RPI Part #RXT077)

All parts also sold separately.



WARNING

- Personal protective equipment, including eye protection and gloves, should be worn.
- Properly dispose of all collected oil as required by law, and according to local, state and federal regulations.