



replacement parts industries, inc.

CALL 800-221-9723 or 818-882-8611 • FAX 818-882-7028
E-MAIL: techsupport@rpiparts.com • WEBSITE: www.rpiparts.com

**RPI Part #SDK004
OEM Part #100797-02
PM-1 KIT (ADIXEN/PFEIFFER PUMP)
FITS STERRAD
Model: NX**

Includes:

- (1) Oil Mist Filter (Red) (RPI Part #SDF002)
- (1) Bottle (1 liter/32 fl. oz.) (RPI Part #RPB939)
- (1) Vacuum Pump Oil (Adixen/Pfeiffer Pump) (RPI Part #RPL926)
SDS sheet available on RPI website, www.rpiparts.com
- (1 pc) Absorbent Pad (RPI Part #RPP958)
- (18" lg.) Tubing (3/8" ID x 1/2" OD Clear) (RPI Part #RPT646)

All parts also sold separately.



WARNING

- Personal protective equipment, including eye protection and gloves, should be worn.
- Properly dispose of all collected oil as required by law, and according to local, state and federal regulations.