



## New Parts and Upgrades to fit A-DEC UTILITY MODULE and REGULATORS

Now Available ... RPI's Clear Valve Body!

Easily locate malfunctioning valve diaphragms and pistons at-a-glance. The new clear valve body is see-through so that you can actually see what's happening inside the valve and troubleshoot from the outside. There's no need to dismantle the entire unit to find the troublesome valve!

### WHITE VALVE BODY

RPI Part #ADB197  
OEM Part #24-0135-00

**REPLACED BY THE  
CLEAR VALVE BODY -  
ADB202!**



### BLACK VALVE BODY

RPI Part #ADB198  
OEM Part #24-0355-00

**REPLACED BY THE  
CLEAR VALVE BODY -  
ADB202!**



**Now Available ...  
THE NEW RPI CLEAR VALVE BODY.  
Replacement for the White and Black Bodies!**

### CLEAR VALVE BODY

RPI Part #ADB202   
OEM Part #24-0135-00 & 24-0355-00

• Clear Plastic


Fits: Air Regulator Valve "White" Bodies

OEM Part #24-0128-00 & #24-0135-00 &


Water Regulator Valve "Black" Bodies

OEM Part #24-0129-00 & #24-0138-00

Models: Cascade, Decade, Excellence & Mini-Trol



Now It's Clear ... So you can troubleshoot  
without taking apart the unit!



Save time and troubleshoot from the outside:

- Easily check for air bubbles
- Diagnose diaphragm malfunction
- Evaluate piston performance
- Examine o-ring and spring integrity
- Identify debris contamination from Air Compressor

RPI Part #ADK153

OEM Part # (No OEM Part # Available)



### WHITE VALVE BODY REPAIR KIT

• This kit has been upgraded to include the Clear Valve Body (RPI Part #ADB202)

Fits: Utility Module OEM Part #24-0352-00

and Air Regulator Valve "White" Bodies

OEM Part #24-0128-00 & #24-0135-00

Models: Cascade, Decade, Excellence & Mini-Trol

RPI Part #ADK173

OEM Part # (No OEM Part # Available)



### BLACK VALVE BODY REPAIR KIT

• This kit has been upgraded to include the Clear Valve Body (RPI Part #ADB202)

Fits: Utility Module OEM Part #24-0352-00

and Water Regulator Valve "Black" Bodies

OEM Part #24-0129-00 & #24-0138-00

Models: Cascade, Decade, Excellence & Mini-Trol

RPI Part #ADK201

OEM Part #24-0228-00/  
24-0234-00/030-019-00



### FILTER KIT

• New round design for enhanced performance


• The RPI Filter Housing is clear so you can see when the filter needs replacement without having to remove the housing

Fits: Utility Module OEM Part #24-0352-00 and Filter Regulators OEM Part #'s 24-0168-00, 24-0170-00 & 24-0175-00

Models: Cascade, Decade, Excellence & Mini-Trol

### CLEAR VALVE BODY!


RPI's valve body is now clear for easy troubleshooting!



All parts in Kit also sold separately

### CLEAR VALVE BODY!

RPI's valve body is now clear for easy troubleshooting!



All parts in Kit also sold separately

### CLEAR HOUSING!

The RPI filter housing is clear so you can see when the filter needs replacement without having to remove the housing. It's the RPI Advantage.

