



New Parts and Upgrades to fit A-DEC UTILITY MODULE and REGULATORS

Now Available ... RPI's Clear Valve Body!

Easily locate malfunctioning valve diaphragms and pistons at-a-glance. The new clear valve body is see-through so that you can actually see what's happening inside the valve and troubleshoot from the outside. There's no need to dismantle the entire unit to find the troublesome valve!

WHITE VALVE BODY

RPI Part #ADB197
OEM Part #24-0135-00

**REPLACED BY THE
CLEAR VALVE BODY -
ADB202!**



BLACK VALVE BODY

RPI Part #ADB198
OEM Part #24-0355-00

**REPLACED BY THE
CLEAR VALVE BODY -
ADB202!**



**Now Available ...
THE NEW RPI CLEAR VALVE BODY.
Replacement for the White and Black Bodies!**

CLEAR VALVE BODY

RPI Part #ADB202
OEM Part #24-0135-00 & 24-0355-00

• Clear Plastic


Fits: Air Regulator Valve "White" Bodies

OEM Part #24-0128-00 & #24-0135-00 &


Water Regulator Valve "Black" Bodies

OEM Part #24-0129-00 & #24-0138-00

Models: Cascade, Decade, Excellence & Mini-Trol



**Now It's Clear ... So you can troubleshoot
without taking apart the unit!**



Save time and troubleshoot from the outside:

- Easily check for air bubbles
- Diagnose diaphragm malfunction
- Evaluate piston performance
- Examine o-ring and spring integrity
- Identify debris contamination from Air Compressor

RPI Part #ADK153

OEM Part # (No OEM Part # Available)



WHITE VALVE BODY REPAIR KIT

• This kit has been upgraded to include the Clear Valve Body (RPI Part #ADB202)

Fits: Utility Module OEM Part #24-0352-00

and Air Regulator Valve "White" Bodies

OEM Part #24-0128-00 & #24-0135-00

Models: Cascade, Decade, Excellence & Mini-Trol

RPI Part #ADK173

OEM Part # (No OEM Part # Available)



BLACK VALVE BODY REPAIR KIT

• This kit has been upgraded to include the Clear Valve Body (RPI Part #ADB202)

Fits: Utility Module OEM Part #24-0352-00

and Water Regulator Valve "Black" Bodies

OEM Part #24-0129-00 & #24-0138-00

Models: Cascade, Decade, Excellence & Mini-Trol

RPI Part #ADK201

OEM Part #24-0228-00/
24-0234-00/030-019-00



FILTER KIT

• New round design for enhanced performance


• The RPI Filter Housing is clear so you can see when the filter needs replacement without having to remove the housing

Fits: Utility Module OEM Part #24-0352-00 and Filter Regulators OEM Part #'s 24-0168-00, 24-0170-00 & 24-0175-00

Models: Cascade, Decade, Excellence & Mini-Trol

CLEAR VALVE BODY!


RPI's valve body is now clear for easy troubleshooting!



All parts in Kit also sold separately

CLEAR VALVE BODY!

RPI's valve body is now clear for easy troubleshooting!



All parts in Kit also sold separately

CLEAR HOUSING!

The RPI filter housing is clear so you can see when the filter needs replacement without having to remove the housing. It's the RPI Advantage.

