








## LEADWIRE SETS FOR VETERINARY USE

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END	NOTES
LVM701	3 LEADWIRE SET - 24" DIN/ALLIGATOR	3	~24"	Din/Unshielded .060" Safety Socket	Alligator Clip 	<ul style="list-style-type: none"> <li>• Reusable</li> <li>• Recommended for veterinary use only</li> </ul> <b>Also available:</b> <ul style="list-style-type: none"> <li>• 3 Lead Organizer (RPI Part #LAR004)</li> </ul>
LVM702	3 LEADWIRE SET - 24" DIN/ALLIGATOR	3	~24"	Din/Unshielded .060" Safety Socket	Alligator Clip 	<ul style="list-style-type: none"> <li>• Reusable</li> <li>• Recommended for veterinary use only</li> </ul> <b>Also available:</b> <ul style="list-style-type: none"> <li>• 3 Lead Organizer (RPI Part #LAR004)</li> </ul>
LVM703	5 LEADWIRE SET - 24" DIN/ALLIGATOR	5	~24"	Din/Unshielded .060" Safety Socket	Alligator Clip 	<ul style="list-style-type: none"> <li>• Reusable</li> <li>• Recommended for veterinary use only</li> </ul> <b>Also available:</b> <ul style="list-style-type: none"> <li>• 5 Lead Organizer (RPI Part #LAR005)</li> </ul>
LVM704	7 LEADWIRE SET - 24" DIN/ALLIGATOR	7	~24"	Din/Unshielded .060" Safety Socket	Alligator Clip 	<ul style="list-style-type: none"> <li>• Reusable</li> <li>• Recommended for veterinary use only</li> </ul> <b>Also available:</b> <ul style="list-style-type: none"> <li>• 5 Lead Organizer (RPI Part #LAR005)</li> <li>• 3 Lead Organizer (RPI Part #LAR004)</li> </ul>
LVM705	3 LEADWIRE SET - 30" DUAL/ALLIGATOR	3	~30"	Dual/Shielded	Alligator Clip 	<ul style="list-style-type: none"> <li>• Reusable</li> <li>• Recommended for veterinary use only</li> </ul> <b>Also available:</b> <ul style="list-style-type: none"> <li>• 3 Lead Shielded Accessory Kit (RPI Part #LAK006)</li> </ul>
LVM706	3 LEADWIRE SET - 30" DUAL/ALLIGATOR	3	~30"	Dual/Shielded	Alligator Clip 	<ul style="list-style-type: none"> <li>• Reusable</li> <li>• Recommended for veterinary use only</li> </ul> <b>Also available:</b> <ul style="list-style-type: none"> <li>• 3 Lead Shielded Accessory Kit (RPI Part #LAK006)</li> </ul>
LVM707	5 LEADWIRE SET - 30" DUAL/ALLIGATOR	5	~30"	Dual/Shielded	Alligator Clip 	<ul style="list-style-type: none"> <li>• Reusable</li> <li>• Recommended for veterinary use only</li> </ul> <b>Also available:</b> <ul style="list-style-type: none"> <li>• 5 Lead Shielded Accessory Kit (RPI Part #LAK007)</li> </ul>