




A-dec 300 & 500 Series Handpiece Tubing Hook-Ups

Silicone Handpiece Tubing Connections



Silicone Handpiece Tubing Colors and Functions


Color	Function
Clear	Drive Air
Red	Exhaust
Blue	Air Coolant
Green	Water Coolant

Handpiece Tubing Types


Handpiece Tubing with Borden Terminal, 3-Hole



Handpiece Tubing with Midwest Terminal, 4-Hole



Handpiece Tubing, Fiber Optic with Bulb, 4-Hole



Handpiece Tubing, Fiber Optic, 6-pin

