



... Proud to be a **Healthcare Components Group** company ...



**HEALTHCARE
Components Group**

A-DEC


New Parts to fit

A-dec Dental Delivery Units

300 & 500 Series Control Blocks!

RPI Part #ADS279

OEM Part #38.1825.00 (Version 2)



STEM (AIR COOLANT)

- Version 2 design
- Material: Stainless steel
- Double flat stem
- Includes: (1pc) O-Ring (RPI Part #RPO484)
- Also available: Knob (Flow Adj.) (RPI Part #ADK281), High Temp Lubricant (RPI Part #RPL090), and Version 1 Stem (Air Coolant) (RPI Part #ADS234)

Fits: Handpiece Control Block Version 2 (3 and 4 Position) Models: 332, 333, 334, 335, 336 & 541

Models: Stems, Handpiece Control Block 300 Series (38.1871.00), Handpiece Control Block 300 Series (38.1970.00) & Handpiece Control Block 500 Series (38.1972.00)

RPI Part #ADS280

OEM Part #38.1824.00 (Version 2)



STEM (WATER FLOW ADJ.)

- Version 2 design
- Material: Stainless steel
- Double flat stem
- Includes: (1 pc) O-Ring (RPI Part #RPO752)
- Also available: Knob (Flow Adj.) (RPI Part #ADK281), High Temp Lubricant (RPI Part #RPL090), and Version 1 Stem (Water Flow Adj.) (RPI Part #ADS233)

Fits: Handpiece Control Block Version 2 (3 and 4 Position) Models: 332, 333, 334, 335, 336 & 541

Models: Stems, Handpiece Control Block 300 Series (38.1871.00), Handpiece Control Block 300 Series (38.1970.00) & Handpiece Control Block 500 Series (38.1972.00)

RPI Part #ADK281

OEM Part #38.1823.00 (Version 2)



KNOB (FLOW ADJ.)

- Version 2 design
- 5 per package
- Material: Gray polyurethane
- Fits: Double flat stem
- Also available: Stem (Air Coolant) (RPI Part #ADS279) and Stem (Water Flow Adj.) (RPI Part #ADS280)

Fits: Handpiece Control Block Version 2 (3 and 4 Position) Models: 332, 333, 334, 335, 336 & 541

Models: Knobs, Handpiece Control Block 300 Series (38.1871.00), Handpiece Control Block 300 Series (38.1970.00) & Handpiece Control Block 500 Series (38.1972.00)



RPI Part #ADR283

OEM Part #38.1800.01

BARB RETAINER

- Version 1 (3 Position) and Version 2 (3 Position) compatible
- Material: White delrin

Fits: Handpiece Control Block Version 1 and 2 (3 Position) Models: 332, 333, 334, 335, and 336

Models: Barbs and Barb Kits, Handpiece Control Block 300 Series (38.1831.00), Handpiece Control Block 300 Series (38.1871.00)



RPI Part #ADD282

OEM Part #38.1833.01

DIAPHRAGM (3 POSITION)

- Version 1 (3 Position) and Version 2 (3 Position) compatible
- 5 per package
- Material: Fairprene

Fits: Handpiece Control Block Version 1 and 2 (3 Position) Models: 332, 333, 334, 335, and 336

Models: Diaphragms, Handpiece Control Block 300 Series (38.1831.00), Handpiece Control Block 300 Series (38.1871.00)

Don't Forget ...

RPI's offers other parts to fit A-dec 300 & 500 Handpiece Control Blocks including:

- Barb Retainer (4 Position)
- Diaphragm (4 Position)
- Control Block Kit (Major)
- Control Block Kit (Minor)
- Helical Compression Springs
- Many more!



replacement parts industries, inc. A Healthcare Components Group Company

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@ripiparts.com • WEBSITE www.ripiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. 2003AD RD1 Show 11/22

SECTION

4

PAGE

554

UPDATED

11/22